


Portfolio of Ophthalmic Surgical Devices



www.hamptonsurgical.co.uk

Table of Contents

Knives and Blades	5
Safety Knives	6
Conventional Knives	8
Conventional Knives	10
Mini-Blades	11
Griddles Mini-Blades	12
Blades	13
Cannulas	14
Anesthesia Cannulas	15
Anterior Chamber Cannulas	16
Hydrodissection Cannulas	17
Cystotome Cannulas	20
Irrigation Cannulas	22
Irrigation and Aspiration Cannulas	23
Polishers	24
Lacrimal Cannulas	26
Refractive and FLACS	27
Vitreoretinal Cannulas	28
Single-Use Instruments	29
Sprung Speculums	30
Pediatric Speculums	31
Fixation Forceps - Toothed	31
Fixation Forceps - Notched	32
Capsulorhexis Forceps	33
General Forceps	35
IOL Forceps	36
Tying Forceps	36
Choppers	37
Manipulators	38
Dissection Scissors	40
Dressing Scissors	40
Sprung Scissors	41
Dressing Forceps	42
Needle Holders and Drivers	42
Flap Elevators and Separators	43
Epithelium Spatulas, Hoes and Wells	45
Markers, Gauges and Calipers	46
Probes and Dilators	47
Clamps and Curettes	48
Fluid Control	
PVA Eye Spears	49
Full Body Drapes	51
Half Body Drapes	52
Mini Drapes	54
Accessories	55
Eye Shields and bandages	56
Custom Pack	57

Cataract Knives

Trusted Performance, Precision and Innovation

Hampton offers a wide range of conventional knives for cataract surgery. Slit, Sideport, Crescent and Micro-sharp knives come in a wide range of size to meet the needs of ophthalmic surgeons.



Knives & Blades

Conventional Knives

Slit Knives

Slit Knives are available in a wide range of sizes designed for smooth penetration and accurate incision width.



Single Bevel

Brand: HAMPTONS

Description	MODEL
1.8 mm Slit, 45 Degrees, Single bevel up	HSE-373718
2.2 mm Slit, 45 Degrees, Single bevel up	HSE-373722
2.4 mm Slit, 45 Degrees, Single bevel up	HSE-373724
2.5 mm Slit, 45 Degrees, Single bevel up	HSE-373725
2.6 mm Slit, 45 Degrees, Single bevel up	HSE-373726
2.75 mm Slit, 45 Degrees, Single bevel up	HSE-373727
2.8 mm Slit, Straight, Matte	HSE-373328
2.8 mm Slit, 45 Degrees, Single bevel up	HSE-373728
3 mm Slit, 45 Degrees, Single bevel up	HSE-373730
3.2 mm Slit, Straight, Matte	HSE-373332
3.2 mm Slit, 45 Degrees, Single bevel up	HSE-373732

Double Bevel

Brand: HAMPTONS

Description	MODEL
2 mm Slit, 45 Degrees, Double bevel up	HSE-370820
2.2 mm Slit, 45 Degrees, Double bevel up	HSE-370822
2.4 mm Slit, 45 Degrees, Double bevel up	HSE-370824
2.5 mm Slit, 45 Degrees, Double bevel up	HSE-370825
2.65 mm Slit, 45 Degrees, Double bevel up	HSE-370826
2.75 mm Slit, 45 Degrees, Double bevel up	HSE-370827
2.8 mm Slit, 45 Degrees, Double bevel up	HSE-370828
3 mm Slit, 45 Degrees, Double bevel up	HSE-370830

Crescent Knives

Crescent Knives offer the flexibility to create optimal tunnel architecture. Available in straight, bevel up or bevel down.



Brand: HAMPTONS

Description	MODEL
2.5 mm, 55 Degrees, Bevel up	HSE-373807
2.5 mm, 55 Degrees, Straight	HSE-373808
2.5 mm, 55 Degrees, Bevel down	HSE-373809

Sideport Knives

Sideport Knives are designed to create precise, consistent and fluid tight stab incisions.



Brand: HAMPTONS

Description	MODEL
1 mm Sideport Knife, Angled	HSE-373710
1.2 mm Sideport Knife, Angled	HSE-373712

Knives & Blades

Conventional Knives

EdgeAhead Sideport and MVR Knives



Sclerotome Knife

The Sclerotome knife has a smaller cutting edge than the Sclerotome blade (375700) and is designed to have a more precise, controlled dissection.



Brand:
HAMPTONS

Description	MODEL
Multi sided sharp all around	HSE-375701

Brand: HAMPTONS

Description	MODEL
0.6 mm Stiletto, 5/Bx	HSE-585232
0.6 mm MVR, Straight, 5/Bx	HSE-585231
1.1 mm MVR, Straight, 5/Bx	HSE-585240
1.15 mm MVR, 20G, Straight, 10/Bx	HSE-376630
1.15 mm Sideport, Angled, 10/Bx	HSE-376620

Mini-Blade Knife



Brand: HAMPTONS

Description	MODEL
Surgical blade, 180 Degree cutting surface	HSE-374769



Knives & Blades

Conventional Knives

Optimum Knives

Optimum Knives are color coded for identification. Offered in a wide range of sizes.



Brand: HAMPTONS

Description	MODEL
30 Degrees, 15 Degree bend left	HSE-374890
15 Degrees, Straight, Green	HSE-374891
22.5 Degrees, Straight, Slate blue	HSE-374892
30 Degrees, Straight, Black	HSE-374893
45 Degrees, Straight, Tan	HSE-374894

Micro-Unitome Knives

Micro-Unitome Knives are manufactured using a proprietary grindless technology to support a more consistent sharpness. Each knife is color coded for identification.



Brand: HAMPTONS

Description	MODEL
1.5 mm, 15 Degrees, Orange	HSE-377512
3 mm, 15 Degrees, Blue	HSE-377514
5 mm, 15 Degrees, Green	HSE-377516
4 mm, 22.5 Degrees, Light blue	HSE-377523
3.5 mm, 30 Degrees, Black	HSE-377531
3 mm, 45 Degrees, Tan	HSE-377546



Knives & Blades

Mini-Blades

Available in a variety of designs to meet all small incision and excision needs.

Mini-Blade



Brand: HAMPTONS

Description	MODEL
Sharp on two sides with pointed tip	HSE-376300

Mini-Blade



Brand: HAMPTONS

Description	MODEL
Pointed tip, Smaller alternative to #11 general blade	HSE-376500

Mini-Blade

Brand: HAMPTONS



Description	MODEL
MIS 6400, Narrower version of Beaver6400 for minimally invasive surgery	HSE-376401

Sclerotome Blade

The Sclerotome is a versatile blade suitable for glaucoma and retinal procedures.



Brand: HAMPTONS

Description	MODEL
Multi sided sharp all around S/Su	HSE-375700



Knives & Blades

Grindless HAMPTONS Mini-Blades

Hampton's Mini-Blades are manufactured using a proprietary grindless technology and are designed to meet all small incision and excision needs. Mini-Blades are preassembled into a unique, patented tip, designed for safeassembly. Each Mini-Blade is color coded for identification.

6400 Mini-Blade



Brand: HAMPTONS

Description	MODEL
Rounded tip, Sharp on one side, Green	HSE-6400

6700 Mini-Blade



Brand: HAMPTONS

Description	MODEL
Smaller alternative to #15 general blade, Orange	HSE-6700

6600 Mini-Blade



Brand: HAMPTONS

Description	MODEL
Lamellar, Angled 60 degrees, Smooth undersurface, Sharp all around, Bevel up, Black	HSE-6600

6900 Mini-Blade



Brand: HAMPTONS

Description	MODEL
Straight, Sharp all around, 180 Degree cutting surface, Double bevel blade, Blue	HSE-6900

6610 Mini-Blade



Brand: HAMPTONS

Description	MODEL
Lamellar, Straight, Smooth undersurface, Sharp all around, Black	HSE-6610

Knives & Blades

Blades

Micro-Sharp® Blades

Micro-Sharp Blades are manufactured using a proprietary grindless technology. Each blade is color coded for identification.



Brand: HAMPTONS

Description	MODEL
1.5 mm, 15 Degrees, Orange	HSE-377511
3 mm, 15 Degrees, Blue	HSE-377513
5 mm, 15 Degrees, Green	HSE-377515
4 mm, 22.5 Degrees, Light blue	HSE-377522
3.5 mm, 30 Degrees, Black	HSE-377530

Optimum Blades

Optimum stab blades are color coded for identification. Offered in a wide range of sizes.



Brand: HAMPTONS

Description	MODEL
15 Degrees, Straight, Green	HSE-374881
22.5 Degrees, Straight, Slate blue	HSE-374882
30 Degrees, Black	HSE-374883



- Full range of single-use cannulas developed to support all ophthalmic surgical procedures
- A history of innovation to meet your diverse surgical needs



Cannulas

Anesthesia Cannulas

Retrobulbar (Atkinson)

Intended for injecting anaesthetic agents into the muscle cone.



QTY: 5/bx
Brand: HAMPTONS

Description	MODEL
● 23G x 1-1/2 in, 0.6 x 38 mm	HSE-581638
● 25G x 1-1/2 in, 0.5 x 38 mm	HSE-581637

Retrobulbar/Peribulbar (Atkinson)



QTY: 5/bx
Brand: HAMPTONS

Description	MODEL
● 23G x 1-1/4 in, 0.6 x 32 mm	HSE-585001
● 25G x 1-1/4 in, 0.5 x 32 mm	HSE-585108
● 27G x 1-1/4 in, 0.4 x 32 mm	HSE-585109

Peribulbar (Atkinson)

Intended for injecting anaesthetic agents around the globe (periocular).



QTY: 5/bx
Brand: HAMPTONS

Description	MODEL
● 23G x 7/8 in, 0.6 x 22 mm	HSE-585144
● 25G x 7/8 in, 0.5 x 22 mm	HSE-585027
● 27G x 7/8 in, 0.4 x 22 mm	HSE-585145

Retrobulbar (Sharp Atkinson)



QTY: 5/bx
Brand: HAMPTONS

Description	MODEL
● 25G x 1-1/2 in, 0.5 x 38 mm	HSE-581641

Retrobulbar (Sharp Point)

Intended for injecting anaesthetic agents into the muscle cone.



QTY: 5/bx
Brand: HAMPTONS

Description	MODEL
● 25G x 1-1/2 in, 0.5 x 38 mm	HSE-581275

Retrobulbar, Curved (Atkinson/Straus)

Designed to course around the globe and enter the anterior region of the muscle cone inferotemporal and parallel to the optic nerve. Straight portion 11 mm, curved portion 22 mm.



QTY: 5/bx
Brand: HAMPTONS

Description	MODEL
● 25G x 1-3/8 in, 0.5 x 34 mm	HSE-585019

Cannulas

Anesthesia Cannulas

Sub-Tenon's Anesthesia Cannula

Intended for administering anaesthetic agents into the posterior sub-tenon's space.



Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
● 19G x 1 in, 1.1 x 25 mm	HSE-585176

Anterior Chamber Cannulas

Intended to be used for irrigation/injection of liquid.

Anterior Chamber (Rycroft)

Angled 45 degrees, 4 mm from end.



Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
● 21G x 7/8 in, 0.8 x 22 mm	HSE-581281
● 23G x 7/8 in, 0.6 x 22 mm	HSE-585025
● 25G x 7/8 in, 0.5 x 22 mm	HSE-581639
● 27G x 7/8 in, 0.4 x 22 mm	HSE-581280
●* 30G x 7/8 in, 0.3 x 22 mm	HSE-581273

Anterior Chamber Cannula

Angled 45 degrees, 6 mm from end.



Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
● 25G x 7/8 in, 0.5 x 22 mm	HSE-581642
● 27G x 7/8 in, 0.4 x 22 mm	HSE-581284
●* 30G x 7/8 in, 0.3 x 22 mm	HSE-581285

Cannulas

Hydrodissection Cannulas

Intended to separate the cortex from the capsule via irrigation (hydrodissection).

Nucleus Hydrodissector (Pearce)

Angled 35 degrees, 8 mm from end.



Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
● 25G x 7/8 in, 0.5 x 22 mm	HSE-585037
● 27G x 7/8 in, 0.4 x 22 mm with 8 mm bend	HSE-585036
● 27G x 7/8 in, 0.4 x 22 mm	HSE-585274

Nucleus Hydrodelineator (P. Koch)

Intended to break up lens contents and facilitate their aspiration via irrigation.



Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
● 25G x 7/8 in, 0.5 x 22 mm	HSE-585015

Cortical Cleaving Hydrodissector (Flat End)

45 degrees, 11 mm from end. For cleaving cortex from capsular bag.



Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
● 27G 7/8 in, 0.4 x 22 mm	HSE-585158

Nucleus Hydrodissector (Sinskey)

35 degrees, 8 mm from end.



Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
● 25G x 7/8 in, 0.5 x 22 mm	HSE-585195

Nucleus Hydrodissector

45 degrees, 7 mm from end.



Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
● 30G x 7/8 in, 0.3 x 22 mm	HSE-581592

Nucleus Hydrodissector (Sauter)

Tip is flattened on the vertical plane.



Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
● 27G x 7/8 in, 0.4 x 22 mm	HSE-585099

Cannulas

Hydrodissection Cannulas

Intended to separate the cortex (tissue) from the capsule via irrigation (hydrodissection).

Nucleus Hydrodissector

Tip is flattened on the horizontal plane.



Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
● 27G 7/8 in, 0.4 x 22 mm	HSE-585155

Nucleus Hydrodissector (Chang)

90 Degree bend, 1 mm flattened tip.



Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
● 27G x 7/10 in, 0.4 x 18 mm, Bevel angle 33 degrees	HSE-585294
● 27G x 7/10 in, 0.4 x 18 mm, Bevel angle 47 degrees	HSE-585295

Tapered Hydrodelineator (Blumenthal)

Tip is tapered to 0.3 mm (30G), 45 degrees, 4 mm from end.



Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
● 25G x 7/8 in, .5 x 22 mm	HSE-585107

Cortical Cleaving Hydrodissector (Round End)

45 degrees, 11 mm from end.



Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
● 27G x 7/8 in, 0.4 x 22 mm	HSE-585157

Nucleus Hydrodissector (Buratto)

Formed, 1 mm bent tip.



Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
● 25G x 7/8 in, 0.5 x 22 mm	HSE-585199

Cortical Cleaving Hydrodissector (Bolger)

Tip is tapered to 0.3 mm (30G), 33 degrees, 1 mm from end.



Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
● 25G x 7/8 in, 0.5 x 22 mm	HSE-585159

Cannulas

Hydrodissection Cannulas

Intended to separate the cortex (tissue) from the capsule via irrigation (hydrodissection).

Hydrodissection Cannula (Akahoshi Style)

Tapered tip, Thin wall, Angled 1 mm from end.



Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
● 27G 7/8 in, 0.4 x 22 mm	HSE-585198

Micro J Shape Cannula

0.25 mm open end port.
Hook is 1.5 mm long x 1.75 mm wide.



Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
● 25G x 7/8 in, 0.5 x 22 mm	HSE-585003
● 27G x 7/8 in, 0.4 x 22 mm	HSE-585018

Micro J Shape Cannula (Hessburg)

Hook is 0.85 mm long x 1.7 mm wide.



Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
● 25G x 3/4 in, 0.5 x 20 mm	HSE-585148

Micro J Shape Hydrodissector (Pearce)

Hook is 2.25 mm long x 1.75 mm wide.



Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
● 27G x 7/8 in, 0.4 x 22 mm	HSE-585029

Cortical Cleaving Hydrodissector (Grabow)

40 degrees, 12 mm from end.



585197

Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
● 30G x 5/8 in, 0.3 x 16 mm	HSE-585197

HydroExpression Cannula (Corydon)

1.8 mm from end.



585191

Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
● 23G x 7/8 in, 0.6 x 22 mm	HSE-585091
● 25G x 7/8 in, 0.5 x 22 mm	HSE-585090
● 27G x 7/8 in, 0.4 x 22 mm	HSE-585191

Cannulas

Cystotome Cannulas

Intended to perform a capsulorhexis. Indicated for use in cataract surgery.

Irrigating Cystotome



Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
● 23G x 5/8 in, 0.6 x 16 mm	HSE-581599
● 25G x 5/8 in, 0.5 x 16 mm	HSE-581610
● 27G x 5/8 in, 0.4 x 16 mm	HSE-581618
● 30G x 5/8 in, 0.3 x 16 mm	HSE-581630

Irrigating Cystotome (Blumenthal)



Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
● 27G x 5/8 in, 0.4 x 16 mm	HSE-585043

Irrigating Cystotome

Formed 8.5 mm from bend to tip. Shortened curve designed for deep set eyes or small pupils.



Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
● 25G 5/8 in, 0.5 x 16 mm	HSE-585004
● 27G 5/8 in, 0.4 x 16 mm	HSE-585008

Irrigating Formed Cystotome

Tip 45 degrees.



Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
● 25G x 5/8 in, 0.5 x 16 mm	HSE-585049

Irrigating Cystotome



Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
● 23G x 5/8 in, 0.6 x 16 mm	HSE-581594
● 25G x 5/8 in, 0.5 x 16 mm	HSE-581617

Irrigating Side Cutting Cystotome



Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
● 25G x 5/8 in, 0.5 x 16 mm	HSE-581602

Cannulas

Cystotome Cannulas

Intended to perform a capsulorhexis. Indicated for use in cataract surgery.

Irrigating Cystotome (Blumenthal)

Formed 8.5 mm from bend to tip.



Qty: 5/Bx
Brand:
HAMPTONS

Description	MODEL
● 25G x 5/8 in, 0.5 x 16 mm	HSE-585113
● 27G x 5/8 in, 0.4 x 16 mm	HSE-585042

Irrigating Cystotome (Simcoe)

45 Degrees and 10.5 mm from tip.



Qty: 5/Bx
Brand:
HAMPTONS

Description	MODEL
● 27G x 1/2 in, 0.4 x 12.5 mm	HSE-585093

Blunt Tip Irrigating Cystotome (Kellan)

Tip angled toward bevel. Unique blunt tip creates a wide capsular flap. Useful for pear shape capsulotomies.



Qty: 5/Bx
Brand:
HAMPTONS

Description	MODEL
● 25G x 7/8 in, 0.5 x 22 mm	HSE-585115

Irrigating Cystotome (Joyce)

Tip 0.5 mm long.



Qty: 5/Bx
Brand:
HAMPTONS

Description	MODEL
● 25G x 5/8 in, 0.5 x 16 mm	HSE-585030
● 27G x 5/8 in, 0.4 x 16 mm	HSE-585039

Irrigating Cystotome (Reverse)

Tip 90 degrees toward bevel.



Qty: 5/Bx
Brand :HAMPTONS

Description	MODEL
● 25G x 5/8 in, 0.5 x 16 mm	HSE-585034

Cannulas

Irrigation Cannulas

Intended for irrigation of liquids and body fluids.

Bishop Harmon Cannula

Intended for irrigation of liquids and body fluids. 40 degrees, 8 mm from end.



Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
● 19G x 1 in, 1.1 x 25 mm	HSE-585040

Anterior Chamber Cannula

30 degrees, 6 mm from end.



Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
○ 16G x 1 in, 1.6 x 25 mm	HSE-581277

Irrigation and Aspiration Cannulas

Intended to be used for irrigation and/or aspiration of liquid and air.

Viscoflow Formed

Formed 11 mm from end, radius is 19 mm.



Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
● 27G x 7/8 in, 0.4 x 22 mm	HSE-585007

Viscoflow Cannula

Angled 45 Degrees, 9 mm from end.



Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
● 23G x 7/8 in, 0.6 x 22 mm	HSE-585009
● 25G x 7/8 in, 0.5 x 22 mm	HSE-585045
● 27G x 7/8 in, 0.4 x 22 mm	HSE-585026
● 30G x 7/8 in, 0.3 x 22 mm	HSE-585046

Single-Use Instruments

Flap Elevators and Separators

LASIK Flap Separator/Elevator (Seibel)



Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
0.16 mm x 0.28 mm x 1 mm Tip, 20 mm Straight shaft, 0.5 mm x 10.5 mm pointed hook, 52 mm Handle	HSE-1172S

LASIK Flap Separator and Elevator, Binder



Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
0.1 mm x 0.3 mm Rounded tip 0.7 mm long, 10 mm Curved 60 degree angled shaft, 75 mm Handle	HSE-1206S

LASIK Flap Elevator Curved 8 mm



Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
0.2 mm x 0.15 mm Tip, Curved 8 mm tapered shaft at 90 degrees, 75 mm Handle	HSE-1177S

LASIK Flap Separator/Elevator Curved



Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
0.15 mm x 0.3 mm x 8 mm Curved separator at 90 degree to shaft, 0.15 mm x 0.3 mm x 1 mm Pocketing hook, 52 mm Handle	HSE-1294S

Y Shape Corneal Marker



Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
2 mm x 8 mm Marking blades, 75 mm Handle	HSE-1134S

Reinstein Lenticule Separator



Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
0.22 mm x 2 mm Pocketing hook, 1.2 mm x 2 mm Spatula tip, Curved separator at 90 degrees to shaft, 52 mm Handle	HSE-1297S

Single-Use Instruments

Flap Elevators and Separators

LASIK Flap Separator/
Elevator Straight



Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
0.15 mm x 0.3 mm x 8 mm Straight separator at 90 degree to shaft, 0.15 mm x 0.3 mm x 1 mm Pocketing hook, 52 mm Handle	HSE-12955

LASIK Flap Roller

Used after realignment. Acts as a squeegee to remove excess fluid from the interface; smooths flap and removes folds.



Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
14 mm sleeve	HSE-585233



Single-Use Instruments

Epithelium Spatulas, Hoes and Wells

Corneal Epithelium Spatula Flat



Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
5 mm Spatula with 0.9 mm flat blade, 0.8 mm x 1.7 mm Round bladed tip, 75 mm Handle	HSE-1113S

Epithelium Well



Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
3 mm x 8.5 mm ID well, Angled 45 degrees, 75 mm Handle	HSE-1197S
3 mm x 9 mm ID well, Angled 45 degrees, 75 mm Handle	HSE-1196S

Corneal Epithelium Spatula Round



Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
2.65 mm x 5.5 mm x 0.35 mm Rounded blade, 75 mm Handle	HSE-1187S

Spatula PRK



Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
Spatula tip and top are slightly sharpened to facilitate more aggressive scraping.	HSE-585213



Markers, Gauges and Calipers

Mendez Gauge



Qty: 5/Bx

Brand: HAMPTONS

Description	MODEL
2.5 mm Polycarbonate ring with 13 mm ID, 0-180 Degree marks in increments of 10, 75 mm Handle	HSE-1152S

Double Ended IVT Guide Pad



Qty: 5/Bx

Brand: HAMPTONS

Description	MODEL
7 mm Pad with 3.5 guide, 0.8 mm Pad with 4 mm guide, 6 mm 45 Degree angled shafts, 90 mm Long	HSE-1285S

Pointed IVT Marker 3.5/4, Blue



Qty: 5/Bx

Brand: HAMPTONS

Description	MODEL
2 x 0.55 mm Tips with 3.5 mm spread, 2 x 0.65 mm Tips with 4 mm spread, Polycarbonate, 108 mm Long	HSE-1277S

Single-Use Instruments

Probes and Dilators

Punctal Dilator and Probe



Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
25 mm (1 in) medium taper, Blunt, Wilder style, dilator mated with a size 0 Bowman lacrimal probe, 0.7 mm x 45 mm (22G x 1 3/4 in)	HSE-585124

Golf Club Spud



Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
0.1 mm x 0.8 mm x 2 mm Oval profile, 18 mm Straight shaft, 75 mm Handle	HSE-14175

Punctum Dilator,
Fine/Medium



Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
0.4 mm x 21.5 mm Tip, 0.6 mm x 18 mm Tip, 77 mm Handle	HSE-11815

Punctum Dilator/
Lacrimal Probe



Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
0.28 mm x 25 mm Dilator, 0.7 mm x 44.5 mm Probe, 70 mm Handle	HSE-14195

Punctum Dilator/Seeker



Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
0.45 mm x 21.5 mm Tip, 5 mm Pointed tip, 18 mm Straight shaft, 77 mm Handle	HSE-11935



Single-Use Instruments

Clamps and Curettes

Chalazion Clamp, Oval



Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
13 mm x 20 mm Oval profile, 10 mm x 17 mm Aperture, 90 mm Long	HSE-1149S

Chalazion Curette, Double, 1 mm/2 mm



Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
1 mm Curette, 1 mm Depth, 2 mm Curette, 1.35 mm Depth, 128 mm Long	HSE-1157S

Chalazion Clamp, Round



Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
12 mm Round profile, 9 mm Round aperture, 90 mm Long	HSE-1150S

Chalazion Curette, Double, 2 mm/3 mm



Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
2 mm Curette, 1.35 mm Depth, 3 mm Curette, 1.6 mm Depth, 128 mm Long	HSE-1158S

Fixation Rings


Swivel Fixation Ring



Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
13 mm Diameter	HSE-372242

Fluid Control



Fluid Control

PVA Eye Spears

Spear overall handle length is 7.16 cm or 2.82 in.

UltraClean™ LASIK PVA Eye Spears

Compressed shape, designed for fast expanding. Made from highly absorbent, virtually lint free PVA material. Designed specifically for rapid recovery of fluids during LASIK and other corneal procedures.



Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
UltraClean™ Eye Spears	HSE-40810
UltraClean™ Eye Spears	HSE-40815

Visitec® LASIK PVA Expanded Eye Spear

Made from highly absorbent PVA material. Pre-expanded for convenience. Designed to provide instantaneous wicking. No need to pre-soak and fan.




Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
LASIK Expanded Spear	HSE-581709

Visitec® LASIK PVA Eye Spears

Made from highly absorbent PVA material. Designed for optimal tissue manipulation and controlled fluid absorption to facilitate flap alignment and adhesion.



Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
LASIK Eye Spear - Orange	HSE-581761

Weck-Cel® Cellulose Eye Spears

Made from highly absorbent, natural, cellulose material. Fast wicking action delivers quick absorption of blood and fluid from the surgical area. The flat tip provides precise application at the surgical incision.



Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
Weck-Cel® Spears, 6 Pack, 144/Bx	HSE-0008680
Weck-Cel® Spears, 10 Pack, 180/Bx	HSE-0008685

Visitec® Cellulose Eye Spears

Made from highly absorbent, natural, cellulose material. Fast wicking action delivers quick absorption of fluid from the surgical area. The flat tip provides precise application at the surgical incision. Magenta handle.



Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
Sponge Eye, 7 cm (2.75 in) overall length.	HSE-581089

Fluid Control

Corneal Light Shields

Visitec® Corneal Light Shield

PVA absorbent pad. Soft, porous structure designed to protect the retina from operating microscope light and contains no cellulose. Upon moistening, shield becomes soft and conforms to shape of cornea. Designed to keep cornea moist.



Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
Visitec® Corneal Light Shield, 8 mm (0.3 in) diameter	HSE-581062



Fluid Control

Full Body Drapes

Visidrape® full body drapes feature wider widths and a three-sided channel bag for maximum fluid capacity. Full body drapes can be utilized in all ophthalmic procedures.

Full Body Incise Drape

Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
Visidrape® Non-Woven Incise Drape with Fluid Collection Bag	HSE-581087
Overall size:	140 cm x 244 cm (55 in x 96 in)
Incise area:	9 cm x 11.5 cm (3.5 in x 4.5 in)
Fluid collection bag:	25.4 cm x 15.2 cm (10 in x 5.98 in)
Capacity:	500 ml
Material Type:	Non-Woven Wetlaid



Full Body Aperture Drape

Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
Visidrape® Full Body Aperture Drape	HSE-581088
Overall size:	140 cm x 244 cm (55 in x 96 in)
Aperture area:	5 cm x 7.5 cm (2 in x 3 in)
Capacity:	500 ml
Material Type:	Non-Woven Wetlaid



Full Body Incise Drape with Plastic Head and Channel Bag

Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
Visidrape® Full Body Incise Drape with Plastic Head and Channel Bag	HSE-581715
Overall size:	165 cm x 254 cm (65 in x 100 in)
Incise Area:	9 cm x 11 cm (3.5 in x 4 in)
Channel Bag:	62 cm x 40 cm (24 in x 15.75 in)
Capacity:	1000 ml
Material Type:	SMS (63 in) and Opaque (37 in)



Fluid Control

Half Body Drapes

Embossed Plastic Incise Drape with Fluid Collection Bag

Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
Visidrape® Embossed Plastic Incise Drape with Fluid Collection Bag	HSE-581071
Overall size:	101.5 cm x 127 cm (40 in x 50 in)
Incise area:	9 cm x 11.5 cm (3.5 in x 4.5 in)
Fluid collection bag:	15.2 cm x 16.5 cm (6 in x 6.5 in)
Capacity:	250 ml
Material Type:	Polyethylene Opaque Film



Embossed Plastic Aperture Drape with Fluid Collection Bag

Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
Visidrape® Embossed Plastic Aperture Drape with Fluid Collection Bag	HSE-581078
Overall size:	101.5 cm x 127 cm (40 in x 50 in)
Aperture area:	5 cm x 7.5 cm (2 in x 3 in)
Fluid collection bag:	15.2 cm x 16.5 cm (6 in x 6 1/2 in)
Capacity:	250 ml
Material Type:	Polyethylene Opaque Film



Non-Woven Incise Drape with Fluid Collection Bag

Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
Visidrape® Non-Woven Incise Drape with Fluid Collection Bag	HSE-581072
Overall size:	101.5 cm x 127 cm (40 in x 50 in)
Incise area:	9 cm x 11.5 cm (3.5 in x 4.5 in)
Fluid collection bag:	15.2 cm x 16.5 cm (6 in x 6.5 in)
Capacity:	250 ml
Material Type:	Non-Woven Wetlaid



Fluid Control

Half Body Drapes

Non-Woven Aperture Drape with Fluid Collection Bag

Qty: 5/Bx
Brand: HAMPTONS


Description	MODEL
Visidrape® Non-Woven Aperture Drape with Fluid Collection Bag	HSE-581079
Overall size:	101.5 cm x 127 cm (40 in x 50 in)
Aperture area:	5 cm x 7.5 cm (2 in x 3 in)
Fluid collection bag:	15.2 cm x 16.5 cm (6 in x 6.5 in)
Capacity:	250 ml
Material Type:	Non-Woven Wetlaid



Non-Woven Incise Drape with Channel Bag

Qty: 5/Bx
Brand: HAMPTONS


Description	MODEL
Visidrape® Non-Woven Incise Drape with Channel Bag	HSE-581084
Overall size:	140 cm x 127 cm (55 in x 50 in)
Incise area:	9 cm x 11.5 cm (3.5 in x 4.5 in)
Fluid collection bag with channels:	14.6 cm x 25.4 cm (5.75 in x 10 in)
Capacity:	500 ml
Material Type:	Non-Woven Wetlaid



Non-Woven Aperture Drape with Channel Bag

Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
Visidrape® Non-Woven Aperture Drape with Channel Bag	HSE-581086
Overall size:	140 cm x 127 cm (55 in x 50 in)
Aperture area:	5 cm x 7.5 cm (2 in x 3 in)
Fluid collection bag with channels:	14.6 cm x 25.4 cm (5.75 in x 10 in)
Capacity:	500 ml
Material Type:	Non-Woven Wetlaid



Fluid Control

Half Body Drapes

Non-Woven Half Body Incise Drape* with Channel Bag

Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
Visidrape® Non-Woven Half Body Incise Drape* with Channel Bag	HSE-581723
Overall size:	165 cm x 183 cm (65 in x 72 in)
Incise area:	9 cm x 11 cm (3.5 in x 4 in)
Channel bag:	62 cm x 40 cm (24 in x 15.75 in)
Capacity:	1000 ml
Material Type:	SMS

*Without wick



Mini Drapes

Ideal for refractive or office procedures.

Mini Plastic Incise Drape

Adhesive backing.

Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
Visidrape® Mini Plastic Incise Drape	HSE-581095
Overall size:	41 cm x 41 cm (16 in x 16 in)
Incise area:	5 cm x 7.5 cm (2 in x 3 in)
Material Type:	Polyethylene Opaque Film



Mini Plastic Aperture Drape

Adhesive backing.

Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
Visidrape® Mini Plastic Aperture Drape	HSE-581096
Overall size:	41 cm x 41 cm (16 in x 16 in)
Aperture area:	5 cm x 7.5 cm (2 in x 3 in)
Material Type:	Polyethylene Opaque Film



Fluid Control


Mini Drapes

Mini Plastic Aperture Drape

Small incise provides complete coverage of lashes, lids and brows. Converts aperture to incise drape.

Qty: 5/Bx
Brand: HAMPTONS

Description		MODEL
Visidrape® Incise Drape		HSE-581700
Overall size:	10.2 cm x 20.4 cm (4 in x 8 in)	
Adhesive area:	11 cm x 16 cm (4 in x 6.25 in)	
Material Type:	Polyethylene Opaque Film	



Dual Refractive Drape

Qty: 5/Bx
Brand: HAMPTONS

Description		MODEL
Visidrape® Dual Refractive Drape		HSE-581710
Overall size:	40 cm x 26 cm (16 in x 10.25 in)	
Aperture area:	2.5 cm x 3.75 cm (1 in x 1.5 in)	
Material Type:	Polyethylene Opaque Film	




Accessories
Quality Accessories to Support
You and Your Patients Needs



Eye Shields and Bandages

Eye Shield (Left/Right)

White lightweight plastic eye shield with ventilation holes. Flexible eye shield, shaped to fit the right or left eye for protection. Smooth edges do not require a cover. Packaged sterile.



Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
Eye Shield, Left/Right 8.5 cm (3.3 in)	HSE-581019

Universal Pediatric Eye Shield

Clear lightweight plastic and features smooth edges and ventilation holes. For use on left or right eye. Packaged sterile.



Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
Pediatric Eye Shield	HSE-581708

Universal Eye Shield

Clear lightweight plastic, for use on left or right eye, that features smooth edges and ventilation holes. Packaged sterile.



Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
Universal Eye Shield, 7.5 cm x 6.5 cm (3 in x 2.5 in)	HSE-581040

Visitec® Guibor Eye Bandage

Disposable self-adhesive eye dressing designed to protect the eye during drying conditions, such as "dry eye". Straightforward to apply. Packaged sterile.



Qty: 5/Bx
Brand: HAMPTONS

Description	MODEL
Eye Bandage, Guibor Bubble Bandage	HSE-8590300

Our Cataract Pack

Hamptons Cataract Pack

Hamptons Cataract Pack

Our Cataract is designed according to research of surgeon preference world wide and recommended by top Surgeon World



Product	Description	Quantity
Swab	Non-Woolen (10x10cm 4 ply)	10 pcs
Tray	Half Shallow (Plastic)	01 pcs
Gallipot	Polypropylene (120ML)	02 pcs
Kidney Dish	Polypropylene (800ML)	01 pcs
Needle Blunt Fill Filter	Size (18g/ 19g x1.5")	01 pcs
Syringe	5 ML Luer lock	01 pcs
Syringe	3 ML Luer lock	03 pcs
Syringe	1 ML Luer lock	01 pcs
Scissors Conjunctival	Stainless Steel Single use	01 pcs
Scissors West-Cott	Curved Round Tip (55)	01 pcs
Drysdale Nucleus Spatula	Paddle Shape(Stainless Steel Single-use)	01 pcs
Speculum Lid	Stainless Steel Single use	01 pcs
Eye Shield	Clear (Plastic)	01 pcs
Drape Trolley Cover	Size (140cm x 190cm)	01 pcs
Paper Crepe white	Size (90cm x 90cm)	01 pcs
Drape Ophthalmic	Size(140cmx127 cm)with Incise& Pouch	01 pcs
Forcep Sponge Holding	Plastic Single Use (19cm 7.5")	01 pcs
Hoskins Forcep	Stainless Steel Single use	01 pcs
Hydrodissector (Akaoshi)	Angled needle (27G x7/8 0.4 x 22mm)	01 pcs
VIS Cannula Anterior	Chamber Rycroft (27G x 7/8")	01 pcs
Cannula Cystotome	Irrigating Blumenthal (27G x5/8")	01pcs
Slit Knife	Single bevel up angled (2.2mm)	01 pcs
MVR Knife	Side-port angled (1.0mm)	01 pcs
Total		37 pcs

Customize Cataract Pack

Customize Cataract Pack

Customize Cataract Pack

Standard Cataract Pack is been designed with nesscary items , In this surgeon can select any instrument needed



Product	Description	Quantity
Swab	Non-Woolen (10x10cm 4 ply)	10 pcs
Tray	Half Shallow (Plastic)	01 pcs
Gallipot	Polypropylene (120ML)	02 pcs
Kidney Dish	Polypropylene (800ML)	01 pcs
Needle Blunt Fill Filter	Size (18g/ 19g x1.5")	01 pcs
Syringe	5 ML Luer lock	01 pcs
Syringe	3 ML Luer lock	03 pcs
Syringe	1 ML Luer lock	01 pcs
Eye Shield	Clear (Plastic)	01 pcs
Drape Trolley Cover	Size (140cm x 190cm)	01 pcs
Paper Crepe white	Size (90cm x 90cm)	01 pcs
Drape Ophthalmic	Size(140cmx127 cm)with Incise& Pouch	01 pcs
Forcep Sponge Holding	Plastic Single Use (19cm 7.5")	01 pcs
If needed : can be provided on demand		