



MT MEDICAL

LED OPERATION LAMP

YDE700/500

A Type



- Streamlined, ultra-thin, enclosed head with reduced air resistance to minimize contamination
- Manual /Automatic focus
- Digital controlled light intensity with 10 grades of brightness and automatic memory for added flexibility
- Built in power supply systems avoiding complicated installation of wall mounted control boxes
- Auxiliary light source included for laparoscopic operations
- Rigid suspension arm with easy mobility and stable positioning
- No infrared or ultraviolet radiation

PARAMETERS:

Diameters:	700mm	500mm
Illuminance	160000Lux.	140000Lux.
Color Temperature	3600--5500K(±500K)variable	
Color Rendering Index	≥ 95	≥ 95
Ultraviolet light irradiance within400nm wave length:	≤0.006W/m2	
Power Of Bulb(W)	1W X N60	1W X N45
50% Spot Diameter d50(mm)	150	120
10% Spot Diameter d10(mm)	280	230
Depth Of Illuminance (cm)	140	120
The Adjusting Range Of Illuminance	25-100%	
Focus System	Electric / Manual (optional)	
Voltage	220V	220V
Frequency	50Hz	50Hz
Input Power	100W	80W
LED Service Life	100,000h	
The Temperature On Surgeon's Head Is Less Than	1°C	
Camera System	Optional	
Battery	Optional	

YDE500



LED OPERATION LAMP



YDE700/500

B Type


- Streamlined, ultra-thin, enclosed head with reduced air resistance to minimize contamination
- Manual /Automatic focus
- Digital controlled light intensity with 10 grades of brightness and automatic memory for added flexibility
- Built in power supply systems avoiding complicated installation of wall mounted control boxes
- Auxiliary light source included for laparoscopic operations
- Rigid suspension arm with easy mobility and stable positioning
- No infrared or ultraviolet radiation

PARAMETERS:

Diameters:	700mm	500mm
Illuminance	160000Lux.	140000Lux.
Color Temperature	3600--5500K(±500K)variable	
Color Rendering Index	≥ 95	≥ 95
Ultraviolet light irradiance within400nm wave length:	≤0.006W/m ²	
Power Of Bulb(W)	1W X N60	1W X N45
50% Spot Diameter d50(mm)	150	120
10% Spot Diameter d10(mm)	280	230
Depth Of Illuminance (cm)	140	120
The Adjusting Range Of Illuminance	25-100%	
Focus System	Electric / Manual (optional)	
Voltage	220V	220V
Frequency	50Hz	50Hz
Input Power	100W	80W
LED Service Life	100,000h	
The Temperature On Surgeon's Head Is Less Than	1°C	
Camera System	Optional	
Battery	Optional	

YDE500D





YDE700/500

Slim Design

- Utilising the most advanced LED technology available, the MT-medical LED70 series of operating theatre lights have been designed to provide a high quality, shadowless beam with extremely low running costs.
- The LED operation lamp has the backlight compensation function, one parts of LEDs over surgeon's head will be closed and Other Light source will be enhanced.
- LEDs (light emitting diodes) have been around for a long time. Recently this technology has become leading -edge in the field of lighting - a breakthrough technological solution to the inherent problems of the incandescent light bulb.
- We have taken the technology and developed it further to provide high definition lighting solutions with no infrared, low environmental impact, high brightness, low power consumption, no toxic chemicals and a long operating life.

PARAMETERS:

Diameters:	700mm	500mm
Illuminance	160000Lux.	140000Lux.
Color Temperature	3600--5500K(±500K)variable	
Color Rendering Index	≥97	≥97
R9 : 95	R13 : 92	
Power Of Bulb(W)	1W X N128	1W X N64
Spot Diameter d50(mm)	130-350	
Depth Of Illuminance (cm)	140	120
The Adjusting Range Of Illuminance	10-100%	
Focus System	Manual	
Minium Installation Height is	2800mm	
Voltage	220V	220V
Frequency	50Hz	50Hz
Input Power	160W	80W
LED Service Life	100,000h	
The Temperature On Surgeon's Head Is Less Than	1°C	
Camera System	Optional	
Battery	Optional	



YDE500

Slim Design



YDZ700/500

- shadowless operating lamps meet the requirements of various surgical operations. It is an ideal lighting facility for surgical operations.
- Special heat-isolated radiation filter reduces heat to the minimum level.
- Imported aluminum sheet squamulose reflector and transparent panel are adopted. Average service life of imported lamp receptacle and bulb can reach 1500 hours.
- The imported aluminum reflecting mirror makes light softer and more well-proportioned and increasing safety. CPU is adopted for the control, capable of increasing or decreasing and boasting brightness memory. When the main bulb malfunctions, the standby bulb will come into operation within 1/2 second. The major and auxiliary bulbs can be replaced easily. The dischargeable handle sheath can be sterilized in steam sterilizer.

PARAMETERS:

Diameters:	700mm	500mm
Illuminance	150000Lux.	120000Lux.
Color Temperature	3000K---6700K	
Color Rendering Index	≥95	
Spot Diameter	120--320mm	
Depth Of Illuminance	≥800mm	
Light-adjusting mode	10 grade dimmer system	
Focus System	Manual	
Temperature At Surgeon's Head	≤2°C	
Temperture At Surgical Field	≤8°C	
Working Distance:	55-170cm	
Power Supply	220V±22V 50Hz±1Hz	
Input Power	200W	
Power Of Bulb	24V/150W	
Bulb average Life	1500h	
Shift Of Major And Stand-by Lamps	YES	
Camera System	Optional	
Battery	Optional	



YDZ700



YDZ500D



MT MEDICAL

LED OPERATION LAMP

YDE700/700

C Type



- Digital controlled light intensity with 10 grades of brightness and automatic memory for added flexibility
- Built in power supply systems avoiding complicated installation of wall mounted control boxes
- Auxiliary light source included for laparoscopic operations
- Rigid suspension arm with easy mobility and stable positioning
- No infrared or ultraviolet radiation

PARAMETERS:

Diameters:	700mm	700mm
Illuminance	160000Lux.	160000Lux.
Color Temperature	3600--5500K(±500K)variable	
Color Rendering Index	≥ 95	≥ 95
R9 : 95 R13 : 92		
Power Of Bulb(W)	1W X N80	1W X N80
50% Spot Diameter d50(mm)	150	150
10% Spot Diameter d10(mm)	280	280
Depth Of Illuminance (cm)	140	140
The Adjusting Range Of Illuminance	25-100%	
Focus System	Electric / Manual (optional)	
Voltage	220V	220V
Frequency	50Hz	50Hz
Min.installation height	2800mm	
LED Service Life	100,000h	
The Temperature On Surgeon's Head	Is Less Than 1°C	
Camera System	Optional	
Battery	Optional	



YDE500 D



YDE700 /700 D Type

- Streamlined, ultra-thin, enclosed head with reduced air resistance to minimize contamination
- Manual /Automatic focus
- Digital controlled light intensity with 10 grades of brightness and automatic memory for added flexibility
- Built in power supply systems avoiding complicated installation of wall mounted control boxes
- Auxiliary light source included for laparoscopic operations
- Rigid suspension arm with easy mobility and stable positioning
- No infrared or ultraviolet radiation

PARAMETERS:

Diameters:	700mm	700mm
Illuminance	160000Lux.	160000Lux.
Color Temperature	3600--5500K(±500K)variable	
Color Rendering Index	≥95	≥95
R9 : 95 R13 : 92		
Power Of Bulb(W)	1W X N64	1W X N64
50% Spot Diameter d50(mm)	150	150
10% Spot Diameter d10(mm)	280	280
Depth Of Illuminance (cm)	140	140
The Adjusting Range Of Illuminance	25-100%	
Focus System	Electric / Manual (optional)	
Voltage	220V	220V
Frequency	50Hz	50Hz
Min.installation height	2800mm	
LED Service Life	100,000h	
The Temperature On Surgeon's Head Is Less Than 1°C		
Camera System	Optional	
Battery	Optional	



YDE700 E Type



YDE300

Wall / Ceiling mounted type

- YDE300 Portable LED examination light is an invaluable and cost effective tool. It has been specifically designed for minor procedures under local anaesthetic in health centres, GP surgeries and outpatient departments.
- FEATURE :
 - The use of multiple LED's provides a high, and even light intensity throughout the entire viewing area to reduce eye strain
 - A concise and self-explanatory 3 button control system on the lightthead itself
 - Balanced arms for one touch adjustment and a great range of movement
 - Ideal for minor surgery and gynaecological procedures
 - Removeable, Autoclavable Handle

PARAMETERS:

Diameters:	350mm	
Illuminance	80000Lux	
Color Temperature	4500K(±500K)variable	
Color Rendering Tndex (Ra)	≥90	
Input Power	50 W	
Power Of Bulb(W)	1W X N20	
Spot Diameter d50(mm)	80	
Depth Of Illuminance (cm)	80	
Battery	Optional	
Focus System	Automatic	
LED Service Life	100,000h	
Voltage	220V	220V
Frequency	50Hz	50Hz
Packing Details	0.3 CBM/unit, G.W.:42 KGS	



YDE300

Portable Type

YDE200D

Deep Irradiation Exam Light

Lamp head diameter(mm)	180
Illuminance(Lux)	≥50,000
Color temperature(K)	4,000K±500
Color rendering index(Ra.%)	≥88
Spot diameter(mm)	160~220
Depth of illuminance(mm)	600
Power of bulb(W)	50
Bulb service life(hours)	1,000
Voltage	220V 50Hz or 110V 60Hz
Backup battery	optional
Packing details	11 CBM/set, G.W.:20KGS



YDE180D

Portable LED Exam Light

Illumination	50000Lux
Color temperature	5500±500K
Voltage	AC90-240V
Bulb Power	LED 7*3W
LED service life	50000h
Power	7*3 W LED
Packing details	0.03 CBM/set,G.W.:8.5KGS



Hole Type Exam Light

YDK04D

Illuminance	≥40,000Lux
Color temperature	3,700K±500K
Color rendering index(Ra.%)	≥85
Spot diameter(mm)	160~220
Depth of illuminance(mm)	900
Power of bulb(W)	25
Bulb service life(hours)	1,000
Voltage	220V 50Hz
Backup battery	optional

YDE200 LED Portable LED exam light

Illuminance	60000Lux.
Color Temperature	4800 ±300 K
Power supply voltage	AC 220V/50-60Hz / optional 110V/60Hz
Power of bulb	3.2V/1W
Diameter of light head	180mm



YDE200C Ceiling LED exam light


Illuminance	40000Lux.
Color Temperature	5000±500 K
Input Power:	36 W
Voltage	AC90~240V
LED Service Life	50,000h



YDE-A LED exam light

YDE-B LED exam light

YDE10 LED exam light



DST-700G

Imaging Orthopedic Operation Table

DST-700G Electric orthopedic navigation image operation table, for the expansion of minimally invasive treatment to improve surgical accuracy with modern medical diagnostic imaging technology, stereotactic designed and product design ergonomic characteristics of human anatomy, application of modern technology and user-friendly design philosophy, can meet the requirements of neurosurgery, spinal surgery, ENT, plastic surgery, radiosurgery and 3D imaging. The operation table adopts the imported micro-motor, with electric gear mesh mechanism structure instead of the traditional hydraulic transmission. The table top adopts the carbon fibre material, with low x-ray absorption, reduces the radiation damage of x-ray to the patients and paramedic. 1150mm length for table top without metal (full carbon fibre), meet the requirement for 3D and 360° X-ray photograph of orthopedics navigation systems.

PARAMETERS:

Bed panel size(Carbon fibre):	(2400*530)±10mm	One key reset function
Max Height :	1000±20mm	Battery Optional
Min Height	700±20mm	Standard accessories
Route of longitudinal movement :	280±20mm(Electric)	Y-001 Mattress
Head tilt angle:	25°±2° (Electric)	Y-002 Controller
Foot tilt angle:	25°±2° (Electric)	Y-003 Lateral support
Left side leaning angle:	21°±2° (Electric)	Y-004 Carbon fibre arm support
Right side leaning angle:	21°±2° (Electric)	Y-005 Fixed belt
Back section up:	65°±2° (Electric)	Y-006 Anaesthesia screen
Back section down:	20°±2° (Electric)	Y-007 Leg support
Leg section up:	35°±5° (manual)	Y-008 Ccarbon fibre head frame
Leg section down:	90°±5° (manual)	Y-008 Detachable edge rail
Head section up:	45°±5° (manual)	P-001 Lock block/clamping block
Head section down:	90°±5° (manual)	
Bed body locking:	Electric	
Power supply (dual purpose) :	AC:220V 50Hz 300VA /DC:24V	
Net weight	300kg	

Bed panel X-ray attenuation coefficient(Optional): ≤0.44mmAL



DST-700

Luxury Electric Operation Table

Optional orthopedics shelf

Application of modern technology, combine with the idea of humanized design, meet the requirements of operation in the process of video observation and diagnosis, comprehensive strength, for each part of the human body surgery hospital operating room, i.e., the brain, heart, blood vessels, extracorporeal circulation, abdominal cavity, the gallbladder excision, renal resection, the neck, cervical vertebra, lumbar spine, spine, back, urology, gynecology, cranial nerve, surgical, rectum, plam traction, upper limb surgery, knee surgery, orthopedics lower limb traction (such as with orthopaedics traction frame) surgery.

PARAMETERS:

Length	2220 mm	Back section down	22°
Total width	600mm(with side rail)	Trendelenburg	24°
Head section length	300mm	Re-trendelenburg	24°
Back section length	490mm	Lateral tilt	22°
Main section length	540mm	Leg section adjustment	Up 15° Down 90°
Leg section length	580mm	Spread angle of leg plates	200°
Supporting bottom size	1030x650x200mm	Head section adjustment	Up 40° Down 90°
Castor size	125mm	N.W. Of overall table	280kgs
Height adjustment range	700-1000mm	Max loading weight	320kgs
Kidney bridge rise	≥110mm	Patient Horizontal max loading weight	300kgs
Sliding movement	350mm	Integrated battery	Yes, for over 150 movements
Back section up	70°	Voltage	220V 50Hz or 110V, 60HZ for choice
Standard accessories: Leg support-1 pairs, Arm supports-1pairs		Shoulder and waist holders- 2pairs, Anesthesia shelf-1pcs	

**DST-500A****Electric Operation Table**

*This type electric operation table has advantages of beautiful appearance, convenient, labor-saving operation, safe and reliable use, has the advantage of smooth action, etc.. The platform lifting, roll, tilt, before and after the lock is realized by the special motor, push-button. With the C type arm, no dead angle, parts can be disassembled, easy to use.

* product is easy to operate, strong versatility, wide application, material selection, workmanship excellent, quality reached the international level.

*Head section: pneumatic auxiliary device for head plate, plate easy to operate.

*Bed section: special composite material plate, X-ray high transmittance, good uniformity, meet the requirements of operation.

*Mattress: Memory foam, it adopts the international popular material, streamlined appearance.

*Base: thin base design, hidden castors, both to ensure enough ground contact, to ensure stable operation platform, but also can be the greatest degree of operation in the C-arm film.

PARAMETERS:

Length	2040mm	Leg board open	≥ 90°, detachable
Max Height	920mm	Back board fold upward	≥ 65°
Min Height	670mm	Back board fold downward	≥ 30°
Width	520mm	Left Tilt	≥ 19°
Reversed trendelenburg	≥ 20°	Right Tilt	≥ 19°
Trendelenburg	≥ 20°	Kidney bridge elevation	350mm
Lateral Tilt	≥ 19°	Power Supply	AC 220±22v, 50HZ±1HZ
Head board fold upward	≥ 20°, detachable	Translation moving	>300mm±20mm
Head board fold downward	≥ 90°, detachable	Battery	Optional
Up and down range	250mm	Standard accessories: Leg support-1 pairs, Arm supports-1pairs	
Leg board fold downward	≥ 90°, detachable	Shoulder and waist holders- 2pairs, Anesthesia shelf-1pcs	



DST-1A

Electric Operation Table

This type operating table is suitable for operation of any parts of human body. Besides all the functions of multi-purpose operating table, it has the unique advantages of its own. For example, the leg panel can be disassembled and easily adjusted which is excellent for urology operation; the table board can be translated by a great degree and the high strength transparent board is adopted on the table panel, which is very convenient for X ray examination; the back panel can be upforward folded, which can overcome the defect caused by translation. All the movement of the table is controlled by electric drive pusher with safety and reliability. This operating table can be matched up with the osteology support especially for C-arm.

PARAMETERS:

Table size	1950x550mm
Height adjustable range	660-910mm
Trendelenburg	$\geq 20^{\circ} \pm 2^{\circ}$
Reversed trendelenburg	$\geq 20^{\circ} \pm 2^{\circ}$
Lateral tilt	left $\geq 20^{\circ} \pm 2^{\circ}$ & right $\geq 20^{\circ} \pm 2^{\circ}$
Table surface Sliding range	350mm
Head Section adjustable range	up $\geq 45^{\circ}$ & down $\geq 90^{\circ}$
Safe bearing load	200kg
Back Section adjustable range	up $\geq 70^{\circ}$ & down $\geq 20^{\circ}$
Leg Section adjustable range	up $\geq 45^{\circ}$ & down $\geq 90^{\circ}$
Kidney bridge elevation	≥ 120 mm
Power supply	220V 50Hz
Standard accessories: Leg support-1 pairs, Arm supports-1pairs	Shoulder and waist holders- 2pairs, Anesthesia shelf-1pcs
Battery	Optional



DST-III

Electric Hydraulic Operation Table

*DST-III , the application of modern and humanism combined with the design concept, using electric imported hydraulic transmission technology, follow the development trend of the clinical operation, meet the image observation and diagnosis of the current operation process requirements.

*product is easy to operate, strong versatility, wide application , material selection, workmanship excellent, quality reached the international level.

*Fully meet the general surgery, heart, kidney, thyroid, extracorporeal circulation, head and neck, gallbladder and pleural, department of urology, department of orthopedics, department of gynaecology, rectum, spine surgical operation requirements.

Features:

- * The use of electric hydraulic transmission structure, stable speed, reliable performance;
- * The bed panel is a transparent material, the mattress adopts artificial wrapped high density memory foam;
- * The Lifting, translation, head and foot tilt operation driven by independent system.
- * Imported precision linear guide rail, ensure smooth running operation bed. 304 stainless steel material, the stainless steel cover, anti glare.
- * The design of large casters, mobile and flexible;
- * There are two sets of independent electric operating system, a remote control, a bed body standby operation control panel system, ensure the operation bed can still operate reliably in the remote controller failure.

PARAMETERS:

Size	2015 x 540 x 650-1000mm
Reversed /trendelenburg	25°/25° electric
V & ^ function	yes electric
Head section angle	up 30° down 90° manual
Leg section	up 28° down 90° manual
Tilt angle	left 22° right 22° electric
Back section	up 70° down 32° electric
Table panel x-ray attenuation coefficient ≤0.76mmAL	
Power Supply	AC 220±22v, 50HZ±1HZ
Translation moving	>350mm±20mm electric
Battery	optional
Standard accessories: Leg support-1 pairs, Arm supports-1pairs	
Shoulder and waist holders- 2pairs, Anesthesia shelf-1pcs	
One key reset function	optional



DST-III(B)

Electric Hydraulic Operation Table

- *.The controller device lock key protecting from misoperation.
- *.The table is made of penetrable X-ray material, there is cassette guide rail under the table which make radiography.
- *.The table is fixed a button which can control the electric wheel to draw in and put out so that the table can be pulled easily.
- *.Accessories and guide rail are made of stainless steel, polyurethane material etc.
- *.The main accessory hydraulic pump and electromagnetic valve are made of imported component from Taiwan.
- *.The mattress is made of form polyurethane, acid and alkali resistance, high temperature resistance and easy clean.
- *.The head board and leg board are adopted special imported gas-spring from Germany, it's convenient operate.

PARAMETERS:

Table size	2100x500mm
Height adjustable range	750-1000mm
Trendelenburg	≥ 22°
Reversed trendelenburg	≥ 22°
Lateral tilt	left ≥ 15° & right ≥ 15°
Table surface Sliding range	350mm
Head Section adjustable range	up ≥ 45° & down ≥ 90°
Safe bearing load	240kg
Back Section adjustable range	up ≥ 45°
Leg Section adjustable range	up ≥ 45° & down ≥ 90°
Kidney bridge elevation	≥ 120mm
Power supply	220V 50Hz
Standard accessories: Leg support-1 pairs, Arm supports-1pairs	Shoulder and waist holders- 2pairs, Anesthesia shelf-1pcs
Battery	Optional

DST-IB

Hydraulic Operation Table



HOT-H multi-purpose operation table is an improvement product developed on the basis of 3008 series head-manipulation type multipurpose operation table and is for use of general operations in chest surgery, department of ophthalmology, department of otolaryngology, department of obstetrics and gynecology, urological department etc. As a unit of non electronic control system comprehensive operation equipment, it is convenient and flexible to operate.

Features:

- * A multifunctional operation table that applies to various surgical operations and examinations.
- * Adoption of a foot hydraulic pump that integrates oil tank of oil pump and drive mechanism.
- * Adoption of the imported cylinder that controls actions of table top and is convenient to operate and reliable to locate.
- * Compact and rational structure and originative model.
- * Built-in waist bridge is convenient to use.
- * Both guide rail and accessories are made of corrosion resistant stainless steel.
- * All mattresses use the new-style materials that are acid and alkali resistant and disinfection durable.

PARAMETERS:

Size	1960 x 480 x 690~930mm
Reversed /trendelenburg	25°/25°
Lateral tilt	left≥18° & right≥18°
Head section angle	up 45° down 50°
Leg section	down 90°
Kidney bridge elevation	≥120mm
Back section	up 75° down 25°
G.W	295KGS
Power Supply	AC 220±22v,50HZ±1HZ
Translation moving	>350mm±20mm
Battery	optional
Standard accessories: Leg support-1 pairs, Arm supports-1pairs	
Shoulder and waist holders- 2pairs, Anesthesia shelf-1pcs	



DST-I

Hydraulic Operation Table

(head operating type)

1.The operating table is supplied for thoracic surgery, abdominal surgery, ophthalmology, obstetrics and gynecology, urinary surgery.

2.The rise and drop of operating table was handled by pedaling the oil pump, other various movements were handled by the mechanical transmission head.




Table size	2100 x 480 X 800 / 1045 mm
Trendelenburg & reversed trendelenburg	≥35° & ≥55°
Lateral tilt	left≥25° & right≥25°
Head Section adjustable range	up:90° & down:90°(p.s.:detachable)
Back Section adjustable range	up:85° & down:10°
Kidney bridge elevation	≥120mm
Leg Section adjustable range	down:90°(p.s.:detachable)

Standard accessories: Leg support-1 pairs, Arm supports-1pairs Shoulder and waist holders- 2pairs, Anesthesia shelf-1pcs

DST-II

Electric Ophthalmology Operation Table



PARAMETERS:

Size:	1950 x 520 mm
Height adjustable range	500 - 800 mm
Head Height adjustable range	550 mm
Power supply	AC220V , 50Hz

Standard accessories: Anaesthesia screen , Arm supports-1pairs , Remote controller



JR-9000

Interventional Imaging Table

JR-9000 Catheterization table,adopt the advanced electric,mechanic transmission technology,full carbon fibre table top,specially design for various movable medium,small C-arm to imaging of the Cardiovascular, peripheral-vascular,catheterization vascular ang Jography and DSA.The product conform the features of logical structure,easy control,graceful outlook and the function is advanced and achieved the quality of Jatarnalional.According to the clinical practice,it is confirmed as the best quality products in China.

Features:

- * The operation table adopts the imported micro-motor,with electric gear mesh mechanism structure instead of the traditional hydraulic transmission.Not need to change oil termly,and worry about the problem of leak oil.Product with fine steel can move equally and smoothly in any directions;also provides accurate positioning,trusty capability and durable lifetime.
- * Big gear feet wheel(the diameter of large truckle)especially use for medical treatment,the whole table is easily to be moved to clear the floor.
- * Test the table according to the CE European standard by Germany TUV company,achieve the European standard,obtain CE certificate.The operation table achieves the quality of imported products.
- * High precision imported beeline navigation rail,ensure the smooth and balanced movement of operation table.
- * Base with the ABS industrial plastic,helps prevent glare,impact and it is easily clear up and listerize.
- * Handle control is easily operating and takes control of the horizontal move of table top , the brake system is safe and sensitive.
- * All mattresses use the new-style materials that are acid and alkali resistant and disinfection durable.

PARAMETERS:

Table board material	cabon fibre
Whole length	2500±10mm
Length of radiolucent board	1750±10mm
Length without metal inspection	1200±5mm
Whole width	685±5mm
Width of radiolucent board	465mm
Thickness	40mm
X-ray absorption coefficient (aluminium equivalent)	≤0.44mmAL
Max.loading	240kg
Horizontal movement	manual
locking device	electromagnetic locking
longitudinal movement	manual
Move route	250±20mm
locking device	electromagnetic locking
Max.head protract	1880±20mm
Max.foot protract	1050±20mm
Vertical movement	Electric
Min.Bed height from the ground	750±20mm
Max.Bed height from the ground	1100±20mm
Move route	350±20mm
average velocity	8mm/s
Bed surface shifting/Horizontal movement/Locking controllor	
average velocity	8mm/s
Table body locking	Electric
Big Four direction wheel	free moving
Inner configurate	DC24V Panasonic lead-acid battery
Single-phase Alternating Current	220VAC 50/60Hz
Input maximum power	300VA
Protective tube	F5AL250V(5*20mm)
Standard accessories:	Memory foam mattress--1pcs
	Arm supports-1pairs



D0015

Orthopedics Traction Shelf



D001

Orthopedics Traction Shelf



D009

Spine Frame



D0010

Arm Traction Table



D0011

Neurosurgery Headframe



D0012

Neurosurgery Headframe

DST-3003

Luxury Electric Obstetric Table



We learn from Japan's advanced technology, the use of human engineering principles, curve design. The utility model is provided with a plurality of auxiliary functions for gynecology, obstetrics, gynecology, diagnosis and examination, including emergency abdomen, and the utility model is especially provided with a unique foot support

PARAMETERS:

Length	1900±20mm
Max Height	890±20mm
Width	680±20mm
Min Height	620±20mm
Up & down lifting range	250±20mm (Electric)
Back Section up	75°up(Electric)
Back Section down	5°down(Electric)
Supply	AC220±10%,50HZ
Seat section	up20° down 0° (electric)

Standard accessories arm supports, leg supports, foot controller



DST-IVA

Electric Obstetric Operation Table

PARAMETERS:

Size	1810 x 600 mm
Height	600-850mm
Reversed /Trendelenburg	≥20°
Back section up&down	≥75°
Auxiliary board	90° (level moving)
Supply	AC220±10%,50HZ
Battery	Optional



This bed is designed base on the concept of family accompanying delivery. It has got 5 national patents including an exterior patent. The driven and control system use imported motors which are safe, reliable, easy operated and with low noise. Both families and nurses can operate it.

- Bed head end and guardrail are made of high strength plastic, which have nice appearance. They can be assembled or detachable easily.
- Bed height is adjustable electrically; the back section upward, trendelenburg and reverse trendelenburg are electric controlled too. Leg section is removable. There is filth basin under the seat, it is easy to adjust the position of lying-in women and help delivery process by combining use of the leg holder, assistant table and foldable footrest.
- The bed is dirt free, easy cleaned. The central brake provides the bed a better stability.

DHC-II(6001)

Electric Delivery Bed

PARAMETERS:



Size	1920x920mm
Height	600-890mmmm
Back section	≥65°up
Reversed Trendelenburg /Trendelenburg	≥5°
Leg board	up 5° down 5°
Power supply	220V 50Hz
Safe loading	200kg

DST-3004B

Hydraulic Obstetric Table




PARAMETERS:

Size	1800 x 600 mm
Height	650-895mm
Max angle of back board	up>75° down<6°
Reversed trendelenburg	forward>10° backward>25°
Back board	600*800mm
Seat board	600*400mm
Leg board	540*500mm

PARAMETERS:

Size	1400 x 580 mm
Height	580-980mm
Back board up	65°
Reversed trendelenburg	82°-180°
Back board angle	45°



DST-V

Electric Obstetric Bed



DELIVERY BED



A44-4

Stainless Steel Obstetric Bed

Size : 1800 x 600 x 800 mm



A44-5

Expoxy Steel Obstetric Bed

Size : 1800 x 800 x 800 mm



A44-6

Stainless Steel Exam Bed

Size : 1800 x 600 x 800 mm



D4

Stainless Steel Exam Bed

Size : 1800 x 600 x 800 mm



D4-3

Stainless Steel Exam Bed

Size : 1950 x 650 x 680 mm



A44-1

Stainless Steel Delivery Bed

Size : 1800 x 600 x 800 mm



A44-3

Expoxy Steel Delivery Bed

Size : 1800 x 800 x 800 mm



A44-2

Stainless Steel Delivery Bed

Size : 1800 x 600 x 800 mm



A44

Stainless Steel Delivery Bed

Size : 1870 x 600 x 765 mm



DST-IV

Luxury Electric Obstetric Table

1. The main body is made of 45 # steel, The appearance parts are made of epoxy steel and 304 stainless steel, The base cover is high quality ABS;
2. Using imported motors with low noise, stable performance, Precise positioning, strong thrust,
3. The unique three - layer structure has achieved thermal insulation and cushioning. The minimum height 550mm platform, The platform is divided into backplane, seat plate and footboard. Through the operation handle, the upper and lower sides of the back plate can be changed, the bed will be tilted back and forth, so that the doctor can get an ideal surgical location and make the patient get comfortable seats.
4. The newborn baby's chassis can be rotated, which can correspond to the side care. The whole bed is running smoothly without noise and vibration.

PARAMETERS:

Size	1080*2058*550mm
Min height and Max height	550mm / 910mm
Translation distance of rotary tray support plate	60mm
Foot brace folding angle	0-90° gas spring
Foot brace rotation angle	Left 33° & right 33° (Electric)
Back Section up & down	60°up & 10°down (Electric)
Seat Section up & down	25°up & 25°down (Electric)
Power supply	AC220,2.5A 50HZ
Foot board folding angle	75° (manual)
Side handle rotation angle	180°(manual)
Tray support plate moving distance	150mm
Tray support plate rotation angle	90°



YD-580

Medical Ceiling Column




Bedhead Unit

We customize the bedhead unit according to your requirements

Nurse Call System

It is used in patient room & ICU room, we customize the nurse call system, you just tell us your hospital's bed details

The medical column is used for ICU, NICU room, the configuration is choiced by your requests.



YD-180E

Luxury ICU Ceiling Bridge


PARAMETERS:

Length	2000-3800mm
Payload gliders	100kg(220p)
Payload shelf	60kg (110p)
Ambient light	2*54w
Reading light	220v, 2*15w LED spotlights
RJ45	1pcs
Column length	700mm
Electrical socket(220V)	8pcs
Gas outlet(Germany)	O2x2,VACx1,AIRx1
Shelf x1,Drawer x1, basket x1, extension arm x1, infusion rack x1	

Optional item:

- 1. 1-4pcs shelves (25x10/30x10mm)
- 2. Drawer layer:1/2

Medical Gas Outlets Standards: optional British, German, French, Metric, Ohmeda, DISS, Chemetron, Purit-Bennett
Electrical Sockets:customized as clients' request



- ICU Bridge Type Ceiling-Mounted system, specially designed for intensive care rooms, the entire process wards, the recovery room. Supply flexible, first-rank nursing satisfy patient's separate treatment and nursing demand. Unique design lies in it all hides the electrical cable and pipeline in the ceiling mounted arm.
- The modular design and multiple configurations.
- As clients' requirements, providing reasonable ICU Combination Ceiling-Mounted pendant (Apart Dry-Wet).
- Advanced spray technology, the intelligent combination of innovative ways to make the bridge pendant with diverse functions.

YD-180E/X

ICU Ceiling Bridge (Dry And Wet Separated)



PARAMETERS:

Length	2000-3800mm
Payload gliders	100kg(220p)
Payload shelf	60kg (110p)
Ambient light	2*54w
Reading light	220v, 2*15w LED spotlights
RJ45	1pcs
Column length	700mm
Electrical socket(220V)	8pcs
Gas outlet(Germany)	O2x2,VACx1,AIRx1
Shelf x1,Drawer x1, basket x1, extension arm x1, infusion rack x1	

Optional item:

- 1. 1-4pcs shelves (25x10/30x10mm)
 - 2. Drawer layer:1/2
- Medical Gas Outlets Standards: optional British, German, French, Metric, Ohmeda, DISS, Chemetron, Purit-Bennett



YD-3600-2

Double Arm Ceiling Pendant

PARAMETERS:

Double arm/Single arm	800mm+800mm/800mm
Column length	1000mm
Max.loading capacity	350kg
Arm rotating angle	340°
Braking system	Optional Merchanical/Pneumatic braking
Electrical power	220V,50Hz
System	Electric
RJ45	1 pcs
Shelf	3 shelves
Electrical socket	6 pcs
Gas outlet(Germany standard)	O2x1,VACx1,AIRx1
Drawer x1,basket x1,extension arm x1,infusion rack x1,IV pole rack x1	
Optional item:	
1-4pcs shelves (25x10/30x10mm)	
Medical Gas Outlets Standards: optional	

British, German, French, Metric, Ohmeda, DISS, Chemetron, Purit-Bennett , Electrical Sockets:customized as clients' requests.
Optional anesthesia,endoscopy and surgery pendant

- Pendant Applicative Area: Surgical Room/ Emergency Room/ICU/Recovery room
- Ceiling Mounted installed.Within 340°revolved,Adjustterminal box location,convenience for medical staff operation.
- The suspension arm is imported by high-strength aluminum alloy extrusion, surface designed by the first level oxidized rust prevention.
- Friction brake + Pneumatic brake,adopted by German Standard Gas terminal(Insert and Pull more than 20,000 times.)
- Provide multiple configure,make you satisfy.
- Electrical up and down,precise positioning.



YD-3200

Single Arm Ceiling Pendant

YD-3600-2D

Double Arm Electric Ceiling Pendant

PARAMETERS:

Double arm/Single arm	800mm+800mm/800mm
Column length	1000mm
Max.loading capacity	350kg
Arm rotating angle	340°
Braking system	Optional Mechanical/Pneumatic braking
Electrical power	220V,50Hz
System	Electric
RJ45	1 pcs
Shelf	3 shelves
Electrical socket	6 pcs
Gas outlet(Germany standard)	O2x1,VACx1,AIRx1
Drawer x1,basket x1,extension arm x1,infusion rack x1,IV pole rack x1	

Optional item:

1-4pcs shelves (25x10/30x10mm)

Medical Gas Outlets Standards: optional

British, German, French, Metric, Ohmeda, DISS, Chemetron, Purit-Bennett , Electrical Sockets:customized as clients' requests.

Optional anesthesia,endoscopy and surgery pendant



Pendant Applicative Area: Surgical Room/ Emergency Room/ICU/Recovery Room

Ceiling Mounted installed.Within 340°revolved,Adjustterminal box location,convenience for medical staff operation.



The suspension arm is imported by high-strength aluminum alloy extrusion, surface designed by the first level oxidized rust prevention.

Friction brake + Pneumatic brake,adopted by German Standard Gas terminal(Insert and Pull more than 20,000 times.)

Provide multiple configure,make you satisfy.

Electrical up and down,precise positioning.

YD-3200D

Single Arm Electric Ceiling Pendant

PM-90



Multi-Parameter Patient Monitor

- 12.1 inches high resolution color TFT screen
- User friendly and configurable interface design
- Intuitive human voice alarms
- Powerful data review:
 - 480 hours graphic and tabular trends of all parameters
 - 1000 pieces of NIBP record storage
 - 700 pieces of alarm review
- Built-in rechargeable Li-ion battery for up to 5 hours of continuous work
- VGA output for external display(optional)
- Wire/wireless central monitoring system solution(optional)
- Standard Parameters:3/5-lead ECG,PR,RESP,NIBP,SpO₂,TEMP
- Optional: 2-IBP,EtCO₂(mainstream,sidestream),central monitoring system

● Respiration EtCO₂ (optional)



CAPNOSTAT5
mainstream module



LoFlo sidestream module



EtCO₂(Respiration)
"Plug and Play"



PM-15

Multi-Parameter Patient Monitor

- 15 inches high resolution color TFT display
- Optional full touch-screen for easy operation
- Attractive and durable silicon buttons with backlight
- ST segment analysis and 23 types of arrhythmia analysis
- Optional simultaneous 12-lead ECG analysis for further diagnostic confirmation
- Maximum 480-hour graphic and tabular trends,1000NIBP records, 700 alarm events and 2-hour ECG waveform storage and review
- ESU proof and defibrillation synchronization for different working environments
- Dual rechargeable Li-ion batteries for up to 10 hours of continuous work
- Dual USB ports enable software upgrade and also supports external USB storage
- Customer-defined color configurations for all parameters display setting (32 colors)
- Nurse call function and bi-directional communication solution
- VAG output for external display(optional)
- Flexible mounting and trolley solution(optional)

● Masimo SpO₂



Masimo M-LNCS™
SpO₂ extended cable



M-LNCS™ YI adult, pediatric,
neonatal reusable sensor



M-LNCS™ DCI adult
reusable finger clip sensor



M-LNCS™ DCIP pediatric
reusable finger clip sensor

- 1. 15-inch high resolution LED monitor
- 2. standard PC LED backlit silicon gel keyboard
- 3. one key regulation, easy operation
- 4. rechargeable built-in battery
- 5. abundant peripherals: DICOM3.0, VGA, video out, USB etc

C-60

Portable Color Doppler Ultrasound Machine

Detailed Product Description

- A. Imaging optimization technology;
- B. Compound enhance technology
- C. Imaging processing technology
- D. Imaging optimization technology
- E. Compound enhance technology
- F. Multi-beam parallel processing technology
- G. Tissue harmonic imaging

Displaying mode

- B, 2B, 4B, B/M, M
- B/C, B/C/D, B/D CFM, PRF
- PW, velocity, power (direction), histogram
- Triplex / Duplex

Probe

- Electronic convex, linear, transvaginal
- Wide band multi-frequency



- 1, 15-inch high resolution LED monitor
- 2, standard PC LED backlit silicon gel keyboard
- 3, one key regulation, easy operation
- 4, rechargeable built-in battery
- 5, abundant peripherals: DICOM3.0, VGA, video out, USB etc

C-60PLUS

High-class Portable Color Doppler Ultrasound Machine



Flexible and convenient operating system

A type of high array, high definition, multi-function full-digital laptop Color Doppler ultrasonic diagnostic apparatus. Light and portable, more clearer, delicate and smooth image, convenient to operate, and powerful electrical power ability supports the machine for consultation in various environment of cities, towns, outdoor. A variety of charging to ensure the examination consultation under different environments. Supernormal battery life 2-3 hours. 15 inch LED shows the images with more flexible ways.

Operating Mode:

- Black and white image : B, 2B, 4B, left B | M, B | D, PW, M, B partial enlarged mode
- Bloodstream image : Spectrum : B / D, B / C / D CFM: B | C | D, B | C | M, B | C pairs in real time , PW, CFM, CPA
- Color flow image adjustment parameters : a Doppler frequency , sampling frame position and size , baseline color gain, the deflection angle , the wall filter, the cumulative frequency.
- Display depth: Continuously adjustable



VET-6

Veterinary Ultrasound Scanner



1.Perfect Ergonomic Design

- Compact, sophisticated, portable
- 5.6 inch LED display
- Lightweight body
- Rechargeable built-in battery
- User-friendly operation interface
- Simple and beautiful operation panel
- Fast, convenient and flexible touch-type keyboard operation

2.Advanced Technology and Function

- Advanced digital signal processing technology
- Leading digital beam former technology
- NC dynamic frequency scanning technology
- Real-time point by point dynamic focusing reception technology
- Real-time dynamic aperture imaging

DESCRIPTIONS :

The C-80 color doppler ultrasound system bringing an entirely new level of diagnostic performance and workflow efficiency. The C-80 system, combines best-in-class image quality and a robust set of features to meet daily clinical needs. It enhances diagnostic confidence with high quality color and power Doppler, 2D-mode imaging. All-digital signal processing and multibeam formation technology provides best-in-class image quality in 2D-mode, Tissue Harmonic Imaging, Dynamic Tissue Optimization Technology and Doppler modes.

Data and information communication can be easily carried out. It also can connect to PACS(Picture Archiving and Communication System)as a terminal.

PACS is a medical picture archiving and communication systems, which greatly facilitates network management system of the hospital. It is convenient for remote diagnosis. The machine is powerful, convenient and easy to operate, supports B, B/B, B/4B, B/M, B/PWD, CFM, CDE, B/CFM/D. Meanwhile, the display can be moved from up and down, forward and backward, turn to left and right. These powerful features have greatly boosted diagnostic confidence.

USER INTERFACE

- Intuitive Windows-based operating principles
 - User-centric control panel with HomeBase layout and control customization
 - On/Off task light and back-lit illumination of control panel
 - Variable brightness indicates active state of function keys
 - Easily accessible, full size QWERTY keyboard for text entry, function keys and system programming
- Cine Playback, Multiple Arrows, Configurable Worksheets, Exam Review, Pictograms (Body Marks), System Setup Menu

C-80

Full Digital 3D Mobile Ultrasound Machine Color Doppler





Mattress holders



Leg plate at both ends of the regulation, adjust the comfort of the leg board



Extend bedboard



TENTE silent castors central control braking



Manual CPR



X-ray cassette



anticollision wheels with IV pole



DHC-III (FM04)

Luxury Multi-functional Electric ICU Bed With Weighing Function

PARAMETERS:

Size: L2170xW1080xH500/850mm

Bed plate: L1950xW910mm

The whole bed can weigh, with error controlled in 1-3%

Back plate adjusting $70^{\circ} \pm 5^{\circ}$, leg plate adjusting $45^{\circ} \pm 5^{\circ}$

Head and foot tilt $16^{\circ} \pm 3^{\circ}$

Optional weighing function



Model: FM03
4 functions



Model: FM04
5 functions



DHC-IV(FP07)

Luxury Electric ICU Bed With 7 Functions



PARAMETERS:

Size	2240 x 1110 x 630/1030 mm	Overall lifting	630 -1030mm
Back adjustment	40°±2°	Lateral tilt	12°±2°
Leg adjustment	13°±2°	Luxury central control castors	
Trendelengurg angle	17°±2°	ReverseTrendelengurg angle	25°±2°

4 pieces of European syle guaedrail, upside down structure

Foot controller, nurse controller, patient controller

Optional mattress, bedside cabinet, overbed food table

Five motors



Weighing operator



Foot control pedal



Nurse control panel



Panel on guardrail



Bedboard soft connection



ABS guardrail



TENTE Castors



Bed panel with punching design

DHC-IV(FP05)

Electric ICU Hospital Bed With 5 Functions



PARAMETERS:

Size	2300 x 1050 x 530mm/710mm	Bedboard	1955 x 900 mm
Net weight	158kg	Loading capacity	250kg
Back board	0-65°	Foot board	0-43°
Trendelenburg/Reversed Trendelenburg			0-10°
"Reset" functions , be can be back to "zero" position quickly. (manual CPR)			
Backup battery for emergency power failure			
Bed board has safe connection which make the bending of the bed board.			
Detachable head&foot board in PP plastic,can be fixed when assembled.			
4 pcs foldable PP siderails			
Central locking castors controlled by electric.			
I.V.pole mount bracked at 4 corners of the bed frame			
With easily operated touching switch,solve the problem of limited service life,mechanical keypad and keypad invalid due to constantly operating			
125mm silent castors			
Optional mattress, bedside cabinet,overbed food table, i.v stand pole etc.			



FO03 Electric Hospital Bed

PARAMETERS:

Bed size 2280 x 1120 x 420/780mm

Height lifting ≥ 200 mm

Backboard up $\geq 65^\circ$

Leg board up $\geq 35^\circ$

Asymmetric European style guardrail ,Turn-over structure

Luxury central control castors,brake the bed with one stamp.

Bed panel is made of one-step molding cold rolled steel, with excellent performance in anti-slipping and air permeability.



FO - 2 Manual Hospital Bed

PARAMETERS:

Bed size 2280 x 1120 x 420/780mm

Height lifting ≥ 200 mm

Backboard up $\geq 75^\circ$

Leg board up $\geq 37^\circ$

Asymmetric European style guardrail ,Turn-over structure

Luxury central control castors,brake the bed with one stamp.

Bed panel is made of one-step molding cold rolled steel, with excellent performance in anti-slipping and air permeability.



FO - 1 Manual Hospital Bed

PARAMETERS:

Bed size 2280 x 1120 x 500mm

Backboard up $\geq 75^\circ$

Leg board up $\geq 37^\circ$

Asymmetric European style guardrail ,Turn-over structure

Luxury central control castors,brake the bed with one stamp.

Bed panel is made of one-step molding cold rolled steel, with excellent performance in anti-slipping and air permeability.



FP04 Electric Hospital Bed



Bed size	2280 x 1020 x 420/780mm
Height lifting	≥200mm
Backboard up	≥65°
Leg board up	≥35°
Luxury central control castors,brake the bed with one stamp	

FP03 Electric Hospital Bed



Bed size	2280 x 1020 x 420/780mm
Height lifting	≥200mm
Backboard up	≥65°
Leg board up	≥35°
Luxury central control castors,brake the bed with one stamp	

FP - 2 Manual Hospital Bed



Bed size	2280 x 1020 x 420/780mm
Height lifting	≥200mm
Backboard up	≥75°
Leg board up	≥37°
Luxury central control castors,brake the bed with one stamp	

FP - 3 Manual Hospital Bed



Bed size	2280 x 1020 x 500mm
Backboard up	≥75°
Leg board up	≥37°
Luxury central control castors,brake the bed with one stamp	
Asymmetric European style guardrail ,Turn-over structure	

DHC-II(FE01)

Electric Medical Bed With 5 Functions



PARAMETERS:

Size	2175 x 1090 x 420/720mm
Back board tilt	65° electric / manual
Leg board tilt	35° electric / manual
Central controlling	Silent castors
Optional mattress, bedside cabinet, overbed food table , IV Stand	Optional 5,3,2 functions or manual function

DHC-II(FE01)

Electric Medical Bed With 3 Functions



FN-2

Manual Hospital Bed With 3 Functions



PARAMETERS:

Size	2200 x 1020 x 420/760mm or height 500mm for two function
Back board tilt	75°
Leg board	37°
Safe loading	250kg
Optional mattress, bedside cabinet, overbed food table , IV Stand	

FN-1

Manual Hospital Bed With 2 Functions



DHC-II (FH05)

Economical Electric Home Care Bed



DHC-II (FH06)

Luxury Electric Home Care Bed



PARAMETERS:



Size	2150x1030x460mm	Height(from floor)	460-860mm
Back plate adjusting	0~70°	Knee section adjusting	0~40°
Trendelenburg/Reversal trendelenburg		13°/13°	
Auto contour position		CPR functions	
Safety load capacity		250kgs	
Optional mattress, bedside cabinet, traction shelf			

DHC-II (FH07)

Home Nursing Bed



DHC-II (FH08)

Home Nursing Bed



Bed size : 2100 x 960 x 540mm (including mattress)
 The hole size: 240*200mm
 The overbed table : 1030 x 300mm
 The castors : 125mm braking

Bed size : 2100 x 960 x 540mm (including mattress)
 The hole size: 240*200mm
 The overbed table : 1030 x 300mm
 The castors : 125mm braking

FA-5

Electric Bed With 3 Functions



Size : 2200 x 960 x 420/720mm

FA-6

Manual Hospital Bed With 3 Functions



Size : 2200 x 960 x 440/790mm

FA-7

Double Cranks Hospital Bed



Size : 2200 x 960 x 520mm

FA-8

Double Cranks Hospital Bed



Size : 2150 x 950 x 500 mm

A25

Double Cranks Hospital Bed



Size : 2150 x 900 x 500mm

A23

Single Crank Hospital Bed



Size : 2150 x 900 x 500 mm

A24

Hospital Bed



Size : 2150 x 900 x 500mm

A26

Single Crank Hospital Bed



Size : 2150 x 900 x 500 mm

A27

Hospital Bed



Size : 2150 x 900 x 500mm

A31

Double Crank Hospital Bed
Stainless steel



Size : 2040 x 900 x 500 mm

A34

Double Cranks Hospital Bed
Steel coating



Size : 2040 x 900 x 500mm

A33

Hospital Bed Stainless steel



Size : 2040 x 900 x 500 mm

A36

Hospital Bed



Size : 2040 x 900 x 500 mm

A37-1

Orthopedic Traction Bed



Size : 2020 x 900 x 500 mm

A39D

Child Bed With 1 Function



Size : 1600 x 850 x 610 mm

A39

Double Cranks Child Bed



Size : 1770 x 900 x 550 mm

A39-1P

Luxury Child Bed With 2 Functions



Size : 1770 x 900 x 550 mm

CR



Size : 200/110cm

B-2H

Luxury Electric Transfer Stretcher



With unique motion design, this cart could transport patient comfortably between hospital bed and transfer cart steadily without need of moving patient

Surface size: 2020mm×680(mm)

Size when connected: 3900mm×680(mm)

lifting range : 590mm~880(mm)

move range of shifting plate : 510mm

move speed of shifting plate : 30mm/s

Max loading: : 135kg

Accessories :

IV pole 1piece hand controller 1set

power cord 1piece connecting cord 1piece

B-2

Medical Connecting Transfer Stretcher



Product composition:

This product has A car, B and C three major components, including car A and car B, each has A screw and four control of castor, senior C car there are two big barrier and ABS engineering plastics of bed plate.

PE docking cart A

Size: 1790x450 mm , Height: 570~850mm

Movements distance: 280±20 mm(manual)

PE docking cart B

Size: 1835x450mm , Height: 570~850mm

Movements distance: 280±20 mm(manual)

PE docking cart C

Size: 1955x585mm

Width(guardrail up): 780±20 mm , Width(guardrail down): 710±20 mm

Height : 220±20 mm

Back adjustment: 90°±5°(manual)

PE docking cart (ABC) Size: 3575x770mm

B-1Y

Hydraulic Emergency Stretcher



PARAMETERS:

This Hydraulic Stretcher has Central-lock castors,diameter 150mm;bed board for X-ray Examination;

Hi-low range:800mm-550mm;

Backrest angle:0-60°

Knee angle:0-40°

Trendelenburg & Rev.Trendelenburg:0-12°

B-1

Manual Emergency Stretcher



PARAMETERS:

Size:L1930*W750*H600/830mm

Specification:

- 1.main frame of the flat cart is jet-molding carbon steel,solid and good looking.
- 2.cart surface and guardrail is structured by one-step molding ABS material,solid and elegant.
- 3.adopting separate design: back plate is adjusted by air spring,could be operated by only one hand,convenient to use.
- 4.the crank of the cart(the screw is equipped with a two-way protection device) can adjust height of the cart surface from 600 to 830mm
- 5.equipped with international advanced central-control braking system,stable and reliable,the oriental castors are easy to operate.
- 6.meeting requirement from overall design concept,more humanized,ideal for transporting patients in hospital

B-1A

Manual Transfer Stretcher



PARAMETERS:

Hospital transfer stretcher trolley

Features:

1. It is a type of manual transfer stretcher. Turn by deasil, the stretcher will be up, and turn by widdershins, the stretcher will be down;
2. Aluminum-Alloy protective rails, can be folded and hidden under the surface on no use;
3. Central-lock castors, diameter 150mm. Matched with direction castor, help the operator controlling the direction easily;
4. Rotate the crand and handle on bedboard, patient on the mattress can be moved parallelly to another bed or flat table;
5. Hi-Low range: 920mm-610mm
6. Weight: 89 kgs

B4-1

Emergency Trolley



It is made of high-strength aluminum alloy thick-wall pipes. The bed surface is made of foamed cushion and the adjustable back making the wounded fell comfortable. This machine is very convenient. One person is enough to send the patient (with the stretcher) to the ambulance car. The mechanical foldaway legs are controlled by the handles on both ends of the stretcher. And the legs can be folded under the stretcher automatically when loading into the vehicle. Also the stretcher can be locked by the fixing devices when on the ambulance car. This type of stretcher is characterized by its light-weighted, durable, safe usage and easy for sterilization. The stretcher adopts stretcher trolley and stretcher to separate structure.

high position: 1900*550*920mm
 low position: 1900*550*260mm
 maximum angle of the bach: 75°
 packing size (1pc): 1980*620*280mm
 loading bearing: ≤159kg
 N.W: 43kg
 G.W: 49kg



MT MEDICAL

MEDICAL TROLLEY SERIES

B4

Stainless Steel Emergency Trolley



Size : 1950 x 600 x 750 mm

B4-2

Portable Emergency Stretcher



Size : 2080 x 550 x 130 mm

MWT-1
Mobile Working Station



Size : 730 x 470 x 910 mm

MWT-2
Mobile Working Trolley



Size : 1050 x 635 x 1110 mm

CQ-08
Emergency Trolley



Size : 810 x 520 x 960 mm

CQ-09
Emergency Trolley



Size : 615 x 550 x 1010 mm

CQ-06
Emergency Trolley



Size : 762 x 520 x 1045 mm

CQ-10
Emergency Trolley



Size : 630 x 470 x 910 mm

CM-01
Anesthesia Vehicle



Size : 615 x 550 x 1010 mm

CM-02
Anesthesia Trolley



Size : 630 x 470 x 1500 mm

CM-03
Anesthesia Trolley



Size : 810 x 520 x 960 mm

CM-07
Anesthesia Vehicle



Size : 762 x 520 x 1045 mm

CZ-08
Treatment Trolley



Size : 775 x 635 x 1080 mm

CZ-04
Treatment Trolley



Size : 785 x 480 x 995 mm

CZ-05
Treatment Trolley



Size : 762 x 520 x 955 mm

CF-05
Medicine Trolley



Size : 775 x 635 x 1080 mm

CF-02
Medicine Trolley



Size : 785 x 480 x 995 mm



MT MEDICAL

MEDICAL TROLLEY SERIES

CF-06

Medicine Trolley



Size : 630 x 490 x 1070 mm

CF-05

Medicine Trolley



Size : 750 x 450 x 1220/1275 mm

B26

Medicine Trolley Stainless Steel



Size : 570 x 400 x 785 mm

B19

Stainless Steel Nursing Trolley



Size : 700 x 420 x 930 mm

CS-05

Infusion Trolley



Size : 630 x 470 x 910 mm

CS-04

Infusion Trolley



Size : 630 x 470 x 910 mm

CB-02 (25)

Medical Files Cart



Size : 510 x 390 x 780 mm

CB-07 (50)

Medical Files Cart



Size : 820 x 390 x 880 mm

CH-02

Luxury Nursing Cart



Size : 1185 x 570 x 950 mm

CH-07

Nursing Cart



Size : 800 x 590 x 890 mm

CH-06

Nursing Trolley



Size : 570 x 430 x 890 mm

B34

Stainless Steel Nursing Trolley



Size : 960 x 550 x 920 mm

B32

Stainless Steel Nursing Trolley



Size : 800 x 600 x 920 mm

B7-2

Luxury Baby Cot



Weighing, Music play and ambient temperature display

B10

Stainless Steel Baby Trolley



Size : 840 x 530 x 840/960 mm

B9

Stainless Steel Baby Cot



Size : 920 x 490 x 1020 mm

B15

Instrument Trolley
Stainless Steel



Size : 750 x 450 x 910 mm

B18

Instrument Trolley
Stainless Steel



Size : 750 x 450 x 910 mm



MT MEDICAL

MEDICAL TROLLEY SERIES

B13

Instrument Trolley



Size : 750 x 450 x 850 mm

B14

Instrument Trolley Stainless Steel



Size : 1400 x 450 x 910 mm

B16

Treatment Trolley Stainless Steel



Size : 750 x 450 x 910 mm

CZ- 05

Nursing Trolley



Size : 610 x 410 x 830 mm

CZ- 06

Luxury Instrument Trolley



Size : 620 x 420 x 860/930mm

CY- 04

Instrument Trolley



Size : 630 x 470 x 910 mm

CY- 01

Instrument Cart



Size : 730x 475 x 1415 mm

B30

Stainless Steel Sending Goods Cart



Size : 1100 x 520 x 930 mm

B35

Stainless Steel Mayo Table



Size : 590 x 450 x 940/1200 mm

MEDICAL CABINET SERIES



B36

Stainless Steel Trolley



Size : 690 x 450 x 940/1200 mm

B40

Power Supply Instrument Vehicle



Size : 650 x 550 x 890 mm

B70

Instrument Trolley



Size : 945 x 450 x 810 mm

B37

Oxygen Bottle Trolley



Size : 810 x 530 x 810 mm

C12

Stainless Steel Working Stable



Size : 1800 x 700 x 810/1750 mm

C32

Medicine Cabinet With Stainless Steel Foundation And Table Surface



Size : 900 x 240/400 x 1750 mm

C35

Medicine Cabinet



Size : 900 x 240/530/400 x 1750 mm

C36

Stainless Steel Instrument Cabinet



Size : 900 x 400 x 1750 mm

C38-1

Stainless Steel Instrument Cabinet



Size : 900 x 400 x 1750 mm



MT MEDICAL

MEDICAL CABINET SERIES

C39

Instrument Cabinet



Size : 900 x 400 x 1750 mm

C47

Steel Disassembly Type Six Door Locker



Size : 960 x 400 x 1750 mm

C42

Steel Disassembly Type Four Door Locker



Size : 900 x 400 x 1750 mm

C49

Steel Disassembly Type Three Door Locker



Size : 960 x 400 x 1750 mm

C48

Steel Detachable Two Door Locker



Size : 900 x 400 x 1750 mm

C50

Shoes Cabinet



Size : 960 x 400 x 1750 mm

C58

Stainless Steel Cabinet



Size : 900 x 250/400 x 1750 mm

C37

Steel Dismounting Type Instrument Cabinet



Size : 900 x 400 x 1750 mm

C38 CR

Stainless Steel Dismounting Type Instrument Cabinet



Size : 900 x 400 x 1750 mm

MEDICAL CABINET SERIES



MT MEDICAL

C38

Steel Dismounting Type
Instrument Cabinet



Size : 900 x 400 x 1750 mm

C40

Steel Dismounting Type
Instrument Cabinet



Size : 900 x 400 x 1750 mm

C42

Sterile Cabinet



Size : 900 x 400 x 1750 mm

C53-1

Stainless Steel Embedded
Instrument Cabinet



Size : 810 x 530 x 810 mm

C7

ABS Bedside Cabinet



Size : 480 x 480 x 750 mm

C8-1

ABS Bedside Cabinet



Size : 480 x 480 x 880 mm

CC-04

Luxury Bedside Cabinet



Size : 520 x 400 x 725 mm

CC-03

Bedside Cabinet



Size : 480 x 470 x 795 mm

CC-05

Luxury Bedside Cabinet



Size : 520 x 470 x 700 mm



MT MEDICAL

MEDICAL CHAIR SERIES

C13 - 1

Stainless Steel Basin
Sink With Sensor



Size : 1600 x 620 x 1900 mm

D11

Stainless Steel Accompany
Chair



Size : 700 x 650 x 880 mm

D10

Steel Coating Accompany
Chair



Size : 700 x 650 x 880 mm

D12 - 1

Luxury Accompany Chair



Size : 780 x 650 x 1060 mm

D5

Stainless Steel Infusion Chair



Size : 700 x 1120 x 1100 mm

D6

Steel Coating Infusion
Chair



Size : 700 x 1120 x 1100 mm

D6 - 1

Transfusion Chair



Size : 700 x 850 x 1100 mm

D6 - 2

Luxury Transfusion Chair



Size : 770 x 680 x 950 mm

D6 - 3

Luxury Transfusion Chair



Size : 710 x 1130 x 1100 mm

X-RAY VIEWING BOX SERIES



D1

New X-ray Film Viewer



Size : 440 × 585 × 110mm

D2

New X-ray Film Viewer



Size : 900 × 585 × 110mm

D3

New X-ray Film Viewer



Size : 1210 × 585 × 110mm

LED - T1

LED X-ray Film View Box



Size : 500 × 500 × 48mm

LED - T2

LED X-ray Film View Box



Size : 1230 × 500 × 48mm

D19

Doctors Stool



Size : Dia340 × 750/920 mm

D25

Doctors Mobile Stool



Size : Dia320 × 460/620 mm

D26

Medical Stool



Size: Dia320 × 430/590 mm

D13

Stainless Steel Folded Screen



Size : 1910 × 1800 mm



MT MEDICAL

MEDICAL SCALE SERIES

D21

Waiting Chair



Size : 1800 x 680 x 800 mm

D22 (A)

Steel Coating Footstep



Size : 620 x 380 x 420 mm

D22 (B)

Stainless Steel Footstep



Size : 540 x 400 x 400 mm

S - 2

Adult Scale



160KG

S - 1

Adult Scale



160KG

S - 4

Baby Scale



20KG

S - 3

Baby Scale



20KG

C55-1

Luxury Overbed Food Table



Size : 850 x 450 x 820/1030 mm

C55

Luxury Overbed Table



Size : 555 x 450 x 940/1200 mm

MEDICAL MATTRESS SERIES



C58

Folding Overbed Table



C57

Overbed Table



D18

Mattress With PU Leather
And High Density Foam



Size : 1950 x 800 x 100 mm

D15

Medical Air Mattress



Size : 2000 x 900 x 70mm

D17

Mattress With Full Foam



Size: 1950 x 800 x 80(mm)

D16

Mattress With Half Of Palm And
Foam



Size : 1950 x 800 x 80(mm)

D8

Expoxy Steel Infusion Chair



Size : 1000 x 600 x 400 mm

D7

Blood Donar Chair

Optional Stainless Steel
Or Steel Coating



Size : 1820 x 520 x 560 mm

B71

Clinical Examination Bed

Optional Stainless Steel
Or Steel Coating



Size : Dia360 x 2 x 760 /1260mm




优业家具

YOUYE 有限公司
FURNITURE CO.,LTD





Company introduction
公司简介




(微信扫一扫)联系我们,我们将竭诚为您服务!



联系人：游小姐



联系人：郭先生



联系人：游先生

广东优业实业有限公司简介

广东优业实业有限公司是一家集设计、开发、生产、销售为一体的、专业性钢制办公家具制造厂商。

位于广东首府广州市白云区，距离广州机场，火车站，省汽车站，市汽车站，黄埔码头不到40分钟的路程，交通极其便利。

优业家具成立于2005年，历经多年，现已发展成为颇具规模办公家具制造商，主要产品包括：文件柜、资料柜、卷门柜、储物柜、更衣柜、电脑台、办公台、密集柜、学生公寓床，课桌椅系列等多种金属结构产品，工厂占地面积20亩（13000平方），其中车间近8000平方，大小机床200余套，采用先进的静电喷涂工艺。3000平方成品仓库，可以提供各常规款式2000余套，足以及时供应国内外代理商的需求。

优业一贯秉承“稳健踏实，勇于创新，追求卓越”的经营理念不断推出满足现代办公需求的新产品，以专业的设计，精湛的工艺为全球的客户带来双赢的共享利益。

WE
HAVE THE
HONOR OF

优业荣誉



追求卓越 务实树典范

回顾过去,创业维艰。

公司坚持以创新为灵魂,以质量为生命
以服务为天职,建立起优质专业的品牌
形象,获得多项殊荣,赢得了社会各界
的一致赞誉。


展望未来,任重道远。

荣誉是肯定,更是鞭策,优业将一如既
往的锐意进取以更优的品质回馈用户。



检测报告

HONOUR QUALIFICATION PAPERS



企业理念

终身学习团队文化

树立终身学习的信念，带领全体员工打造成学习型团队，为企业发展提供智力支持。

和谐创业的人本文化

坚持企业发展的主体是员工，企业发展的动力靠员工，企业发展成果与员工共享。

永无止境的创新文化

不断更新观念、拓宽思路、创新技术、用人性化多功能的产品引领金属家具市场新潮流。

争比奉献的感恩文化

感恩是一种生活态度，做人的起码品格。感恩也是一种责任，用感恩之心构建和谐社会。

企业观点

人才观：企业即人，企业为人，企业靠人。

实施“用事业凝聚人才，用实践造就人才，用机制激励人才，用制度保障人才”的人才战略。

人才是企业最宝贵的财富。企业即人，在优业这个大家庭里，每个人的聪明才智和辛勤劳动都将获得全体成员的一致认可；企业为人，优业不仅重视人才，更注重为人创造良好的工作环境，使员工都能用其所长；企业靠人，优业的每一份子紧密团结，团结一心，努力将自己的知识技能凝聚成为企业的集体智慧。

服务观：您的满意是我们最大的追求

企业因客户而存在。只有实现了客户的利益，企业才能求得自身的发展。优业始终将客户利益摆在第一位，努力缔造高品位办公文化内涵

真诚对待客户，让客户百分百满意。优业的品质与优质的服务是优业永远的承诺。

市场观：没有饱和的市场，只有饱和的产品。

谁了解市场，洞察先机，果断决策，谁就能赢得市场。优业将始终不渝地坚持以市场为导向，精益求精，以不断创新的“零缺陷”产品质量服务引领中国“金属行业”，铸就金属行业“长城”。

企业精神

团队精神

优业强调的并不是表面的和谐，而是深层的凝聚。任何时候优业都强调团队的力量，依靠团队的力量，提升团队的力量。有人说，优业是一家“操正步走”的公司，也正是由于公司上下齐心协力，步调一致才走到今日的成就。

敬业精神

敬业是一种付出，代表员工对企业的忠诚；敬业也是一种收货，表示企业给员工理想的发挥空间。优业相信，每位员工身上都蕴藏着无比的潜能。为了实现人生的理想和信念，每位员工都严格要求自己，兢兢业业，不断为公司奉献着自己的聪明才智。

资质证书

QUALIFICATION CERTIFICATE



企业证书

QUALIFICATION CERTIFICATE



GLASS DOOR FILING CABINET

玻璃文件柜系列



追求品质永远是我们的信仰

每一个精品都需要时间的雕刻，每一个经得起鉴定的艺术品，都必须承受历史的打磨，每一个优业用户，均乐享优业的创造，优业竭诚为用户提供优良的服务。

采用优质冷轧钢板模压成型，表面经除油、去锈、磷化、酸洗等十道工序前处理，国际最新流行色亚光静电自动喷粉，高温固化而成，防锈蚀性能显著。

YOUYE Group file cabinet products, a large fo a shape, phosphate anti-corrosion treatment, automatic spray technology, reasonable design, fashionable, high quality, great to meet the needs of users.

1, available on demand free combination; 2, new hinges, door open light; simple number registered, convenient access to data; quality locks open and flexible, each opening rate is extremely low; 5, blue point deduction hand, novel personality, ABS is not easily broken; double slide, drawer sports facilities, and pull out the 90%; 7, the partition can be freely adjusted up or down or decrease up or down or decrease the full use of cabinet space; 8, skim, phosphating twelve stations before processing to ensure that the coating is not easy to fall off, and weather resistance.





TG-008-1 掩门全玻璃文件柜
Full Glass Open Door Cabinet
900*400*1800mm



TG-023 玻璃移门地柜
Short Glass Cabinet
900*400*900mm



TG-013 单门六抽半玻璃文件柜
Six Drawers Half Glass Cabinet
900*400*1800mm



TG-023-1 玻璃移门文件柜
Short Glass Cabinet
1180*400*900mm



TG-008 移门全玻文件柜
Full Glass Sliding Door Cabinet
900*400*1800mm



TG-009-1 上玻璃下掩门文件柜
Middle Drawer Glass Filing
Cabinet
900*400*1800mm



TG-063-1 上玻下铁移门文件柜
Half Glass Cabinet
900*400*1800mm



TG-063 上铁下铁掩门文件柜
Half Glass Cabinet
900*400*1800mm



TG-066 上移门下掩门中两抽文件柜
Middle Drawer Glass Filing Cabinet
900*400*1800mm



TG-090 上移门下掩门中六抽文件柜
Middle Drawer Glass Filing Cabinet
1180*400*1800mm

STEEL MADE
FILING CABINET

钢制文件柜系列



中国家具协会金属专业委员会会员单位

优业家具文件柜系列产品，采用大面积一次成型，磷化防腐蚀处理，喷塑工艺等先进技术，设计合理，款式新颖，质量上乘，极大满足了用户的需求。

产品特点：

- 1、可按需自由组合；
- 2、新型合页，柜门开启轻便；
- 3、简易号册，存取资料方便快捷；
- 4、优质锁具，开启灵活，互开率极低；
- 5、蓝点扣手，新颖个性，ABS不易破碎；
- 6、复式滑轨，抽屉运动便利，并能拉出90%；
- 7、隔板可自由下调整或增减充分利用柜体空间；
- 8、脱脂、磷化等处理，确保涂膜不易脱落，且耐候性强。

YOUYE Group file cabinet products, a large fo a shape, phosphate anti-corrosion treatment, automatic spray technology, reasonable design, fashionable, high quality, great to meet the needs of users.

1, available on demand free combination; 2, new hinges, door open light; simple number registered, convenient access to data; quality locks open and flexible, each opening rate is extremely low; 5, blue point deduction hand, novel personality, ABS is not easily broken; double slide, drawer sports facilities, and pull out the 90%; 7, the partition can be freely adjusted up or down or decrease up or down or decrease the full use of cabinet space; 8, skim, phosphating twelve stations before processing to ensure that the coating is not easy to fall off, and weather resistance.





TG-065 四门无抽文件柜
Four Doors Cabinet
900*400*1800mm



TG-011 双门内两抽文件柜
Double Doors Inside
Two Drawers Cabinet
900*400*1800mm



TG-062 双门文件柜
Double Doors Cabinet
900*400*1800mm



TG-038 四层档案柜
Four Floors Archive Cabinet
900*400*1800mm



TG-039 五层档案柜
Five Floors Archive Cabinet
900*400*1800mm



TG-065-1 双门地柜
Double Doors Short Cabinet
900*400*1800mm

ROLLING
DOOR CABINET

卷帘门文件柜



FASION DESIGN

AND EXCELLENT QUALITY

Bring You Comfortable and Gladdergladdest Mood

潮流的设计，精湛的工艺，
带给你舒适愉悦的工作心情





卷帘门文件柜-01
01 ROLLING DOOR CABINET
1000*420*1980mm



卷帘门文件柜-02
02 ROLLING DOOR CABINET
1000*420*1980mm



卷帘门文件柜-03
03 ROLLING DOOR CABINET
1000*420*1020mm



卷帘门文件柜-04
04 ROLLING DOOR CABINET
1000*420*1500mm



卷帘门组合地柜
05 ROLLING DOOR CABINET
1000*500*650mm



卷帘门文件柜-06
06 ROLLING DOOR CABINET
1000*420*1980mm

COMBINATION
LOCK

保密文件柜系列



采用国保密码锁具，国家检测认证，使用更放心

GB-01 四门中两抽密码锁柜
900*400*1800mm

人性化的设计不是一句空话，每个解决方案，我们看到的不是生硬的，材料，而是未来使用的主人。

Humanized desing is not a slogan Every design comes from experience of its future user.




锁具型号：M103D



锁具型号：D153



GB-01 上下掩门中两抽密码文件柜
W900*D400*H1800



GB-02 两门无抽密码文件柜
W900*D400*H1800




GB-04 双门地柜密码锁柜
W900*D400*H900





GB-05 上下掩门中两抽密码文件柜
W900*D400*H1800



GB-06 两门无抽密码文件柜
W900*D400*H1800



GB-07 双门地柜密码锁柜
W900*D400*H900



MOBILE
CABINET

快捞资料柜



FASION DESIGN

AND EXCELLENT QUALITY

Bring You Comfortable and Gladderladdest Mood

潮流的设计，精湛的工艺，
带给你舒适愉悦的工作心情





F4快捞柜-01
01 QUICKLY ARK
900*450*(1310-1010-700)mm



F4快捞柜-02
02 QUICKLY ARK
450*620*(1310-1010-700)mm



F4快捞柜-03
01 QUICKLY ARK
450*620*(1310-1010-700)mm



F4快捞柜-04
04 QUICKLY ARK
450*620*(1310-1010-700)mm



A4活动快捞柜-05
05 QUICKLY ARK
410*540*655mm



A4 活动快捞柜-06
06 QUICKLY ARK
410*540*655mm

DOOR
LOCKER

多门储物柜




FASION DESIGN

AND EXCELLENT QUALITY

Bring You Comfortable and Gladdergladdest Mood

潮流的设计，精湛的工艺，
带给你舒适愉悦的工作心情





CG-03 三门储物柜
Three Doors Locker
400*400*1800mm



CG-04 四门储物柜
Four Doors Locker
400*400*1800mm



CG-05 五门储物柜
Five Doors Locker
400*400*1800mm



TG-015 九门储物柜
Nine Doors Locker
980*400*1800mm



TG-051 十二门储物柜
Twelve Doors Locker
980*400*1800mm



TG-049 十五门储物柜
Fifteen Doors Locker
980*400*1800mm



TG-052 十八门储物柜
Eighteen Doors Locker
980*400*1800mm



TG-053 二十四门储物柜
Twenty Four Doors Locker
980*400*1800mm



TG-019 三十二门储物柜
Thirty Two Doors Locker
980*400*1800mm

DOORS
WARDROBE

钢制更衣柜



FASION DESIGN

AND EXCELLENT QUALITY

Bring You Comfortable and Gladderladdest Mood

潮流的设计，精湛的工艺，
带给你舒适愉悦的工作心情





YG-02 上下两门衣柜
Double Doors Wardrobe
400*480*1800mm



YG-03 两门衣柜
Two Doors Wardrobe
980*480*1800mm



YG-035 三门衣柜
Three Doors Wardrobe
980*480*1800mm



YG-16 四门衣柜
Four Doors Wardrobe
980*480*1800mm




YG-22 六门衣柜
Six Doors Wardrobe
980*480*1800mm



YG-21 八门衣柜
Eight Doors Wardrobe
980*480*1800mm

SWIMMING POOL
DOOR LOCKER

感应锁储物柜



(37EM感应锁系列)



(198EM感应锁系列)

GY-01 18门感应锁储物柜
980*400*1800mm

感应开锁 更方便快捷
OFFICE PASSWORD CABINET LOCK

防潮、耐用、易管理

Moistureproof、 Enduring、 Simple management

在努力改善产品质量的同时，我们不忘初心的充分考虑产品的安全性和耐用性。


- 1、 锌合金锁体，更耐用；
- 2、 独立宾客卡开门，万把锁不重复；
- 3、 客户自行更改宾客卡，方便管理；
- 4、 隐藏式外接电源（其它外接电源无法接入）
- 5、 防止暴力开锁，自带警报系统。



GY-01 9门感应锁储物柜
980*400*1800mm



GY-02 12门感应锁储物柜
980*400*1800mm



GY-03 24门感应锁储物柜
980*400*1800mm


相似产品



GY-04 9门感应锁更衣柜
980*400*1800mm



GY-05 12门感应锁更衣柜
980*400*1800mm



GY-06 24门感应锁储物柜
980*400*1800mm

DOORS

WARDROBE

感应锁更衣柜



感应开锁 更方便快捷

OFFICE PASSWORD CABINET LOCK

防潮、耐用、易管理


Moistureproof、Enduring、Simple management

在努力改善产品质量的同时，我们不忘初心的充分考虑产品的安全性和耐用性。


- 1、锌合金锁体，更耐用；
- 2、独立宾客卡开门，万把锁不重复；
- 3、客户自行更改宾客卡，方便管理；
- 4、隐藏式外接电源（其它外接电源无法接入）
- 5、防止暴力开锁，自带警报系统。

6门感应锁更衣柜

980*480*1800mm



GY-07 4门感应锁更衣柜
980*480*1800mm




GY-08 6门感应锁更衣柜
980*480*1800mm




GY-09 8门感应锁更衣柜
980*480*1800mm

相似产品



GY-10 4门感应锁更衣柜
980*480*1800mm



GY-11 6门感应锁更衣柜
980*480*1800mm



GY-12 8门感应锁更衣柜
980*480*1800mm

STEEL OFFICE
DESK

钢制办公台



FASION DESIGN

AND EXCELLENT QUALITY

Bring You Comfortable and Gladderladdest Mood

潮流的设计，精湛的工艺，
带给你舒适愉悦的工作心情





ZT-07 吊抽办公台

Hanging Drawers Office Desk

1000*500*750mm
1000*600*750mm
1200*600*750mm
1400*700*750mm



ZT-012N-1 单抽办公台

Single Side Office Desk

1000*500*750mm
1000*600*750mm
1200*600*750mm
1400*700*750mm



ZT-010N 双抽办公台

Double Side Office Desk

1400*700*750mm
1500*700*750mm



ZT014H 双人电脑台
DOUBLE PERSON COMPUTER DESK
1600*700*750MM



ZT01H 单人电脑台
SINGLE PERSON COMPUTER DESK
1000*600*750MM
1200*700*750MM



ZT02H 单人电脑台
SINGLE PERSON COMPUTER DESK
1500*700*750MM



ZT019H 单人电脑台
DOUBLE PEDESTAL COMPUTER DESK
1400*700*750MM



钢制办公桌-01

01 single side office desk

1200*600*750mm

1400*700*750mm

1600*700*750mm



钢制办公桌-02

02 pedestal office double

1200*600*750mm

1400*700*750mm

1600*700*750mm



钢制办公桌-03

03 single side office desk

1200*600*750mm

1400*700*750mm

1600*700*750mm



钢制办公桌-04

04 single side office desk

1200*600*750mm

1400*700*750mm

1600*700*750mm



钢制办公桌-05

05 single side office desk

1200*600*750mm

1400*700*750mm

1600*700*750mm



钢制办公桌-06

06 single side office desk

1200*600*750mm

1400*700*750mm

1600*700*750mm



L型钢制办公桌-01

01 single side office desk

1600W*700D*750H (主桌) 1000W*430D*750H (侧桌)



L型钢制办公桌-02

02 single side office desk

1600W*700D*750H (主桌) 1000W*430D*700H (侧桌)



L型钢制办公桌-03

03 single side office desk

1600W*700D*750H (主桌) 1000W*430D*750H (侧桌)



SERIES PRODUCT LIST 图书设备系列产品一览

1
SERIES



钢制书架/STEEL SHELF

产品特征：书架组和层高可任意组合，设计合理，兼容性好可安装多种规格插标、工艺精良，压痕均匀，外观好，其效果超越传统的螺钉链接书架强度和稳定性。拆装速度快，运输便捷，从而能达到节约成本和提高拆装速度以及安全适用的目的。

2
SERIES



全钢书架/ALL STEEL SHELVES

产品特征：全钢书架由立柱，侧板，搁板，层板，底座，挡棒，顶板组成。产品各部位的安装牢固、可靠、无松动、承载能力强、外光光滑平整，色泽均匀一致，线条流畅。组装结合严密、无缝隙产品各同类零件，组件之间能进行互换。整体效果整齐美观，端庄气派。

3
SERIES



防火板钢制书架/OAK GUARD A BOOKCASA

产品特征：优业防火板钢制书架采用组合装备式设计，搬迁和拆装均不受损失，木护板采用高档防火板精致而成，板面木纹清晰，彰显绿色自然深邃，整体效果端庄气派尽显古典文化和现代文明之内涵，更显品位与风范，倍增细腻美观的视觉效果。



4
SERIES



期刊架/ANCIENT BOOKCASE

产品特征：展板上可展示供阅览的期刊，展示板后可存放图书资料，有效利用空间，期刊架表面光滑平整，牢固可靠、美观实用、是学校、图书馆、档案馆存放展示图书资料的理想装具。

5
SERIES



阅览桌/CARREL

产品特征：桌架采用静电粉末喷涂，高温固化。先进的九工位表面处理工艺，优质的塑粉使产品的漆膜附着力、硬度、耐冲击性能、光滑度、平整度等符合国家标准。桌面为防火板贴面（进口优质无苯无毒无味白乳胶），遇到潮湿或温度较高的情况不鼓起，不变形。

6
SERIES



图书馆书梯/MUTE BOOK LADDER

产品特征：优业图书馆书梯采用0.8MM优质冷轧钢板，且采用静电粉末喷涂，高温固化。先进的九工位表面处理工艺，优质的塑粉使产品的漆膜附着力、硬度、耐冲击性能、光滑度、平整度等符合国家标准。无甲醛、无有毒物质。



SJ-01 双面五层全钢书架 W1800*D450*H2000



SJ-02 双面六层钢制书架 W1800*D450*H2000



SJ-03 双面六层钢制套色书架 W1800*D450*H2000




SJ-02 单面五层钢制书架 W2700*D450*H2000



SJ-03 单面钢木五层图书架 W2700*D450*H2000



SERIES PRODUCT LIST 图书设备系列产品一览



QK-01 单面期刊架
W900*D360*H1800MM



QK-02 仓储期刊架
W900*D360*H1800MM



QK-03 双面钢木六层活动期刊架
W1800*D650*H2000MM



YL-01 钢木阅览桌
W1500*D900*H7500MM (四人)
W1800*D900*H7500MM (六人)



YL-02 钢木阅览桌
W1500*D900*H7500MM (四人)
W1800*D900*H7500MM (六人)



YL-03 钢木阅览桌
W1500*D900*H7500MM (四人)
W1800*D900*H7500MM (六人)



PJ-01 平板书车
W750*D360*H900MM



PJ-02 书梯
W450*D650*H1300MM



PJ-03 标准书踏、书立
大号：160*235MM 中号：135*175MM



INTELLIGENT
ELECTRIC

智能密集架

产品概述

- 1、密集架为框架结构：由底盘、双柱、隔板及支承板、顶板、护板、门面、传动机构、道轨、防震装置及防尘、防鼠装置等部分组成。
- 2、节型、架高、架宽、层数和列数均可按用户需要选定。
- 3、轨道采用20*20方钢制作，有镀锌的，不锈钢管套的，可根据客户需求进行选择。
- 4、底架一般为分段组合式，加工精度高，具有对接互换性，便于安全运输和安装。
- 5、传动机构采用链条和精密平面轴，传动灵活平稳。
- 6、隔板与扣勾的支承板，组合后平稳牢固，层与层的间距根据需要可自由调整，每层隔板均受重，单面为50KG，双面可达到100KG，如需增大负重，可专门设计。
- 7、每列均装有手刹车制动装置，边列装有锁具，用于整体锁闭。
- 8、每列的接触面均有缓冲和密封装置，顶部有防尘板，底部有防鼠装置，具有良好的防尘、防鼠、防潮、防火功能。



系统优化设计 密集架存储

大容量，科学划分，更大限度提升储物效率



上门安装



支持定制



免费测量

INTELLIGENT ELECTRIC 智能密集架



MJ-01 开放式移动密集架

产品特性：此款密集架最经济实惠，能够满足所有办公环境下档案存储需求。可将空间利用率提升80%，层板可调节高度，摇动过程中轻便省力，可快速存取档案，且结构灵活多变可满足除档案外不同物品的存放要求。



MJ-02 全封闭式移动密集架

产品特性：此款密集架适用于对档案保密性有一定要求的客户，即便密集架制动锁具开启，还需要对应柜体的钥匙才能取出物品。



MJ-03 拉架式快捞密集架

产品特性：此款产品适用于文件需要通过快捞挂起的客户。



MJ-04 底图移动密集架

产品特性：此款产品适用于存放各类型设计图纸，可根据图纸大小定制，如：0号图纸、1号图纸等。

材 料	规 格	1	2
		抽屉内尺寸	1200*900*100
单体型	外形尺寸	1300*940*1160	960*740*1200
密集型	架 宽	940	740
	架 长	1300*节数+135	1300*节数+135
	架 高	1800 2000	2200



MJ-05 全封闭式移动密集架

产品特性：此款密集架适用于对档案保密性有一定要求的客户，即便密集架制动锁具开启，还需要对应柜体的钥匙才能取出物品。



MJ-06 全封闭式移动密集架

产品特性：此款密集架适用于对档案保密性有一定要求的客户，即便密集架制动锁具开启，还需要对应柜体的钥匙才能取出物品。

TECHNICAL PARAMETERS AND
STRUCTURAL ANALYSIS
密集架主要技术参数及结构解析

规 格						规 格					
结构型式	架宽	节柜	架高	层数	层高	节数	架长	道数轨	隔板承重		
双面活动架 (中列)	500	900	2400	5	398	2	1935	2	资料架 >60KG		
	550			6	330	3	2835	2			
	600			7	280	4	3735	3			
	650			8	205	5	4635	3			
	700			12	155	6	5535	4			
双面活动架 (中列)	800							7	6435	4	货架 >80KG
	350							8	7335	5	
	400							9	8235	5	
	450							10	9135	6	



WAREHOUSE
EQUIPMENT
SERIES

仓储设备系列



以灵感挑战传统 以细节成就未来

INSPIRATION CHALLENGE TRADITIONAL DETAILS OF ACHIEVEMENTS OF THE FUTURE


产品特点：

结构合理，拆装方便，可根据用户和物件大小增减托层，层高可任意调节，每层横梁均有锁定装置，防止操作不当导致层板脱落，引发安全事故；

适应行业：广泛用于商场、超市及企业，事业单位；

服务内容：

适用于存储样少量多，批量作业且选择性较小的货物。货物采用托盘一个接一个按深度方向，存放在支撑导轨上，增大存储密度和空间利用率。叉车可直接驶入巷道内存储货物，极为方便。





HJ-01 轻型货架
W2000*D500*H2000



HJ-02 轻型移动货架
W2000*D500*H2000



HJ-03 重型货架
W2000*D500*H2800



HJ-04 重型货架
W1500*D500*H2000

形于效率 而备于精细

规格：长X宽 (mm)		高 (mm)		承载kg/层			
1600*600		1800	2000	可根据客户需求定制			
1600*900							
1600*1180							
2000*600							
2000*900							
1600*1180							
重量型及中量型货架立柱及横梁规格示意图							
立柱横梁	立柱横梁						
	名称	型材	C材	P材	矩管	槽钢	横梁锁
	规格	70*90	40*65	60*40*20	40*60	40*80	30*70

STUDENT
APARTMENT
BUNK BED

组合公寓床



本产品采用高强度钢管制作，表面烤漆处理，不可接触带酸性的腐蚀液体，请用于干燥软布清洁铁床，请勿用酒精或其他清洁剂清洁铁床，安装过程中尽量避免使用影舞锤击。

EVERY
DETAIL REVEALS A FINE

没有复杂的元素
每个细节都透着精致





ZH-01 单人组合床
尺寸：2000W*980D*2100H mm



ZH-02 单人组合床
尺寸：2000W*980D*2100H mm



ZH-03 单人组合床
尺寸：2000W*980D*2100H mm



ZH-04 单人组合床
尺寸：2000W*980D*2100H mm



ZH-05 单人组合床
尺寸：2000W*980D*2100H mm



ZH-06 单人组合床
尺寸：2000W*980D*2100H mm

APARTMENT COMBINED BED
双人连体组合床



SQ-01 双人组合床
尺寸：4000W*980D*2000H mm



SQ-02 双人组合床
尺寸：4000W*980D*2000H mm



SQ-03 双人组合床
尺寸：4000W*980D*2000H mm



SQ-04 双人组合床
尺寸：4500W*980D*2000H mm



SQ-05 双人组合床
尺寸：4000W*980D*2000H mm



SQ-06 三人组合床
尺寸：4500W*980D*2000H mm

STUDENT DESK
AND CHAIR

课桌椅系列



双人课桌椅

规格 桌面：W1200*D500*H700mm

规格 椅面：W380*D375*H440mm

材质 桌面E0级环保板，一体全注塑封边
椅面PP材质，中空吹塑，钢架静电喷塑





KZ-01 单人升降注塑课桌椅

规格：桌子：650*450*700mm
椅子：400*360*700mm



KZ-02 单人升降注塑课桌椅

规格：650*450*(680-760)mm



KZ-03 单人胶合板升降课桌椅

规格：600*400*(65-75)mm



KZ-04 单人胶合板升降课桌椅(无靠背)

规格：600*400*(65-75)mm



KZ-05 单人注塑升降课桌椅

规格：桌子：400*600*(630-760)mm
椅面：380*360*(360-420)mm



KZ-06 单人注塑升降课桌椅(无靠背)

规格：桌子：400*600*(630-760)mm
椅面：340*240*(360-420)mm



KZ-07 双人注塑封边课桌椅

规格：桌子：1200*500*700mm

椅子：380*375*440mm



KZ-08 双人升降注塑课桌椅

规格：桌子：1200*420*730mm

椅子：350*410*740mm



KZ-09 双人注塑升降课桌椅（无靠背）

规格：桌子：1100*400*（630-760）mm

椅面：340*240*（360-420）mm



KZ-10 培训折叠双人连体桌

规格：1050*400*750mm



KZ-11 双人胶合板升降课桌椅

规格：1050*400*780mm



KZ-12 双人高档注塑课桌椅

规格：桌子：1200*420*730mm

椅子：350*410*740mm



CLASSROOM SEATING

阶梯教室排椅

基因决定一切

从经久耐用的品质，到简约实用的设计，优业家具的名字已经与力量、耐用性、实用设计理念连接在一起。流金岁月中的优业，一直在改变，却又宛如初见。

我们不断征服影响家具寿命和高成本的因素，完美实现了耐用，使用以及高性价比的结合，让我们赢得更多的信任。

优业因专业而领先
注定为卓越品质而诞生
DESTINED BIRTH FOR
EXCELLENT QUALITY



AMPHITHEATRE CHAIR 阶梯教室排椅



JT-01 多媒体阶梯教室椅
规格：1040*840*840MM



JT-02 多媒体阶梯教室椅
规格：1040*840*840MM



JT-03 多媒体阶梯教室椅
规格：1040*840*840MM



JT-04 多媒体阶梯教室椅
规格：1040*840*840MM



JT-05 多媒体阶梯教室椅
规格：1040*840*840MM



JT-06 多媒体阶梯教室椅
规格：1040*840*840MM

部分项目案例





DINING TABLE
CHAIR

食堂餐桌椅

COLORFUL LIEF STYLE 多彩生活·简约格调



- 1、本餐桌分为4人位、6人位、8人位、可拆装餐桌，是公司、工厂食堂、员工、快餐店、学校、部队等公共就餐区的理想选择！
- 2、产品采用对不同使用环境分别采用户外玻璃钢和室内防火板，使产品使用寿命有效提高。
- 3、耐热、耐磨、耐划，独特的张力，超强承重，不易变形，稳重兼具美感。
- 4、主材料用不饱和树脂和玻璃纤维布经固化等化学反应复合成玻璃钢中间夹一层15mm厚高密度木板，表面层涂两次耐磨胶衣，成品厚度达到17~18mm，有耐磨、耐腐蚀、抗冲击、防水，不裂缝等优点。
- 5、材 质：支架采用47*47*1.1mm厚优质铁方通，烤漆，台面与凳面为玻璃钢材质。
- 6、本产品具有耐冲击、抗老化、抗氧化等特点，方便清洁、色调明快符合国际环保要求。



CZ-01 长条四人不锈钢餐桌
规格：1430*750*1100MM



CZ-02 长条四人不锈钢餐桌
规格：1430*750*1100MM



CZ-03 长条四人玻璃钢餐桌
规格：1430*750*1100MM



CZ-04 四人玻璃钢圆凳餐桌
规格：1430*750*1100MM



CZ-05 四人分体曲木餐桌
规格：1650*750*1200MM



CZ-06 四人分体曲木餐桌
规格：1650*750*1200MM



CZ-07 双人曲木餐桌
规格：1650*750*600MM



CZ-06 四人分体曲木餐桌
规格：1650*750*1200MM



TRAINING
CHAIR

培训椅系列

高靠背紧密贴合脊背曲线，人体工程学设计缓解久坐背部压力，集美观与舒适为一体

HIGH BACK FIT CLOSELY BACKBONE CURVE, EASE SITTING BACK PRESSURE OF HUMAN BODY
ENGINEERING DESIGN, BEAUTIFUL AND COMFORTABLE AS A WHOLE





ZD-01 折叠培训椅
规格：830*550*485MM



ZD-02 折叠培训椅
规格：830*550*485MM



ZD-03 培训椅
规格：755*530*520MM



ZD-04 培训椅
规格：750*580*680MM



ZD-05 培训椅
规格：785*540*520MM



ZD-06 培训椅
规格：825*440*580MM

用心服务

尊重客户，理解客户

持续提供超越客户期望的产品与服务

做客户永远的伙伴。

这是我们一直坚持和倡导的服务理念



HONESTY DECLARATION

诚信宣言

诚信就是真心实意，诚实守信。它是立世之本，处世之道，创业之基。

诚信就是鱼人交往，要诚信相处，诚信相交，诚心相待，不搞虚情假意，不搞使用主义。

诚信就是诚恳，守信，言必言，行必果，说话算数，兑现承诺。

诚信要坚守职业道德，宁可失去金钱，决不失去信誉。

用“一切为了用户，为了用户一切，为了一切用户”作为座右铭。

“把产品当作人品做，把产品当作良心做”。

切实为用户提供规范、标准、优质、可靠的全程无忧服务。用诚信这张名片，领取进入市场的通行证。

YOUYE

ENTERPRISE SERVICES
SERIES

企业服务篇

成就如星 实力筑辉煌

天道酬勤，成绩斐然，
力求市场、艺术与技术的完美结合，
以精品呈现实力，以实力铸就辉煌！

15年来积极进取，成就了今天的基业；8年来求索展翅，优业豪情满怀。在15年来的发展中，优业创下了无数的精品工程，尤其是进8年以来，更是以飞速发展的姿态，进行着一次又一次的辉煌蝶变。

God rewards the diligent , Have made brilliant achievements,
Strive to market,the perfect combination of art and technology,
Strength boutique rendering,
Stength Brillianteness!

15 years of positive progress and achievement of today's base industry;
In the past eight years, we have sought to spread our wings. In theOver
the past 15 years, the company has set up numerous boutique projects,
especially in the past eight yearsThe rapid development of the posture,
carry on the brilliant butterfly change again and again.

实力铸就辉煌
成就如星

部分长期合作单位



长隆集团员工更衣储物柜供应商



白云机场档案密集柜、候机椅供应商



东风日产集团员工更衣柜、密集柜供应商



中国光大银行
CHINA EVERBRIGHT BANK
广州光大银行档案密集柜供应商



广东岭南职业技术学院

岭南职业技术学院学生宿舍公寓床供应商



中山大学
SUN YAT-SEN UNIVERSITY

中山大学学生宿舍公寓床、图书馆图书架供应商



广州地铁文件柜、储物柜供应商



广州塔游客休息椅、文件柜供应商



广东财经大学学生宿舍公寓床供应商



广州大学附属中学课桌椅、学生宿舍公寓床供应商



美的集团广州区域员工储物柜文件柜供应商



广州市中医院档案密集柜、钢制中药柜供应商



携手成就未来

CREATE THE FUTURE HAND IN HAND



广东优业实业有限公司

[Http://www.youyef.com](http://www.youyef.com)

座机：020-87451308 传真：020-87450530

地址：广东省广州市白云区钟落潭镇金盆村金埔南路23-25号

(图片仅供参考，以实物为准)

Brief of Youye Furniture

Youye Furniture Co., LTD is a professional furniture manufacturer, which is located in Baiyun District, Guangzhou, only 40 minutes course to Guangzhou airport, Train Station, Province Bus Station, City Bus Station and Huangpu Port.

Youye, established in 2005, has been developed into a large-scale furniture manufacturer, which main product include: School furniture, Student desk and chair, Apartment Furniture, Bunk Beds, Wardrobe and Writing table. Office working station, Filing Cabinet, Drawer Filing Cabinet, Rolling Door Filing Cabinet , Locker, Wardrobe, Computer Desk, Office Desk, Filing Mobile Archive Cabinet as well as various kinds of metal structure series product. The factory covers 20 mu (13000 square meters), the workshop is about 10000 square meters, in which there are 200 equipment for stamping, bending, cutting, and welding etc. The advanced powder-coating was used in coating workshop to prevent products stained and rusted. 3000 square meter warehouse stores 2000 sets products of all kinds of regular model, enough to meet demands for supplying domestic and foreign distributors in time.

Youye adheres to the business philosophy of “Practical, Creative, Flexible, Excellent” constantly developing the trend-leading new products in line with demand of modern office. We share with you with win-win benefit base on our professional design and perfect technique process.

公司地址地图导航

