

## CATARACT PACK

<b>Product CHECK LIST</b>	<b>Description / Status</b>	<b>Quantity</b>
Swab	Non-Woolen (10x10cm 4 ply)	10 pcs
Tray	Half Shallow (Plastic)	01 pcs
Gallipot	Polypropylene (120ML)	02 pcs
Kidney Dish	Polypropylene (800ML)	01 pcs
Needle Blunt Fill Filter	Size (18g/ 19g x1.5")	01 pcs
Syringe	5 ML Luer lock	01 pcs
Syringe	3 ML Luer lock	03 pcs
Syringe	1 ML Luer lock	01 pcs
Scissors Conjunctival	Stainless Steel Single use	01 pcs
Scissors West-Cott	Curved Round Tip (SS)	01 pcs
Drysdale Nucleus Spatula	Paddle Shape(Stainless Steel Single-use)	01 pcs
Speculum Lid	Stainless Steel Single use	01 pcs
Eye Shield	Clear (Plastic)	01 pcs
Drape Trolley Cover	Size (140cm x 190cm)	01 pcs
Paper Crepe white	Size (90cm x 90cm)	01 pcs
Drape Ophthalmic	Size(140cmx127 cm)with Incise& Pouch	01 pcs
Forcep Sponge Holding	Plastic Single Use (19cm 7.5")	01 pcs
Hoskins Forcep	Stainless Steel Single use	01 pcs
Hydrodissector (Akaoshi)	Angled needle (27G x7/8 0.4 x 22mm)	01 pcs
VIS Cannula Anterior	Chamber Rycroft ( 27G x 7/8")	01 pcs
Cannula Cystotome	Irrigating Blumenthal (27G x5/8")	01pcs
Slit Knife	Single bevel up angled (2.2mm)	01 pcs
MVR Knife	Side-port angled (1.0mm)	01 pcs
<b>Total</b>		<b>37 pcs</b>

