

Led Operation Lamp

YDE700/500

A Type



- Streamlined, ultra-thin, enclosed head with reduced air resistance to minimize contamination
- Manual /Automatic focus
- Digital controlled light intensity with 10 grades of brightness and automatic memory for added flexibility
- Built in power supply systems avoiding complicated installation of wall mounted control boxes
- Auxiliary light source included for laparoscopic operations
- Rigid suspension arm with easy mobility and stable positioning
- No infrared or ultraviolet radiation

PARAMETERS:

Diameters:	700mm	500mm
Illuminance	160000Lux.	140000Lux.
Color Temperature	3600--5500K(±500K)variable	
Color Rendering Index	> 95	> 95
Ultraviolet light irradiance within400nm wave length:	<0.006W/m2	
Power Of Bulb(W)	1W X N60	1W X N45
50% Spot Diameter d50(mm)	150	120
10% Spot Diameter d10(mm)	280	230
Depth Of Illuminance (cm)	140	120
The Adjusting Range Of Illuminance	25-100%	
Focus System	Electric / Manual (optional)	
Voltage	220V	220V
Frequency	50Hz	50Hz
Input Power	100W	80W
LED Service Life	100,000h	
The Temperature On Surgeon's Head	Is Less Than 1°C	
Camera System	Optional	
Battery	Optional	

YDE500



Page : 1

Led Operation Lamp



YDE700/500
B Type

- Streamlined, ultra-thin, enclosed head with reduced air resistance to minimize contamination
- Manual /Automatic focus
- Digital controlled light intensity with 10 grades of brightness and automatic memory for added flexibility
- Built in power supply systems avoiding complicated installation of wall mounted control boxes
- Auxiliary light source included for laparoscopic operations
- Rigid suspension arm with easy mobility and stable positioning
- No infrared or ultraviolet radiation

PARAMETERS:

Diameters:	700mm	500mm
Illuminance	160000Lux.	140000Lux.
Color Temperature	3600—5500K(±500K)variable	
Color Rendering Index	> 95	> 95
Ultraviolet light irradiance within400nm wave length:	<0.006W/m2	
Power Of Bulb(W)	1W X N60	1W X N45
50% Spot Diameter d50(mm)	150	120
10% Spot Diameter d10(mm)	280	230
Depth Of Illuminance (cm)	140	120
The Adjusting Range Of Illuminance	25-100%	
Focus System	Electric / Manual (optional)	
Voltage	220V	220V
Frequency	50Hz	50Hz
Input Power	100W	80W
LED Service Life	100,000h	
The Temperature On Surgeon's Head Is Less Than 1°C		
Camera System	Optional	
Battery	Optional	

YDE500D



Led Operation Lamp



YDE700/500

Slim Design

- Utilising the most advanced LED technology available, the MT-medical LED70 series of operating theatre lights have been designed to provide a high quality, shadowless beam with extremely low running costs.
- The LED operation lamp has the backlight compensation function, one parts of LEDs over surgeon's head will be closed and Other Light source will be enhanced.
- LEDs (light emitting diodes) have been around for a long time. Recently this technology has become leading -edge in the field of lighting - a breakthrough technological solution to the inherent problems of the incandescent light bulb.
- We have taken the technology and developed it further to provide high definition lighting solutions with no infrared, low environmental impact, high brightness, low power consumption, no toxic chemicals and a long operating life.

PARAMETERS:

Diameters:	700mm	500mm
Illuminance	160000Lux.	140000Lux.
Color Temperature	3600---5500K(±500K)variable	
Color Rendering Index	≥97	≥97
R9 : 95	R13 : 92	
Power Of Bulb(W)	1W X N128	1W X N64
Spot Diameter d50(mm)	130-350	
Depth Of Illuminance (cm)	140	120
The Adjusting Range Of Illuminance	10-100%	
Focus System	Manual	
Minium Installation Height is	2800mm	
Voltage	220V	220V
Frequency	50Hz	50Hz
Input Power	160W	80W
LED Service Life	100,000h	
The Temperature On Surgeon's Head Is Less Than	1°C	
Camera System	Optional	
Battery	Optional	



YDE500

Slim Design

Halogen Operation Lamp



YDZ700/500

- shadowless operating lamps meet the requirements of various surgical operations. It is an ideal lighting facility for surgical operations.
- Special heat-isolated radiation filter reduces heat to the minimum level.
- Imported aluminum sheet squamulose reflector and transparent panel are adopted. Average service life of imported lamp receptacle and bulb can reach 1500 hours.
- The imported aluminum reflecting mirror makes light softer and more well-proportioned and increasing safety. CPU is adopted for the control, capable of increasing or decreasing and boasting brightness memory. When the main bulb malfunctions, the standby bulb will come into operation within 1/2 second. The major and auxiliary bulbs can be replaced easily. The dischargeable handle sheath can be sterilized in steam sterilizer.

PARAMETERS:

Diameters:	700mm	500mm
Illuminance	150000Lux.	120000Lux.
Color Temperature	3000K---6700K	
Color Rendering Index	≥95	
Spot Diameter	120--320mm	
Depth Of Illuminance	≥800mm	
Light-adjusting mode	10 grade dimmer system	
Focus System	Manual	
Temperature At Surgeon's Head	<2°C	
Temperture At Surgical Field	<8°C	
Working Distance:	55-170cm	
Power Supply	220V±22V 50Hz±1Hz	
Input Power	200W	
Power Of Bulb	24V/150W	
Bulb average Life	1500h	
Shift Of Major And Stand-by Lamps	YES	
Camera System	Optional	
Battery	Optional	



YDZ700



YDZ500D

Led Operation Lamp

YDE700/700

C Type



- Digital controlled light intensity with 10 grades of brightness and automatic memory for added flexibility
- Built in power supply systems avoiding complicated installation of wall mounted control boxes
- Auxiliary light source included for laparoscopic operations
- Rigid suspension arm with easy mobility and stable positioning
- No infrared or ultraviolet radiation

PARAMETERS:

Diameters:	700mm	700mm
Illuminance	160000Lux.	160000Lux.
Color Temperature	3600—5500K(±500K)variable	
Color Rendering Index	>95	>95
R9: 95 R13: 92		
Power Of Bulb(W)	1W X N80	1W X N80
50% Spot Diameter d50(mm)	150	150
10% Spot Diameter d10(mm)	280	280
Depth Of Illuminance (cm)	140	140
The Adjusting Range Of Illuminance	25-100%	
Focus System	Electric / Manual (optional)	
Voltage	220V	220V
Frequency	50Hz	50Hz
Min.installation height	2800mm	
LED Service Life	100,000h	
The Temperature On Surgeon's Head	Is Less Than 1°C	
Camera System	Optional	
Battery	Optional	



YDE500 D

Led Operation Lamp



YDE700 /700 D Type

- Streamlined, ultra-thin, enclosed head with reduced air resistance to minimize contamination
- Manual /Automatic focus
- Digital controlled light intensity with 10 grades of brightness and automatic memory for added flexibility
- Built in power supply systems avoiding complicated installation of wall mounted control boxes
- Auxiliary light source included for laparoscopic operations
- Rigid suspension arm with easy mobility and stable positioning
- No infrared or ultraviolet radiation

PARAMETERS:

Diameters:	700mm	700mm
Illuminance	160000Lux.	160000Lux.
Color Temperature	3600—5500K(±500K)variable	
Color Rendering Index	≥95	≥95
R9 : 95 R13 : 92		
Power Of Bulb(W)	1W X N64	1W X N64
50% Spot Diameter d50(mm)	150	150
10% Spot Diameter d10(mm)	280	280
Depth Of Illuminance (cm)	140	140
The Adjusting Range Of Illuminance	25-100%	
Focus System	Electric / Manual (optional)	
Voltage	220V	220V
Frequency	50Hz	50Hz
Min.installation height	2800mm	
LED Service Life	100,000h	
The Temperature On Surgeon's Head	Is Less Than 1°C	
Camera System	Optional	
Battery	Optional	



YDE700 E Type

Portable Led Examination Light



YDE300

Wall / Ceiling mounted type

- YDE300 Portable LED examination light is an invaluable and cost effective tool. It has been specifically designed for minor procedures under local anaesthetic in health centres, GP surgeries and outpatient departments.
- FEATURE :
 - The use of multiple LED's provides a high, and even light intensity throughout the entire viewing area to reduce eye strain
 - A concise and self-explanatory 3 button control system on the lightthead itself
 - Balanced arms for one touch adjustment and a great range of movement
 - Ideal for minor surgery and gynaecological procedures
 - Removeable, Autoclavable Handle

PARAMETERS:

Diameters:	350mm	
Illuminance	80000Lux	
Color Temperature	4500K(±500K)variable	
Color Rendering Tndex (Ra)	≥ 90	
Input Power	50 W	
Power Of Bulb(W)	1W X N20	
Spot Diameter d50(mm)	80	
Depth Of Illuminance (cm)	80	
Battery	Optional	
Focus System	Automatic	
LED Service Life	100,000h	
Voltage	220V	220V
Frequency	50Hz	50Hz
Packing Details	0.3 CBM/unit, G.W:42 KGS	



YDE300

Portable Type

Portable Examination Light

YDE200D Deep Irradiation Exam Light

Lamp head diameter(mm)	180
Illuminance(Lux)	> 50,000
Color temperature(K)	4,000K±500
Color rendering index(Ra.%)	> 88
Spot diameter(mm)	160-220
Depth of illuminance(mm)	600
Power of bulb(W)	50
Bulb service life(hours)	1,000
Voltage	220V 50Hz or 110V 60Hz
Backup battery	optional
Packing details	11 CBM/set, G.W.:20KGS



YDE180D Portable LED Exam Light

Illumination	50000Lux
Color temperature	5500±500K
Voltage	AC90-240V
Bulb Power	LED 7*3W
LED service life	50000h
Power	7*3 W LED
Packing details	0.03 CBM/set,G.W.:85KGS

Hole Type Exam Light

Illuminance	> 40,000Lux
Color temperature	3,700K±500K
Color rendering index(Ra.%)	> 85
Spot diameter(mm)	160-220
Depth of illuminance(mm)	900
Power of bulb(W)	25
Bulb service life(hours)	1,000
Voltage	220V 50Hz
Backup battery	optional

YDK04D



Portable Led Examination Light

YDE200 LED Portable LED exam light

Illuminance	60000Lux.
Color Temperature	4800 ±300 K
Power supply voltage	AC 220V/50-60Hz / optional 110V/60Hz
Power of bulb	3.2V/1W
Diameter of light head	180mm



YDE200C Ceiling LED exam light

Illuminance	40000Lux.
Color Temperature	5000±500 K
Input Power:	36 W
Voltage	AC90-240V
LED Service Life	50,000h



YDE-A LED exam light



YDE-B LED exam light



YDE10 LED exam light

Imaging Operation Table


DST-700G
Imaging Orthopedic Operation Table

DST-700G Electric orthopedic navigation image operation table, for the expansion of minimally invasive treatment to improve surgical accuracy with modern medical diagnostic imaging technology, stereotactic designed and product design ergonomic characteristics of human anatomy, application of modern technology and user-friendly design philosophy, can meet the requirements of neurosurgery, spinal surgery, ENT, plastic surgery, radiosurgery and 3D imaging. The operation table adopts the imported micro-motor, with electric gear mesh mechanism structure instead of the traditional hydraulic transmission. The table top adopts the carbon fibre material, with low x-ray absorption, reduces the radiation damage of x-ray to the patients and paramedic. 1150mm length for table top without metal (full carbon fibre), meet the requirement for 3D and 360° X-ray photograph of orthopedics navigation systems.

PARAMETERS:

Bed panel size(Carbon fibre):	(2400*530)±10mm	One key reset function
Max Height :	1000±20mm	Battery Optional
Min Height:	700±20mm	Standard accessories
Route of longitudinal movement :	280±20mm(Electric)	Y-001 Mattress
Head tilt angle:	25°±2° (Electric)	Y-002 Controller
Foot tilt angle:	25°±2° (Electric)	Y-003 Lateral support
Left side leaning angle:	21°±2° (Electric)	Y-004 Carbon fibre arm support
Right side leaning angle:	21°±2° (Electric)	Y-005 Fixed belt
Back section up:	65°±2° (Electric)	Y-006 Anaesthesia screen
Back section down:	20°±2° (Electric)	Y-007 Leg support
Leg section up:	35°±5° (manual)	Y-008 Carbon fibre head frame
Leg section down:	90°±5° (manual)	Y-008 Detachable edge rail
Head section up:	45°±5° (manual)	P-001 Lock block/clamping block
Head section down:	90°±5° (manual)	
Bed body locking:	Electric	
Power supply (dual purpose) :	AC:220V 50Hz 300VA /DC:24V	
Net weight:	300kg	

Bed panel X-ray attenuation coefficient(Optional): ≤0.44mmAL

Electric Operation Table



DST-700

Luxury Electric Operation Table

Application of modern technology, combine with the idea of humanized design,meet the requirements of operation in the process of video observation and diagnosis,comprehensive strength,for each part of the human body surgery hospital operating room,i.e.,the brain,heart,blood vessels,extracorporeal circulation,abdominal cavity,the gallbladder excision,renal resection,the neck, cervical vertebra,lumbar spine,spine,back,urology,gynecology,cranial nerve,surgical,rectum,plam traction,upper limb surgery,knee surgery,orthopedics lower limb traction (such as with orthopaedics traction frame) surgery.

PARAMETERS:

Length	2220 mm	Back section down	22°
Total width	600mm(with side rail)	Trendelenburg	24°
Head section length	300mm	Re-trendelenburg	24°
Back section length	490mm	Lateral tilt	22°
Main section length	540mm	Leg section adjustment	Up 15° Down 90°
Leg section length	580mm	Spread angle of leg plates	200°
Supporting bottom size	1030x650x200mm	Head section adjustment	Up 40° Down 90°
Castor size	125mm	N.W. Of overall table	280kgs
Height adjustment range	700-1000mm	Max loading weight	320kgs
Kidney bridge rise	≧110mm	Patient Horizontal max loading weight	300kgs
Sliding movement	350mm	Integrated battery	Yes,for over 150 movements
Back section up	70°	Voltage	220V 50Hz or 110V,60HZ for choice
Standard accessories: Leg support-1 pairs, Arm supports-1pairs		Shoulder and waist holders- 2pairs, Anesthesia shelf-1pcs	

Electric Operation Table



DST-500A

Electric Operation Table

*This type electric operation table has advantages of beautiful appearance, convenient, labor-saving operation, safe and reliable use, has the advantage of smooth action, etc.. The platform lifting, roll, tilt, before and after the lock is realized by the special motor, push-button. With the C type arm, no dead angle, parts can be disassembled, easy to use.

* product is easy to operate, strong versatility, wide application, material selection, workmanship excellent, quality reached the international level.

*Head section: pneumatic auxiliary device for head plate, plate easy to operate.

*Bed section: special composite material plate, X-ray high transmittance, good uniformity, meet the requirements of operation.

*Mattress: Memory foam, it adopts the international popular material, streamlined appearance.

*Base: thin base design, hidden castors, both to ensure enough ground contact, to ensure stable operation platform, but also can be the greatest degree of operation in the C-arm film.

PARAMETERS:

Length	2040mm
Max Height	920mm
Min Height	670mm
Width	520mm
Reversed trendelenburg	> 20°
Trendelenburg	> 20°
Lateral Tilt	> 19°
Head board fold upward	> 20°, detachable
Head board fold downward	> 90°, detachable
Up and down range	250mm
Leg board fold downward	> 90°, detachable

Leg board open	> 90°, detachable
Back board fold upward	> 65°
Back board fold downward	> 30°
Left Tilt	> 19°
Right Tilt	> 19°
Kidney bridge elevation	350mm
Power Supply	AC 220±22v, 50HZ±1HZ
Translation moving	> 300mm±20mm
Battery	Optional
Standard accessories: Leg support-1 pairs, Arm supports-1pairs	
Shoulder and waist holders- 2pairs, Anesthesia shelf-1pcs	

Electric Operation Table



DST-1A
Electric Operation Table

This type operating table is suitable for operation of any parts of human body. Besides all the functions of multi-purpose operating table, it has the unique advantages of its own. For example, the leg panel can be disassembled and easily adjusted which is excellent for urology operation; the table board can be translated by a great degree and the high strength transparent board is adopted on the table panel, which is very convenient for X ray examination; the back panel can be upforward folded, which can overcome the defect caused by translation. All the movement of the table is controlled by electric drive pusher with safety and reliability. This operating table can be matched up with the osteology support especially for C-arm.

PARAMETERS:

Table size	1950x550mm
Height adjustable range	660-910mm
Trendelenburg	$\geq 20^{\circ} \pm 2^{\circ}$
Reversed trendelenburg	$\geq 20^{\circ} \pm 2^{\circ}$
Lateral tilt	left $\geq 20^{\circ} \pm 2^{\circ}$ & right $\geq 20^{\circ} \pm 2^{\circ}$
Table surface Sliding range	350mm
Head Section adjustable range	up $\geq 45^{\circ}$ & down $\geq 90^{\circ}$
Safe bearing load	200kg
Back Section adjustable range	up $\geq 70^{\circ}$ & down $\geq 20^{\circ}$
Leg Section adjustable range	up $\geq 45^{\circ}$ & down $\geq 90^{\circ}$
Kidney bridge elevation	≥ 120 mm
Power supply	220V 50Hz
Standard accessories: Leg support-1 pairs, Arm supports-1pairs	Shoulder and waist holders- 2pairs, Anesthesia shelf-1pcs
Battery	Optional

Electric Hydraulic Operation Table



DST-III
Electric Hydraulic Operation Table

*DST-III, the application of modern and humanism combined with the design concept, using electric imported hydraulic transmission technology, follow the development trend of the clinical operation, meet the image observation and diagnosis of the current operation process requirements.

*product is easy to operate, strong versatility, wide application, material selection, workmanship excellent, quality reached the international level.

*Fully meet the general surgery, heart, kidney, thyroid, extracorporeal circulation, head and neck, gallbladder and pleural, department of urology, department of orthopedics, department of gynaecology, rectum, spine surgical operation requirements.

Features:

- * The use of electric hydraulic transmission structure, stable speed, reliable performance;
- * The bed panel is a transparent material, the mattress adopts artificial wrapped high density memory foam;
- * The Lifting, translation, head and foot tilt operation driven by independent system.
- * Imported precision linear guide rail, ensure smooth running operation bed. 304 stainless steel material, the stainless steel cover, anti glare.
- * The design of large casters, mobile and flexible;
- * There are two sets of independent electric operating system, a remote control, a bed body standby operation control panel system, ensure the operation bed can still operate reliably in the remote controller failure.

PARAMETERS:

Size	2015 x 540 x 650-1000mm
Reversed /trendelenburg	25°/25° electric
V & ^ function	yes electric
Head section angle	up 30° down 90° manual
Leg section	up 28° down 90° manual
Tilt angle	left 22° right 22° electric
Back section	up 70° down 32° electric
Table panel x-ray attenuation coefficient < 0.76mmAL	
Power Supply	AC 220±22v, 50HZ±1HZ
Translation moving	>350mm±20mm electric
Battery	optional
Standard accessories: Leg support-1 pairs, Arm supports-1pairs	
Shoulder and waist holders- 2pairs, Anesthesia shelf-1pcs	
One key reset function	optional

Electric Hydraulic Operation Table



DST-III(B)

Electric Hydraulic Operation Table

- *.The controller device lock key protecting from misoperation.
- *.The table is made of penetrable X-ray material, there is cassette guide rail under the table which make radiography.
- *.The table is fixed a button which can control the electric wheel to draw in and put out so that the table can be pulled easily.
- *.Accessories and guide rail are made of stainless steel, polyurethane material etc.
- *.The main accessory hydraulic pump and electromagnetic valve are made of imported component from Taiwan.
- *.The mattress is made of form polyurethane, acid and alkali resistance, hight emperature resistance and easy clean.
- *.The head board and leg board are adopted special imported gas-spring from Germany, it's convenient operate.

PARAMETERS:

Table size	2100x500mm
Height adjustable range	750-1000mm
Trendelenburg	> 22°
Reversed trendelenburg	> 22°
Lateral tilt	left> 15°& right> 15°
Table surface Sliding range	350mm
Head Section adjustable range	up> 45°& down> 90°
Safe bearing load	240kg
Back Section adjustable range	up> 45°
Leg Section adjustable range	up> 45°& down> 90°
Kidney bridge elevation	> 120mm
Power supply	220V 50Hz
Standard accessories: Leg support-1 pairs, Arm supports-1pairs	Shoulder and waist holders- 2pairs, Anesthesia shelf-1pcs
Battery	Optional

Hydraulic Operation Table

DST-IB

Hydraulic Operation Table



HOT-H multi-purpose operation table is an improvement product developed on the basis of 3008 series head-manipulation type multipurpose operation table and is for use of general operations in chest surgery, department of ophthalmology, department of otolaryngology, department of obstetrics and gynecology, urological department etc. As a unit of non electronic control system comprehensive operation equipment, it is convenient and flexible to operate.

Features:

- * A multifunctional operation table that applies to various surgical operations and examinations.
- * Adoption of a foot hydraulic pump that integrates oil tank of oil pump and drive mechanism.
- * Adoption of the imported cylinder that controls actions of table top and is convenient to operate and reliable to locate.
- * Compact and rational structure and originative model.
- * Built-in waist bridge is convenient to use.
- * Both guide rail and accessories are made of corrosion resistant stainless steel.
- * All mattresses use the new-style materials that are acid and alkali resistant and disinfection durable.

PARAMETERS:

Size	1960 x 480 x 690-930mm
Reversed /trendelenburg	25°/25°
Lateral tilt	left ≥ 18° & right ≥ 18°
Head section angle	up 45° down 50°
Leg section	down 90°
Kidney bridge elevation	≥ 120mm
Back section	up 75° down 25°
G.W	295KGS
Power Supply	AC 220±22v,50HZ±1HZ
Translation moving	>350mm±20mm
Battery	optional
Standard accessories: Leg support-1 pairs, Arm supports-1pairs	
Shoulder and waist holders- 2pairs, Anesthesia shelf-1pcs	

Ophthalmology Operation Table

DST- I

Hydraulic Operation Table (head operating type)

- 1.The operating table is supplied for thoracic surgery, abdominal surgery, ophthalmology, obstetrics and gynecology, urinary surgery.
- 2.The rise and drop of operating table was handled by pedaling the oil pump, other various movements were handled by the mechanical transmission head.



Table size	2100 x 480 X 800 / 1045 mm
Trendelenburg & reversed trendelenburg	> 35° & > 55°
Lateral tilt	left > 25° & right > 25°
Head Section adjustable range	up:90° & down:90°(p.s.detachable)
Back Section adjustable range	up:85° & down:10°
Kidney bridge elevation	≥ 120mm
Leg Section adjustable range	down:90°(p.s.detachable)
Standard accessories: Leg support-1 pairs, Arm supports-1pairs	Shoulder and waist holders- 2pairs, Anesthesia shelf-1pcs

DST-II

Electric Ophthalmology Operation Table



PARAMETERS:

Size:	1950 x 520 mm
Height adjustable range	500 - 800 mm
Head Height adjustable range	550 mm
Power supply	AC220V , 50Hz
Standard accessories: Anaesthesia screen , Arm supports-1pairs , Remote controller	

Imaging Operation Table

JR-9000

Interventional Imaging Table



JR-9000 Catheterization table, adopt the advanced electric, mechanic transmission technology, full carbon fibre table top, specially design for various movable medium, small C-arm to imaging of the Cardiovascular, peripheral-vascular, catheterization vascular ang Jography and DSA. The product conform the features of logical structure, easy control, graceful outlook and the function is advanced and achieved the quality of Jatamalional. According to the clinical practice, it is confirmed as the best quality products in China.

Features:

- * The operation table adopts the imported micro-motor, with electric gear mesh mechanism structure instead of the traditional hydraulic transmission. Not need to change oil termly, and worry about the problem of leak oil. Product with fine steel can move equally and smoothly in any directions; also provides accurate positioning, trusty capability and durable lifetime.
- * Big gear feet wheel (the diameter of large truckle) especially use for medical treatment, the whole table is easily to be moved to clear the floor.
- * Test the table according to the CE European standard by Germany TUV company, achieve the European standard, obtain CE certificate. The operation table achieves the quality of imported products.
- * High precision imported beeline navigation rail, ensure the smooth and balanced movement of operation table.
- * Base with the ABS industrial plastic, helps prevent glare, impact and it is easily clear up and listerize.
- * Handle control is easily operating and takes control of the horizontal move of table top, the brake system is safe and sensitive.
- * All mattresses use the new-style materials that are acid and alkali resistant and disinfection durable.

PARAMETERS:

Table board material	cabon fibre
Whole length	2500±10mm
Length of radiolucent board	1750±10mm
Length without metal inspection	1200±5mm
Whole width	685±5mm
Width of radiolucent board	465mm
Thickness	40mm
X-ray absorption coefficient (aluminium equivalent)	≤0.44mmAL
Max.loading	240kg
Horizontal movement	manual
locking device	electromagnetic locking
longitudinal movement	manual
Move route	250±20mm
locking device	electromagnetic locking
Max.head protract	1880±20mm
Max.foot protract	1050±20mm
Vertical movement	Electric
Min.Bed height from the ground	750±20mm
Max.Bed height from the ground	1100±20mm
Move route	350±20mm
average velocity	8mm/s
Bed surface shifting/Horizontal movement/Locking controller	
average velocity	8mm/s
Table body locking	Electric
Big Four direction wheel	free moving
Inner configurate DC24V	Panasonic lead-acid battery
Single-phase Alternating Current	220VAC 50/60Hz
Input maximum power	300VA
Protective tube	F5AL250V(5*20mm)
Standard accessories:	Memory foam mattress--1pcs
	Arm supports-1pairs

Optional Operation Table Accessories




D0015

Orthopedics Traction Shelf



D001

Orthopedics Traction Shelf



D009

Spine Frame



D0010

Arm Traction Table



D0011

Neurosurgery Headframe



D0012

Neurosurgery Headframe

Obstetric Operation Table

DST-3003

Luxury Electric Obstetric Table



We learn from Japan's advanced technology, the use of human engineering principles, curve design. The utility model is provided with a plurality of auxiliary functions for gynecology, obstetrics, gynecology, diagnosis and examination, including emergency abdomen, and the utility model is especially provided with a unique foot support

PARAMETERS:

Length	1900±20mm
Max Height	890±20mm
Width	680±20mm
Min Height	620±20mm
Up & down lifting range	250±20mm (Electric)
Back Section up	75°up(Electric)
Back Section down	5°down(Electric)
Supply	AC220±10%,50HZ
Seat section	up20° down 0° (electric)

Standard accessories: arm supports, leg supports, foot controller

Obstetric Operation Table

DST-IVA

Electric Obstetric Operation Table

PARAMETERS:

Size	1810 x 600 mm
Height	600-850mm
Reversed /Trendelenburg	≥20°
Back section up&down	≥75°
Auxiliary board	90° (level moving)
Supply	AC220±10%,50HZ
Battery	Optional



This bed is designed base on the concept of family accompanying delivery. It has got 5 national patents including an exterior patent. The driven and control system use imported motors which are safe, reliable, easy operated and with low noise. Both families and nurses can operate it.

- Bed head end and guardrail are made of high strength plastic, which have nice appearance. They can be assembled or detachable easily.
- Bed height is adjustable electrically; the back section upward, trendelenburg and reverse trendelenburg are electric controlled too. Leg section is removable. There is filth basin under the seat, it is easy to adjust the position of lying-in women and help delivery process by combining use of the leg holder, assistant table and foldable footrest.
- The bed is dirt free, easy cleaned. The central brake provides the bed a better stability.

DHC-II(6001)

Electric Delivery Bed



PARAMETERS:

Size	1920x920mm
Height	600-890mmmm
Back section	≥65°up
Reversed Trendelenburg /Trendelenburg	≥5°
Leg board	up 5° down 5°
Power supply	220V 50Hz
Safe loading	200kg

Obstetric Operation Table

DST-3004B

Hydraulic Obstetric Table



PARAMETERS:

Size	1800 x 600 mm
Height	650–895mm
Max angle of back board	up>75° down<6°
Reversed trendelenburg	forward>10° backward>25°
Back board	600*800mm
Seat board	600*400mm
Leg board	540*500mm

PARAMETERS:

Size	1400 x 580 mm
Height	580–980mm
Back board up	65°
Reversed trendelenburg	82°-180°
Back board angle	45°



DST-V

Electric Obstetric Bed



Delivery Bed



A44-4

Stainless Steel Obstetric Bed
Size : 1800 x 600 x 800 mm



A44-5

Epoxy Steel Obstetric Bed
Size : 1800 x 800 x 800 mm



A44-6

Stainless Steel Exam Bed
Size : 1800 x 600 x 800 mm



D4

Stainless Steel Exam Bed
Size : 1800 x 600 x 800 mm



D4-3

Stainless Steel Exam Bed
Size : 1950 x 650 x 680 mm



A44-1

Stainless Steel Delivery Bed
Size : 1800 x 600 x 800 mm



A44-3

Epoxy Steel Delivery Bed
Size : 1800 x 800 x 800 mm



A44-2

Stainless Steel Delivery Bed
Size : 1800 x 600 x 800 mm



A44

Stainless Steel Delivery Bed
Size : 1870 x 600 x 765 mm

Obstetric Operation Table



DST-IV

Luxury Electric Obstetric Table

1. The main body is made of 45 # steel, The appearance parts are made of epoxy steel and 304 stainless steel, The base cover is high quality ABS;
2. Using imported motors with low noise, stable performance, Precise positioning, strong thrust,
3. The unique three - layer structure has achieved thermal insulation and cushioning. The minimum height 550mm platform, The platform is divided into backplane, seat plate and footboard. Through the operation handle, the upper and lower sides of the back plate can be changed, the bed will be tilted back and forth, so that the doctor can get an ideal surgical location and make the patient get comfortable seats.
4. The newborn baby's chassis can be rotated, which can correspond to the side care. The whole bed is running smoothly without noise and vibration.

PARAMETERS:

Size	1080*2058*550mm
Min height and Max height	550mm / 910mm
Translation distance of rotary tray support plate	60mm
Foot brace folding angle	0-90° gas spring
Foot brace rotation angle	Left 33° & right 33° (Electric)
Back Section up & down	60°up & 10°down (Electric)
Seat Section up & down	25°up & 25°down (Electric)
Power supply	AC220,2.5A 50HZ
Foot board folding angle	75° (manual)
Side handle rotation angle	180°(manual)
Tray support plate moving distance	150mm
Tray support plate rotation angle	90°

Medical Ceiling Pendant Series



YD-580

Medical Ceiling Column




Bedhead Unit

We customize the bedhead unit according to your requirements

Nurse Call System

It is used in patient room & ICU room, we customize the nurse call system, you just tell us your hospital's bed details

The medical column is used for ICU, NICU room, the configuration is choiced by your requests.



Ice Ceiling Bridge

YD-180E

Luxury ICU Ceiling Bridge

PARAMETERS:

Length	2000-3800mm
Payload gliders	100kg(220p)
Payload shelf	60kg (110p)
Ambient light	2*54w
Reading light	220v, 2*15w LED spotlights
RJ45	1pcs
Column length	700mm
Electrical socket(220V)	8pcs
Gas outlet(Germany)	O2x2,VACx1,AIRx1
Shelf x1,Drawer x1, basket x1, extension arm x1, infusion rack x1	

Optional item:

1. 1-4pcs shelves (25x10/30x10mm) 2. Drawer layer:1/2

Medical Gas Outlets Standards: optional British, German, French, Metric, Ohmeda, DISS, Chemetron, Purit-Bennett Electrical Sockets:customized as clients' request



☒ ICU Bridge Type Ceiling-Mounted system, specially designed for intensive care rooms, the entire process wards, the recovery room. Supply flexible, first-rank nursing satisfy patient's separate treatment and nursing demand. Unique design lies in it all hides the electrical cable and pipeline in the ceiling mounted arm.


☒ The modular design and multiple configurations.

☒ As clients' requirements, providing reasonable ICU Combination Ceiling-Mounted pendant (Apart Dry-Wet).

☒ Advanced spray technology, the intelligent combination of innovative ways to make the bridge pendant with diverse functions.

YD-180E/X

ICU Ceiling Bridge (Dry And Wet Separated)



PARAMETERS:

Length	2000-3800mm
Payload gliders	100kg(220p)
Payload shelf	60kg (110p)
Ambient light	2*54w
Reading light	220v, 2*15w LED spotlights
RJ45	1pcs
Column length	700mm
Electrical socket(220V)	8pcs
Gas outlet(Germany)	O2x2,VACx1,AIRx1
Shelf x1,Drawer x1, basket x1, extension arm x1, infusion rack x1	

Optional item:

1. 1-4pcs shelves (25x10/30x10mm)

2. Drawer layer:1/2

Medical Gas Outlets Standards: optional British, German, French, Metric, Ohmeda, DISS, Chemetron, Purit-Bennett

Medical Ceiling Pendant



YD-3600-2

Double Arm Ceiling Pendant

PARAMETERS:

Double arm/Single arm	800mm+800mm/800mm
Column length	1000mm
Max.loading capacity	350kg
Arm rotating angle	340°
Braking system	Optional Merchanical/Pneumatic braking
Electrical power	220V,50Hz
System	Electric
RJ45	1 pcs
Shelf	3 shelves
Electrical socket	6 pcs
Gas outlet(Germany standard)	O2x1,VACx1,AIRx1
Drawer x1,basket x1,extension arm x1,infusion rack x1,IV pole rack x1	
Optional item:	
1-4pcs shelves (25x10/30x10mm)	
Medical Gas Outlets Standards: optional	

British, German, French, Metric, Ohmeda, DISS, Chemetron, Purit-Bennett , Electrical Sockets:customized as clients' requests.
Optional anesthesia,endoscopy and surgery pendant

- Pendant Applicative Area: Surgical Room/ Emergency Room/ICU/Recovery room
- Ceiling Mounted installed.Within 340°revolved,Adjustterminal box location,convenience for medical staff operation.
- The suspension arm is imported by high-strength aluminum alloy extrusion, surface designed by the first level oxidized rust prevention.
- Friction brake + Pneumatic brake,adopted by German Standard Gas terminal(Insert and Pull more than 20,000 times.)
- Provide multiple configure,make you satisfy.
- Electrical up and down,precise positioning.



YD-3200

Single Arm Ceiling Pendant

Medical Ceiling Pendant

YD-3600-2D

Double Arm Electric Ceiling Pendant

PARAMETERS:

Double arm/Single arm	800mm+800mm/800mm
Column length	1000mm
Max.loading capacity	350kg
Arm rotating angle	340°
Braking system	Optional Merchanical/Pneumatic braking
Electrical power	220V,50Hz
System	Electric
RJ45	1 pcs
Shelf	3 shelves
Electrical socket	6 pcs
Gas outlet(Germany standard)	O2x1,VACx1,AIRx1
Drawer x1,basket x1,extension arm x1,infusion rack x1,IV pole rack x1	

Optional item:

1-4pcs shelves (25x10/30x10mm)

Medical Gas Outlets Standards: optional

British, German, French, Metric, Ohmeda, DISS, Chemetron, Purit-Bennett , Electrical Sockets:customized as clients' requests.

Optional anesthesia,endoscopy and surgery pendant



Pendant Applicative Area: Surgical Room/ Emergency Room/ICU/Recovery Room

Ceiling Mounted installed.Within 340°revolved,Adjustterminal box location,convenience for medical staff operation.



The suspension arm is imported by high-strength aluminum alloy extrusion, surface designed by the first level oxidized rust prevention.

Friction brake + Pneumatic brake,adopted by German Standard Gas terminal(Insert and Pull more than 20,000 times.)

Provide multiple configure,make you satisfy.

Electrical up and down,precise positioning.

YD-3200D

Single Arm Electric Ceiling Pendant

Page : 28

Patient Monitor



PM-90

Multi-Parameter Patient Monitor

- 12.1 inches high resolution color TFT screen
- User friendly and configurable interface design
- Intuitive human voice alarms
- Powerful data review:
 - 480 hours graphic and tabular trends of all parameters
 - 1000 pieces of NIBP record storage
 - 700 pieces of alarm review
- Built-in rechargeable Li-ion battery for up to 5 hours of continuous work
- VGA output for external display(optional)
- Wire/wireless central monitoring system solution(optional)
- Standard Parameters:3/5-lead ECG,PR,RESP,NIBP,SpO2,TEMP
- Optional: 2-IBP,EtCO2(mainstream,sidestream),central monitoring system

● **Respironics EtCO₂ (optional)**



CAPNOSTATS
mainstream module



LoFlo sidestream module



EtCO₂ (Respironics)
"Plug and Play"



PM-15

Multi-Parameter Patient Monitor

- 15 inches high resolution color TFT display
- Optional full touch-screen for easy operation
- Attractive and durable silicon buttons with backlight
- ST segment analysis and 23 types of arrhythmia analysis
- Optional simultaneous 12-lead ECG analysis for further diagnostic confirmation
- Maximum 480-hour graphic and tabular trends,1000NIBP records, 700 alarm events and 2-hour ECG waveform storage and review
- ESU proof and defibrillation synchronization for different working environments
- Dual rechargeable Li-ion batteries for up to 10 hours of continuous work
- Dual USB ports enable software upgrade and also supports external USB storage
- Customer-defined color configurations for all parameters display setting (32 colors)
- Nurse call function and bi-directional communication solution
- VAG output for external display(optional)
- Flexible mounting and trolley solution(optional)

● **Masimo SpO₂**



Masimo M-LNCS™
SpO₂ extended cable



M-LNCS™ Y1 adult, pediatric,
neonatal reusable sensor



M-LNCS™ DC1 adult
reusable finger clip sensor



M-LNCS™ DCIP pediatric
reusable finger clip sensor

Doppler Ultrasound Equipment

- 1. 15-inch high resolution LED monitor
- 2. standard PC LED backlit silicon gel keyboard
- 3. one key regulation, easy operation
- 4. rechargeable built-in battery
- 5. abundant peripherals: DICOM3.0, VGA, video out, USB etc

Detailed Product Description

- A. Imaging optimization technology
- B. Compound enhance technology
- C. Imaging processing technology
- D. Imaging optimization technology
- E. Compound enhance technology
- F. Multi-beam parallel processing technology
- G. Tissue harmonic imaging

Displaying mode

- B, 2B, 4B, B/M, M
- B/C, B/C/D, B/D CFM, PRF
- PW, velocity, power (direction), histogram
- Triplex / Duplex

Probe

- Electronic convex, linear, transvaginal
- Wide band multi-frequency



C-60

Portable Color Doppler Ultrasound Machine



- 1. 15-inch high resolution LED monitor
- 2. standard PC LED backlit silicon gel keyboard
- 3. one key regulation, easy operation
- 4. rechargeable built-in battery
- 5. abundant peripherals: DICOM3.0, VGA, video out, USB etc

C-60PLUS

High-class Portable Color Doppler Ultrasound Machine



Flexible and convenient operating system

A type of high array, high definition, multi-function full-digital laptop Color Doppler ultrasonic diagnostic apparatus. Light and portable, more clearer, delicate and smooth image, convenient to operate, and powerful electrical power ability supports the machine for consultation in various environment of cities, towns, outdoor. A variety of charging to ensure the examination consultation under different environments. Supernormal battery life 2-3 hours. 15 inch LED shows the images with more flexible ways.

Operating Mode:

Black and white image : B, 2B, 4B, left B | M, B | D, PW, M, B partial enlarged mode
 Bloodstream image : Spectrum : B / D, B / C / D CFM: B | C | D, B | C | M, B | C pairs in real time , PW, CFM, CPA
 Color flow image adjustment parameters : a Doppler frequency , sampling frame position and size , baseline color gain, the deflection angle , the wall filter, the cumulative frequency.
 Display depth: Continuously adjustable



Doppler Ultrasound Equipment

VET-6

Veterinary Ultrasound Scanner



1. Perfect Ergonomic Design

- Compact, sophisticated, portable
- 5.6 inch LED display
- Lightweight body
- Rechargeable built-in battery
- User-friendly operation interface
- Simple and beautiful operation panel
- Fast, convenient and flexible touch-type keyboard operation

2. Advanced Technology and Function

- Advanced digital signal processing technology
- Leading digital beam former technology
- NC dynamic frequency scanning technology
- Real-time point by point dynamic focusing reception technology
- Real-time dynamic aperture imaging

DESCRIPTIONS :

The C-80 color doppler ultrasound system bringing an entirely new level of diagnostic performance and workflow efficiency. The C-80 system, combines best-in-class image quality and a robust set of features to meet daily clinical needs. It enhances diagnostic confidence with high quality color and power Doppler, 2D-mode imaging. All-digital signal processing and multibeam formation technology provides best-in-class image quality in 2D-mode, Tissue Harmonic Imaging, Dynamic Tissue Optimization Technology and Doppler modes.

Data and information communication can be easily carried out. It also can connect to PACS (Picture Archiving and Communication System) as a terminal.

PACS is a medical picture archiving and communication systems, which greatly facilitates network management system of the hospital. It is convenient for remote diagnosis. The machine is powerful, convenient and easy to operate, supports B, B/B, B/4B, B/M, B/PWD, CFM, CDE, B/CFM/D. Meanwhile, the display can be moved from up and down, forward and backward, turn to left and right. These powerful features have greatly boosted diagnostic confidence.

USER INTERFACE

- Intuitive Windows-based operating principles
- User-centric control panel with HomeBase layout and control customization
- On/Off task light and back-lit illumination of control panel
- Variable brightness indicates active state of function keys
- Easily accessible, full size QWERTY keyboard for text entry, function keys and system programming

- Cine Playback, Multiple Arrows, Configurable Worksheets, Exam Review, Pictograms (Body Marks), System Setup Menu

C-80

Full Digital 3D Mobile Ultrasound Machine Color Doppler



Hospital Bed Series



- 
Mattress holders
- 
Leg plate at both ends of the regulation, adjust the comfort of the leg board
- 
Extend bedboard
- 
TENTE silent castors central control braking
- 
Manual CPR
- 
X-ray cassette
- 
anticollision wheels with IV pole



DHC-III (FM04)

Luxury Multi-functional Electric ICU Bed With Weighing Function

PARAMETERS:

- Size: L2170xW1080xH500/850mm
- Bed plate: L1950xW910mm
- The whole bed can weigh, with error controlled in 1-3‰
- Back plate adjusting 70°±5°, leg plate adjusting 45°±5°
- Head and foot tilt 16°±3°

Optional weighing function



Model: FM03
4 functions



Model: FM04
5 functions



Hospital Bed Series

DHC-IV(FP07)

Luxury Electric ICU Bed With 7 Functions



PARAMETERS:

Size	2240 x 1110 x 630/1030 mm	Overall lifting	630 -1030mm
Back adjustment	40°±2°	Lateral tilt	12°±2°
Leg adjustment	13°±2°	Luxury central control castors	
Trendelengurg angle	17°±2°	ReverseTrendelengurg angle	25°±2°
4 pieces of European syle guaedrail, upside down structure			
Foot controller, nurse controller, patient controller			
Optional mattress, bedside cabinet, overbed food table			

Five motors



Weighing operator



Foot control pedal



Nurse control panel



Panel on guardrail



Bedboard soft connection



ABS guardrail



TENTE Castors



Bed panel with punching design

Hospital Bed Series

DHC-IV(FP05)

Electric ICU Hospital Bed With 5 Functions



PARAMETERS:

Size	2300 x 1050 x 530mm/710mm	Bedboard	1955 x 900 mm
Net weight	158kg	Loading capacity	250kg
Back board	0-65°	Foot board	0-43°
Trendelenburg/Reversed Trendelenburg			0-10°
"Reset " functions , be can be back to "zero" position quickly. (manual CPR)			
Backup battery for emergency power failure			
Bed board has safe connection which make the bending of the bed board.			
Detachable head&food board in PP plastic,can be fixed when assembled.			
4 pcs foldable PP siderails			
Central locking castors controlled by electric.			
IV.pole mount bracked at 4 corners of the bed frame			
With easily operated touching switch,solve the problem of limited service life,mechanical keypad and keypad invalid due to constantly operating			
125mm silent castors			
Optional mattress, bedside cabinet,overbed food table, iv stand pole etc.			

DHC-II Series Hospital Bed



FO03 Electric Hospital Bed

PARAMETERS:

Bed size	2280 x 1120 x 420/780mm
Height lifting	>200mm
Backboard up	>65°
Leg board up	>35°

Asymmetric European style guardrail ,Turn-over structure
 Luxury central control castors,brake the bed with one stamp.
 Bed panel is made of one-step molding cold rolled steel,
 with excellent performance in anti-slipping and air permeability.



FO - 2 Manual Hospital Bed

PARAMETERS:

Bed size	2280 x 1120 x 420/780mm
Height lifting	>200mm
Backboard up	>75°
Leg board up	>37°

Asymmetric European style guardrail ,Turn-over structure
 Luxury central control castors,brake the bed with one stamp.
 Bed panel is made of one-step molding cold rolled steel,
 with excellent performance in anti-slipping and air permeability.



FO - 1 Manual Hospital Bed

PARAMETERS:

Bed size	2280 x 1120 x 500mm
Backboard up	>75°
Leg board up	>37°

Asymmetric European style guardrail ,Turn-over structure
 Luxury central control castors,brake the bed with one stamp.
 Bed panel is made of one-step molding cold rolled steel,
 with excellent performance in anti-slipping and air permeability.



DHC-III Series Hospital Bed

FP04 Electric Hospital Bed



Bed size	2280 x 1020 x 420/780mm
Height lifting	≥200mm
Backboard up	≥65°
Leg board up	≥35°
Luxury central control castors,brake the bed with one stamp	

FP03 Electric Hospital Bed



Bed size	2280 x 1020 x 420/780mm
Height lifting	≥200mm
Backboard up	≥65°
Leg board up	≥35°
Luxury central control castors,brake the bed with one stamp	

FP - 2 Manual Hospital Bed



Bed size	2280 x 1020 x 420/780mm
Height lifting	≥200mm
Backboard up	≥75°
Leg board up	≥37°
Luxury central control castors,brake the bed with one stamp	

FP - 3 Manual Hospital Bed



Bed size	2280 x 1020 x 500mm
Backboard up	≥75°
Leg board up	≥37°
Luxury central control castors,brake the bed with one stamp	
Asymmetric European style guardrail ,Turn-over structure	

Hospital Bed Series

DHC-II(FE01)

Electric Medical Bed With 5 Functions



PARAMETERS:



Size	2175 x 1090 x 420/720mm
Back board tilt	65° electric / manual
Leg board tilt	35° electric / manual
Central controlling	Silent castors
Optional mattress, bedside cabinet, overbed food table, IV Stand	Optional 5,3,2 functions or manual function

DHC-II(FE01)

Electric Medical Bed With 3 Functions



FN-2

Manual Hospital Bed With 3 Functions



PARAMETERS:



Size	2200 x 1020 x 420/760mm or height 500mm for two function
Back board tilt	75°
Leg board	37°
Safe loading	250kg
Optional mattress, bedside cabinet, overbed food table, IV Stand	

FN-1

Manual Hospital Bed With 2 Functions



Hospital Bed Series

DHC-II (FH05)

Economical Electric Home Care Bed



DHC-II (FH06)

Luxury Electric Home Care Bed



PARAMETERS:



Size	2150x1030x460mm	Height(from floor)	460-860mm
Back plate adjusting	0-70°	Knee section adjusting	0-40°
Trendelenburg/Reversal trendelenburg		13°/13°	
Auto contour position		CPR functions	
Safety load capacity		250kgs	
Optional mattress, bedside cabinet, traction shelf			

DHC-II (FH07)

Home Nursing Bed



Bed size : 2100 x 960 x 540mm (including mattress)
 The hole size : 240*200mm
 The overbed table : 1030 x 300mm
 The castors : 125mm braking

DHC-II (FH08)

Home Nursing Bed



Bed size : 2100 x 960 x 540mm (including mattress)
 The hole size : 240*200mm
 The overbed table : 1030 x 300mm
 The castors : 125mm braking

Hospital Bed Series

FA-5

Electric Bed With 3 Functions



Size : 2200 x 960 x 420/720mm

FA-6

Manual Hospital Bed With 3 Functions



Size : 2200 x 960 x 440/790mm

FA-7

Double Cranks Hospital Bed



Size : 2200 x 960 x 520mm

FA-8

Double Cranks Hospital Bed



Size : 2150 x 950 x 500 mm

A25

Double Cranks Hospital Bed



Size : 2150 x 900 x 500mm

A23

Single Crank Hospital Bed



Size : 2150 x 900 x 500 mm

Hospital Bed Series

A24

Hospital Bed



Size : 2150 x 900 x 500mm

A26

Single Crank Hospital Bed



Size : 2150 x 900 x 500 mm

A27

Hospital Bed



Size : 2150 x 900 x 500mm

A31

Double Crank Hospital Bed
Stainless steel



Size : 2040 x 900 x 500 mm

A34

Double Cranks Hospital Bed
Steel coating



Size : 2040 x 900 x 500mm

A33

Hospital Bed Stainless steel



Size : 2040 x 900 x 500 mm

Hospital Bed Series

A36

Hospital Bed



Size : 2040 x 900 x 500 mm

A37-1

Orthopedic Traction Bed



Size : 2020 x 900 x 500 mm

A39D

Child Bed With 1 Function



Size : 1600 x 850 x 610 mm

A39

Double Cranks Child Bed



Size : 1770 x 900 x 550 mm

A39-1P

Luxury Child Bed With 2 Functions



Size : 1770 x 900 x 550 mm

CR



Size : 200/110cm

Medical Stretcher Series

B-2H

Luxury Electric Transfer Stretcher



With unique motion design, this cart could transport patient comfortably between hospital bed and transfer cart steadily without need of moving patient

Surface size: 2020mm x 680mm

Size when connected: 3900mm x 680mm

lifting range : 590mm - 880mm

move range of shifting plate : 510mm

move speed of shifting plate : 30mm/s

Max loading: : 135kg

Accessories :

IV pole 1piece hand controller 1set

power cord 1piece connecting cord 1piece

B-2

Medical Connecting Transfer Stretcher



Product composition:

This product has A car, B and C three major components, including car A and car B, each has A screw and four control of castor, senior C car there are two big barrier and ABS engineering plastics of bed plate.

PE docking cart A

Size: 1790 x 450 mm , Height 570 - 850mm

Movements distance: 280 ± 20 mm (manual)

PE docking cart B

Size: 1835 x 450 mm , Height 570 - 850mm

Movements distance: 280 ± 20 mm (manual)

PE docking cart C

Size: 1955 x 585 mm

Width (guardrail up): 780 ± 20 mm , Width (guardrail down): 710 ± 20 mm

Height : 220 ± 20 mm

Back adjustment 90° ± 5° (manual)

PE docking cart (ABC)

Size: 3575 x 770 mm

Medical Stretcher Series

B-1Y

Hydraulic Emergency Stretcher



PARAMETERS:

This Hydraulic Stretcher has Central-lock castors,diameter 150mm;bed board for X-ray Examination;

Hi-low range:800mm-550mm;

Backrest angle:0-60°

Knee angle:0-40°

Trendelenburg & Rev.Trendelenburg:0-12°

B-1

Manual Emergency Stretcher



PARAMETERS:

Size:L1930*W750*H600/830mm

Specification:

- 1.main frame of the flat cart is jet-molding carbon steel,solid and good looking.
- 2.cart surface and guardrail is structured by one-step molding ABS material,solid and elegant.
- 3.adopting separate design: back plate is adjusted by air spring,could be operated by only one hand,convenient to use.
- 4.the crank of the cart(the screw is equipped with a two-way protection device) can adjust height of the cart surface from 600 to 830mm
- 5.equipped with international advanced central -control braking system,stable and reliable,the oriental castors are easy to operate.
- 6.meeting requirement from overall design concept,more humanized,ideal for transporting patients in hospital

Medical Stretcher Series

B-1A

Manual Transfer Stretcher



PARAMETERS:

Hospital transfer stretcher trolley

Features:

1. It is a type of manual transfer stretcher. Turn by deasil, the stretcher will be up, and turn by widdershins, the stretcher will be down;
2. Aluminum-Alloy protective rails, can be folded and hidden under the surface on no use;
3. Central-lock castors, diameter 150mm. Matched with direction castor, help the operator controlling the direction easily;
4. Rotate the crand and handle on bedboard, patient on the mattress can be moved parallelly to another bed or flat table;
5. Hi-Low range: 920mm-610mm
6. Weight: 89 kgs

B4-1

Emergency Trolley



It is made of high-strength aluminum alloy thick-wall pipes. The bed surface is made of foamed cushion and the adjustable back making the wounded fell comfortable. This machine is very convenient. One person is enough to send the patient (with the stretcher) to the ambulance car. The mechanical foldaway legs are controlled by the handles on both ends of the stretcher. And the legs can be folded under the stretcher automatically when loading into the vehicle. Also the stretcher can be locked by the fixing devices when on the ambulance car. This type of stretcher is characterized by its light-weighted, durable, safe usage and easy for sterilization. The stretcher adopts stretcher trolley and stretcher to separate structure.

high position: 1900*550*920mm
low position: 1900*550*260mm
maximum angle of the back: 75°
packing size (1pc): 1980*620*280mm
loading bearing: <= 159kg
N.W: 43kg

Medical Trolley Series

B4
Stainless Steel Emergency Trolley



Size : 1950 x 600 x 750 mm

B4-2
Portable Emergency Stretcher



Size : 2080 x 550 x 130 mm

MWT-1
Mobile Working Station



Size : 730 x 470 x 910 mm

MWT-2
Mobile Working Trolley



Size : 1050 x 635 x 1110 mm

CQ-08
Emergency Trolley



Size : 810 x 520 x 960 mm

CQ-09
Emergency Trolley



Size : 615 x 550 x 1010 mm

CQ-06
Emergency Trolley



Size : 762 x 520 x 1045 mm

CQ-10
Emergency Trolley



Size : 630 x 470 x 910 mm

Medical Trolley Series

CM-01
Anesthesia Vehicle



Size : 615 x 550 x 1010 mm

CM-02
Anesthesia Trolley



Size : 630 x 470 x 1500 mm

CM-03
Anesthesia Trolley



Size : 810 x 520 x 960 mm

CM-07
Anesthesia Vehicle



Size : 762 x 520 x 1045 mm

CZ-08
Treatment Trolley



Size : 775 x 635 x 1080 mm

CZ-04
Treatment Trolley



Size : 785 x 480 x 995 mm

CZ-05
Treatment Trolley



Size : 762 x 520 x 955 mm

CF-05
Medicine Trolley



Size : 775 x 635 x 1080 mm

CF-02
Medicine Trolley



Size : 785 x 480 x 995 mm

Medical Trolley Series

CF-06

Medicine Trolley



Size : 630 x 490 x 1070 mm

CF-05

Medicine Trolley



Size : 750 x 450 x 1220/1275 mm

B26

Medicine Trolley Stainless Steel



Size : 570 x 400 x 785 mm

B19

Stainless Steel Nursing Trolley



Size : 700 x 420 x 930 mm

CS-05

Infusion Trolley



Size : 630 x 470 x 910 mm

CS-04

Infusion Trolley



Size : 630 x 470 x 910 mm

CB-02 (25)

Medical Files Cart



Size : 510 x 390 x 780 mm

CB-07 (50)

Medical Files Cart



Size : 820 x 390 x 880 mm

CH-02

Luxury Nursing Cart



Size : 1185 x 570 x 950 mm

Medical Trolley Series

CH-07
Nursing Cart



Size : 800 x 590 x 890 mm

CH-06
NursingTrolley



Size : 570 x 430 x 890 mm

B34
Stainless Steel Nursing Trolley



Size : 960 x 550 x 920 mm

B32
Stainless Steel Nursing Trolley



Size : 800 x 600 x 920 mm

B7-2
Luxury Baby Cot



Weighing, Music play and ambient temperature display

B10
Stainless Steel Baby Trolley



Size : 840 x 530 x 840/960 mm

B9
Stainless Steel Baby Cot



Size : 920 x 490 x 1020 mm

B15
Instrument Trolley
Stainless Steel



Size : 750 x 450 x 910 mm

B18
Instrument Trolley
Stainless Steel



Size : 750 x 450 x 910 mm

Medical Trolley Series

B13

Instrument Trolley



Size : 750 x 450 x 850 mm

B14

Instrument Trolley Stainless Steel



Size : 1400 x 450 x 910 mm

B16

Treatment Trolley Stainless Steel



Size : 750 x 450 x 910 mm

CZ- 05

Nursing Trolley



Size : 610 x 410 x 830 mm

CZ- 06

Luxury Instrument Trolley



Size : 620 x 420 x 860/930mm

CY- 04

Instrument Trolley



Size : 630 x 470 x 910 mm

CY- 01

Instrument Cart



Size : 730x 475 x 1415 mm

B30

Stainless Steel Sending Goods Cart



Size : 1100 x 520 x 930 mm

B35

Stainless Steel Mayo Table



Size : 590 x 450 x 940/1200 mm

Medical Cabinet Series

B36
Stainless Steel Trolley



Size : 690 x 450 x 940/1200 mm

B40
Power Supply Instrument Vehicle



Size : 650 x 550 x 890 mm

B70
Instrument Trolley



Size : 945 x 450 x 810 mm

B37
Oxygen Bottle Trolley



Size : 810 x 530 x 810 mm

C12
Stainless Steel Working Stable



Size : 1800 x 700 x 810/1750 mm

C32
Medicine Cabinet With Stainless Steel Foundation And Table Surface



Size : 900 x 240/400 x 1750 mm

C35
Medicine Cabinet



Size : 900 x 240/530/400 x 1750 mm

C36
Stainless Steel Instrument Cabinet



Size : 900 x 400 x 1750 mm

C38-1
Stainless Steel Instrument Cabinet



Size : 900 x 400 x 1750 mm

Medical Cabinet Series

C39
Instrument Cabinet



Size : 900 x 400 x 1750 mm

C47
Steel Disassembly Type Six Door Locker



Size : 960 x 400 x 1750 mm

C42
Steel Disassembly Type Four Door Locker



Size : 900 x 400 x 1750 mm

C49
Steel Disassembly Type Three Door Locker



Size : 960 x 400 x 1750 mm

C48
Steel Detachable Two Door Locker



Size : 900 x 400 x 1750 mm

C50
Shoes Cabinet



Size : 960 x 400 x 1750 mm

C58
Stainless Steel Cabinet



Size : 900 x 250/400 x 1750 mm

C37
Steel Dismounting Type Instrument Cabinet



Size : 900 x 400 x 1750 mm

C38 CR
Stainless Steel Dismounting Type Instrument Cabinet



Size : 900 x 400 x 1750 mm

Medical Cabinet Series

C38
Steel Dismounting Type
Instrument Cabinet



Size : 900 x 400 x 1750 mm

C40
Steel Dismounting Type
Instrument Cabinet



Size : 900 x 400 x 1750 mm

C42
Sterile Cabinet



Size : 900 x 400 x 1750 mm

C53-1
Stainless Steel Embedded
Instrument Cabinet



Size : 810 x 530 x 810 mm

C7
ABS Bedside Cabinet



Size : 480 x 480 x 750 mm

C8-1
ABS Bedside Cabinet



Size : 480 x 480 x 880 mm

CC-04
Luxury Bedside Cabinet



Size : 520 x 400 x 725 mm

CC-03
Bedside Cabinet



Size : 480 x 470 x 795 mm

CC-05
Luxury Bedside Cabinet



Size : 520 x 470 x 700 mm

Medical Chair Series

C13 - 1
Stainless Steel Basin
Sink With Sensor



Size : 1600 x 620 x 1900 mm

D11
Stainless Steel Accompany
Chair



Size : 700 x 650 x 880 mm

D10
Steel Coating Accompany
Chair



Size : 700 x 650 x 880 mm

D12 - 1
Luxury Accompany Chair



Size : 780 x 650 x 1060 mm

D5
Stainless Steel Infusion Chair



Size : 700 x 1120 x 1100 mm

D6
Steel Coating Infusion
Chair



Size : 700 x 1120 x 1100 mm

D6 - 1
Transfusion Chair



Size : 700 x 850 x 1100 mm

D6 - 2
Luxury Transfusion Chair



Size : 770 x 680 x 950 mm

D6 - 3
Luxury Transfusion Chair



Size : 710 x 1130 x 1100 mm

X-Ray Viewing Box Series

D1
New X-ray Film Viewer



Size : 440 x 585 x 110mm

D2
New X-ray Film Viewer



Size : 900 x 585 x 110mm

D3
New X-ray Film Viewer



Size : 1210 x 585 x 110mm

LED - T1
LED X-ray Film View Box



Size : 500 x 500 x 48mm

LED - T2
LED X-ray Film View Box



Size : 1230 x 500 x 48mm

D19
Doctors Stool



Size : Dia340 x 750/920 mm

D25
Doctors Mobile Stool



Size : Dia320 x 460/620 mm

D26
Medical Stool



Size : Dia320 x 430/590 mm

D13
Stainless Steel Folded Screen



Size : 1910 x 1800 mm

Medical Scale Series

D21
Waiting Chair



Size : 1800 x 680 x 800 mm

D22 (A)
Steel Coating Footstep



Size : 620 x 380 x 420 mm

D22 (B)
Stainless Steel Footstep



Size : 540 x 400 x 400 mm

S - 2
Adult Scale



160KG

S - 1
Adult Scale



160KG

S - 4
Baby Scale



20KG

S - 3
Baby Scale



20KG

C55-1
Luxury Overbed Food Table



Size : 850 x 450 x 820/1030 mm

C55
Luxury Overbed Table



Size : 555 x 450 x 940/1200 mm

Medical Mattress Series

C58
Folding Overbed Table



C57
Overbed Table



D18
Mattress With PU Leather
And High Density Foam



Size : 1950 x 800 x 100 mm

D15
Medical Air Mattress



Size : 2000 x 900 x 70mm

D17
Mattress With Full Foam



Size: 1950 x 800 x 80(mm)

D16
Mattress With Half Of Palm And
Foam



Size : 1950 x 800 x 80(mm)

D8
Epoxy Steel Infusion Chair



Size : 1000 x 600 x 400 mm

D7
Blood Donor Chair
Optional Stainless Steel
Or Steel Coating



Size : 1820 x 520 x 560 mm

B71
Clinical Examination Bed
Optional Stainless Steel
Or Steel Coating



Size : Dia360 x 2 x 760 /1260mm



HOSPITAL FURNITURE

www.hamptonsurgical.co.uk



Glass Door Filing Cabinet



YOUYE Group file cabinet products, a large to a shape, phosphate anti-corrosion treatment, automatic spray technology, reasonable design, fashionable, high quality, great to meet the needs of users.

1, available on demand free combination; 2, new hinges, door open light; simple number registered, convenient access to data; quality locks open and flexible, each opening rate is extremely low; 5, blue point deduction hand, novel personality, ABS is not easily broken; double slide, drawer sports facilities, and pull out the 90%; 7, the partition can be freely adjusted up or down or decrease up or down or decrease the full use of cabinet space; 8, skim, phosphating twelve stations before processing to ensure that the coating is not easy to fall off, and weather resistance.

Glass Door Filing Cabinet



TG-008-1
Full Glass Open Door Cabinet
900*400*1800mm



TG-013
Six Drawers Half Glass Cabinet
900*400*1800mm



TG-023
Short Glass Cabinet
900*400*900mm



TG-023-1
Short Glass Cabinet
1180*400*900mm

Glass Door Filing Cabinet



TG-008
Full Glass Sliding Door Cabinet
900*400*1800mm



TG-009-1
Middle Drawer Glass Filing Cabinet
900*400*1800mm



TG-063-1
Half Glass Cabinet
900*400*1800mm



TG-063
Half Glass Cabinet
900*400*1800mm



TG-066
Middle Drawer Glass Filing Cabinet
900*400*1800mm



TG-090
Middle Drawer Glass Filing Cabinet
1180*400*1800mm

Steel Made Filing Cabinet



YOUYE Group file cabinet products,a large to a shape,phosphate anti-corrosion treatment,automatic spray technology,reasonable design,fashionable,high quality,great to meet the needs of users.
1,available on demand free combination;2,new hinges,door open light:simple number registered,convenient access to data;quality locks open and flexible, each opening rate is extremely low; 5,blue point deduction hand, novel personality. ABS is not easily broken: double slide,drawer sports facilities,and pull out the90%.;7,the partition can be freely adjusted up or down or decrease up or down or decrease the full use of cabinet space;8,skim,phosphating twelve stations before processing to ensure that the coating is not easy to fall off,and weather resistance.

Steel Made Filing Cabinet



TG-065
Four Doors Cabinet
900*400*1800mm



TG-011
Double Doors Inside
Two Drawers Cabinet
900*400*1800mm



TG-062
Double Doors Cabinet
900*400*1800mm



TG-038
Four Floors Archive Cabinet
900*400*1800mm



TG-039
Five Floors Archive Cabinet
900*400*1800mm



TG-065-1
Double Doors Short Cabinet
900*400*1800mm

Rolling Door Cabinet



FASION DESIGN

AND EXCELLENT QUALITY

Bring You Comfortable and Gladdergladdest Mood



Rolling Door Cabinet



01 ROLLING DOOR CABINET
1000*420*1980mm



03 ROLLING DOOR CABINET
1000*420*1020mm



05 ROLLING DOOR CABINET
1000*500*650mm



02 ROLLING DOOR CABINET
1000*420*1980mm



04 ROLLING DOOR CABINET
1000*420*1500mm



06 ROLLING DOOR CABINET
1000*420*1980mm

Combination Lock



GB-01
900*400*1800mm


Humanized desing is not a slogan Every design comes from experience of its future user.



Combination Lock



GB-01
W900*D400*H1800



GB-02
W900*D400*H1800




GB-04
W900*D400*H900



Combination Lock



GB-05
W900*D400*H1800



GB-06
W900*D400*H1800



GB-07
W900*D400*H900



Mobile Cabinet



FASION DESIGN

AND EXCELLENT QUALITY

Bring You Comfortable and Gladdergraddest Mood



Mobile Cabinet



F4 -01
01 QUICKLY ARK
900*450*(1310-1010-700)mm



F4 -03
01 QUICKLY ARK
450*620*(1310-1010-700)mm



A4 -05
05 QUICKLY ARK
410*540*655mm



F4 -02
02 QUICKLY ARK
450*620*(1310-1010-700)mm



F4 -04
04 QUICKLY ARK
450*620*(1310-1010-700)mm



A4 -06
06 QUICKLY ARK
410*540*655mm


Door Locker



FASION DESIGN

AND EXCELLENT QUALITY

Bring You Comfortable and Gladdergladdest Mood



Door Locker



CG-03
Three Doors Locker
400*400*1800mm



CG-04
Four Doors Locker
400*400*1800mm



CG-05
Five Doors Locker
400*400*1800mm



TG-015
Nine Doors Locker
980*400*1800mm



TG-051
Twelve Doors Locker
980*400*1800mm



TG-049
Fifteen Doors Locker
980*400*1800mm



TG-052
Eighteen Doors Locker
980*400*1800mm



TG-053
Twenty Four Doors Locker
980*400*1800mm



TG-019
Thirty Two Doors Locker
980*400*1800mm

Doors Wardrobe



FASION DESIGN

AND EXCELLENT QUALITY

Bring You Comfortable and Gladdergraddest Mood



Doors Wardrobe



YG-02 上下两门衣柜
Double Doors Wardrobe
400*480*1800mm



YG-035 三门衣柜
Three Doors Wardrobe
980*480*1800mm



YG-22 六门衣柜
Six Doors Wardrobe
980*480*1800mm



YG-03 两门衣柜
Two Doors Wardrobe
980*480*1800mm



YG-16 四门衣柜
Four Doors Wardrobe
980*480*1800mm



YG-21 八门衣柜
Eight Doors Wardrobe
980*480*1800mm

Swimming Pool Door Locker



GY-01 18
980*400*1800mm

OFFICE PASSWORD CABINET LOCK

Moistureproof、 Enduring、 Simple management

Swimming Pool Door Locker



GY-01
980*400*1800mm




GY-02
980*400*1800mm



GY-03
980*400*1800mm



GY-04
980*400*1800mm



GY-05
980*400*1800mm



GY-06
980*400*1800mm

Doors Wardrobe



OFFICE PASSWORD CABINET LOCK

Moistureproof、 Enduring、 Simple management




980*480*1800mm

Page : 76

www.hamptonsurgical.co.uk

Doors Wardrobe




GY-07
980*480*1800mm



GY-08
980*480*1800mm



GY-09
980*480*1800mm



GY-10
980*480*1800mm



GY-11
980*480*1800mm



GY-12
980*480*1800mm

Steel Office Desk



FASION DESIGN

AND EXCELLENT QUALITY

Bring You Comfortable and Gladderghaddest Mood



Steel Office Desk



ZT-07
Hanging Drawers Office Desk
1000*500*750mm
1000*600*750mm
1200*600*750mm
1400*700*750mm



ZT-012N-1
Single Side Office Desk
1000*500*750mm
1000*600*750mm
1200*600*750mm
1400*700*750mm



ZT-010N
Double Side Office Desk
1400*700*750mm
1500*700*750mm

Steel Office Desk



ZT014H
DOUBLE PERSON COMPUTER DESK
1600*700*750MM



ZT01H
SINGLE PERSON COMPUTER DESK
1000*600*750MM
1200*700*750MM



ZT02H
SINGLE PERSON COMPUTER DESK
1500*700*750MM



ZT019H
DOUBLE PEDESTAL COMPUTER DESK
1400*700*750MM

Steel Office Desk



01 single side office desk
1200*600*750mm
1400*700*750mm
1600*700*750mm



02 pedestal office double
1200*600*750mm
1400*700*750mm
1600*700*750mm



03 single side office desk
1200*600*750mm
1400*700*750mm
1600*700*750mm



04 single side office desk
1200*600*750mm
1400*700*750mm
1600*700*750mm



05 single side office desk
1200*600*750mm
1400*700*750mm
1600*700*750mm



06 single side office desk
1200*600*750mm
1400*700*750mm
1600*700*750mm

Steel Office Desk



01 single side office desk
1600W*700D*750H

1000W*430D*750H



02 single side office desk
1600W*700D*750H

1000W*430D*700H



03 single side office desk
1600W*700D*750H

1000W*430D*750H

Library Equipment Series



SERIES PRODUCT LIST

1
SERIES



STEEL SHELF

2
SERIES



ALL STEEL SHELVES

3
SERIES



OAK GUARD A BOOKCASA

Library Equipment Series



4
SERIES



ANCIENT BOOKCASE

5
SERIES



CARREL

6
SERIES



MUTE BOOK LADDER

Library Equipment Series




SJ-01 W1800*D450*H2000



SJ-02 W1800*D450*H2000



SJ-03 W1800*D450*H2000



SJ-02 W2700*D450*H2000



SJ-03 W2700*D450*H2000



Library Equipment Series

SERIES PRODUCT LIST



QK-01
W900*D360*H1800MM



QK-02
W900*D360*H1800MM



QK-03
W1800*D650*H2000MM



YL-01
W1500*D900*H7500MM
W1800*D900*H7500MM



YL-02
W1500*D900*H7500MM
W1800*D900*H7500MM



YL-03
W1500*D900*H7500MM
W1800*D900*H7500MM



PJ-01
W750*D360*H900MM



PJ-02
W450*D650*H1300MM



PJ-03
: 160*235MM : 135*175MM

Intelligent Electric



Intelligent Electric

INTELLIGENT ELECTRIC



MJ-01



MJ-02



MJ-03



MJ-04

1		2	
1200*900*100		860*640*100	
1300*940*1160		960*740*1200	
940		740	
1300*	+135	1300*	+135
1800	2000	2200	

Intelligent Electric

INTELLIGENT ELECTRIC



MJ-05



MJ-06

TECHNICAL PARAMETERS AND STRUCTURAL ANALYSIS

500	900	2400	5	398	2	1935	2	>60KG	
550			6	330	3	2835	2		
600			7	280	4	3735	3		
650			8	205	5	4635	3		
700			12	155	6	5535	4		
800					7	6435	4		
350						8	7335	5	>80KG
400						9	8235	5	
450						10	9135	6	



Warehouse Equipment Series



INSPIRATION CHALLENGE TRADITIONAL DETAILS OF ACHIEVEMENTS OF THE FUTURE



Warehouse Equipment Series



[mm]		[mm]	
1600*600		1800	
1600*900		2000	
1600*1180		2400	
2000*600		2700	
2000*900			
1600*1180			

70*90	40*65	60*40*20	40*60	40*80	30*70

Student Apartment Bunk Bed



EVERY
DETAIL REVEALS A FINE



Student Apartment Bunk Bed

APARTMENT COMBINED BED



ZH-01
: 2000W*980D*2100H mm



ZH-02
: 2000W*980D*2100H mm



ZH-03
: 2000W*980D*2100H mm



ZH-04
: 2000W*980D*2100H mm



ZH-05
: 2000W*980D*2100H mm



ZH-06
: 2000W*980D*2100H mm

Student Apartment Bunk Bed


APARTMENT COMBINED BED



SQ-01
: 4000W*980D*2000H mm



SQ-02
: 4000W*980D*2000H mm



SQ-03
: 4000W*980D*2000H mm



SQ-04
: 4500W*980D*2000H mm



SQ-05
: 4000W*980D*2000H mm



SQ-06
: 4500W*980D*2000H mm

Student Desk & Chair



: W1200*D500*H700mm

: W380*D375*H440mm



Student Desk & Chair



KZ-01

: 650*450*700mm
: 400*360*700mm



KZ-02

: 650*450* (680-760) mm



KZ-03

: 600*400* (65-75) mm



KZ-04

: 600*400* (65-75) mm



KZ-05

: 400*600* (630-760) mm
: 380*360* (360-420) mm



KZ-06

: 400*600* (630-760) mm
: 340*240* (360-420) mm

Student Desk & Chair



KZ-07
: 1200*500*700mm
: 380*375*440mm



KZ-08
: 1200*420*730mm
: 350*410*740mm



KZ-09
: 1100*400* (630-760) mm
: 340*240* (360-420) mm



KZ-10
: 1050*400*750mm



KZ-11
: 1050*400*780mm



KZ-12
: 1200*420*730mm
: 350*410*740mm

Classroom Seating



Classroom Seating

AMPHITHEATRE CHAIR



JT-01
: 1040*840*840MM



JT-02
: 1040*840*840MM



JT-03
: 1040*840*840MM



JT-04
: 1040*840*840MM



JT-05
: 1040*840*840MM



JT-06
: 1040*840*840MM



Dining Table Chair



COLORFUL LIEF STYLE



Dining Table Chair



CZ-01
: 1430*750*1100MM



CZ-02
: 1430*750*1100MM



CZ-03
1430*750*1100MM



CZ-04
: 1430*750*1100MM



CZ-05
: 1650*750*1200MM



CZ-06
: 1650*750*1200MM



CZ-07
: 1650*750*600MM



CZ-08
: 1650*750*1200MM

Training Chair



HIGH BACK FIT CLOSELY BACKBONE CURVE, EASE SITTING BACK PRESSURE OF HUMAN BODY
ENGINEERING DESIGN, BEAUTIFUL AND COMFORTABLE AS A WHOLE



Training Chair



ZD-01
: 830*550*485MM



ZD-03
: 755*530*520MM



ZD-05
: 785*540*520MM



ZD-02
: 830*550*485MM



ZD-04
: 750*580*680MM



ZD-06
: 825*440*580MM